




JUN 30, 2023

DNA/RNA extraction from fresh-frozen tissue, AllPrep DNA/RNA/miRNA Universal Kit

Annika Fendler¹

¹Charite



Annika Fendler

OPEN ACCESS

Protocol Citation: Annika Fendler 2023. DNA/RNA extraction from fresh-frozen tissue, AllPrep DNA/RNA/miRNA Universal Kit. **protocols.io** <https://protocols.io/view/dna-rna-extraction-from-fresh-frozen-tissue-allpre-cwcqxavw>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development
We are still developing and optimizing this protocol

Created: Jun 27, 2023

Last Modified: Jun 30, 2023

PROTOCOL integer ID: 84080

Keywords: DNA, RNA, Fresh-frozen tissue, Qiagen AllPrep

ABSTRACT

Protocol for combined RNA and DNA extraction from fresh-frozen tissue using the AllPrep DNA/RNA/miRNA Universal Kit.

MATERIALS

⊗ AllPrep DNA/RNA/miRNA Universal Kit (50) Qiagen Catalog #80224

⊗ Genomic DNA ScreenTape Agilent Technologies Catalog #5067-5365

⊗ Qubit™ dsDNA BR Assay Kit Thermo Fisher Scientific Catalog #Q32853

⊗ Qubit RNA BR Assay Kit Thermo Fisher Scientific Catalog #Q10211

⊗ RNA ScreenTape and Reagents Agilent Technologies

⊗ EB buffer Qiagen Catalog #19086

β-ME

EtOH

Isoprop

1.5 and 2 ml LoBind tubes

TissueLyser Beads 5 mm

SAFETY WARNINGS

! Dnase I stocks can be used 4 weeks after being thawed but should not be frozen again

BEFORE START INSTRUCTIONS

Preparations:

FRN buffer: Add 42 ml Isoprop to new bottle

RPE buffer: Add 44 ml EtOH to new bottle

AW1 buffer: Add 25 ml EtOH to new bottle

AW2 buffer: Add 30 ml EtOH


DNase I stocks: 550 µl RNase-free water to lyophilised DNase I, aliquot and store at -20°C for 9 months)

Immediately

DNase I: 70 µl RDD + 10 µl DNase I per sample


Proteinase K: 60 µl AW1 + 20 µl Proteinase K per sample

Tissue preparation

1 This protocol is for  Sample


Optional: Weigh tissue

2 Transfer tissue in  350 µL for up to 10 mg of tissue

 600 µL for 10 to 30 mg or stabilised tissue

RLT + β-ME in 2 ml DNA LoBind tube.

Make sure that the tissue does not defrost.

3 Add a 5mm bead to each tube and lyse tissue in TissueLyser for  00:02:00 @ 20 Hz

2m

Equipment

TissueLyser II

NAME

Bead Mill

TYPE

QIAGEN









BRAND

85300



SKU
















<https://www.qiagen.com/us/products/human-id-and-forensics/automation/tissuelyser-ii/#orderinginformation>

LINK

- 4 Spin down for  00:01:00 @  9000 x g 1m
- 5 Turn tube rack and lyse tissue for  00:02:00 @ 20 Hz 2m
- 6 Spin down for  00:01:00 @  9000 x g 1m
- 7 Add lysed product to DNA Mini Spin Column
- 8 Spin  00:00:30 @  max x g, repeat if any liquid remains on column 30s
- 9 Transfer column to a new collection tube and store tube at  4 °C until DNA extraction
- 10 Transfer flow-through to new 2 ml LoBind tube

RNA extraction



- 11 Add  50 µL for up to 10 mg tissue  80 µL for 10 to 30 mg or stabilised tissue Proteinase K, mix by pipetting

- 12 Add  200 µL for up to 10 mg tissue  350 µL for 10 to 30 mg or stabilised tissue EtOH abs.,
mix by inverting, spin down liquid
- 13 Incubate  00:10:00 @  Room temperature 10m
- 14 Add  400 µL for up to 10 mg tissue  750 µL for 10 to 30 mg or stabilised tissue EtOH abs. ,
mix by pipetting
- 15 Add  700 µL to RNeasy spin column and spin for  00:00:30 @  max x g and discard 30s
flow-through
- 16 Repeat until all liquid passed through the column
- 17 Add  500 µL RPE
- 18 Spin  00:00:30 @  max x g and discard flow-through 30s
- 19 Add  80 µL DNase I working solution directly onto the membrane and incubate 15m
 00:15:00 at  Room temperature
- 20 Add  500 µL FRN

- 21 Spin  00:00:30 @  max x g and *don't* discard flow-through 30s
- 22 Add flow-through again to column
- 23 Spin  00:00:30 @  max x g and discard flow-through 30s
- 24 Add  500 μ L RPE
- 25 Spin  00:00:30 @  max x g and discard flow-through 30s
- 26 Add  500 μ L EtOH abs
- 27 Spin  00:02:00 @  max x g and place column in new collection tube 2m
- 28 Spin  00:02:00 @  max x g and place column in new 1.5 ml collection tube 2m
- 29 Add  30 μ L of RNase-free H₂O

30 Incubate  00:01:00 @  Room temperature 1m



31 Spin down  00:01:00 @  8000 x g 1m


32 QC: Qubit  Qubit RNA BR Assay Kit Thermo Fisher Scientific Catalog #Q10211 and tape
station  RNA ScreenTape and Reagents Agilent Technologies

DNA extraction

11m
















33  350 µL AW1 auf DNA Mini Spin Column geben

34 Spin  00:00:30 @  max x g and discard flow-through 30s

35  80 µL Proteinase K working solution directly onto membrane

36 Incubate  00:05:00 @  Room temperature 5m

37 Add  350 µL AW1

- 38 Spin  00:00:30 @  max x g and discard flow-through 30s
- 39 Add  350 μ L AW2
- 40 Spin  00:02:00 @  max x g and place column in new 2 ml collection tube 2m
- 41 Spin  00:01:00 @  max x g and place column in new 1.5 ml collection tube 1m
- 42 Add  50 μ L  EB buffer Qiagen Catalog #19086 directly onto membrane
- 43 Incubate  00:01:00 @  Room temperature 1m
- 44 Spin down  00:01:00 @  8000 x g 1m
- 45 QC: Qubit DNA  Qubit™ dsDNA BR Assay Kit Thermo Fisher Scientific Catalog #Q32853
and tape station  Genomic DNA ScreenTape Agilent Technologies Catalog #5067-5365