




Sep 20, 2022

# 🌐 Perfusion-fixation procedure for adult rhesus monkeys

Yoland Smith<sup>1,2</sup><sup>1</sup>Emory University; <sup>2</sup>ASAP1 *Works for me* Share[dx.doi.org/10.17504/protocols.io.5jyl8jn67g2w/v1](https://dx.doi.org/10.17504/protocols.io.5jyl8jn67g2w/v1) johnson.agniswamy

## ABSTRACT

This protocol details Perfusion-fixation procedure for adult rhesus monkeys.

## ATTACHMENTS

[527-1107.pdf](#)

## DOI

[dx.doi.org/10.17504/protocols.io.5jyl8jn67g2w/v1](https://dx.doi.org/10.17504/protocols.io.5jyl8jn67g2w/v1)

## PROTOCOL CITATION

Yoland Smith 2022. Perfusion-fixation procedure for adult rhesus monkeys.  
**protocols.io**  
<https://protocols.io/view/perfusion-fixation-procedure-for-adult-rhesus-monk-cgkutuww>



## KEYWORDS

Perfusion-fixation , Adult rhesus monkeys

## LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Sep 14, 2022

## LAST MODIFIED

Sep 20, 2022

## OWNERSHIP HISTORY

Sep 14, 2022  maria.sSep 16, 2022  johnson.agniswamy

## PROTOCOL INTEGER ID

70004

## MATERIALS TEXT

**Solutions to prepare:**

- Pentobarbital (  **25 mg/mL** solution).

- **Ringer's solution**

Dissolve the following reagents, in order, in  **1000 mL** of distilled water.

| A                               | B        |
|---------------------------------|----------|
| Reagent                         | Quantity |
| HEPES                           | 0.60 g   |
| NaCl                            | 11.86 g  |
| KCl                             | 0.223 g  |
| CaCl <sub>2</sub>               | 0.353 g  |
| NaHCO <sub>3</sub>              | 2.18 g   |
| KH <sub>2</sub> PO <sub>4</sub> | 0.177 g  |
| MgSO <sub>4</sub>               | 0.32 g   |
| D-glucose                       | 1.8 g    |

- Fixative  **2.5 L**

 **Paraformaldehyde 96% Fisher**

1. **Scientific Catalog #AC416785000** )

 **Glutaraldehyde Electron Microscopy**

2. 0.1% **Sciences Catalog #16220**

3.  **0.1 Molarity (M)** phosphate buffer,  **7.4**

Glutaraldehyde must be added to the fixative solution only at the time of

Glutaraldehyde must be added to the fixative solution only at the time of perfusion.

#### FIXATIVE SOLUTION:



PREPARE THIS SOLUTION UNDER A HOOD.

- Final volume: **1 L**
- Concentration:

| A                | B    |
|------------------|------|
| Paraformaldehyde | 4%   |
| Glutaraldehyde   | 0.1% |

- Heat **500 mL** distilled water until the temperature reaches **60 °C**.
- Dissolve **40 g** paraformaldehyde with strong stirring action for **00:20:00**.
- Add NaOH drop by drop until the solution is almost clear.
- Let the solution cool for **01:00:00**.
- Filter the solution in another beaker and add **500 mL** PB (**0.2 Molarity (M)** **pH 7.4**)

KEEP IN THE FRIDGE.

- Just before the perfusion, substitute **4 mL** of fixative by **4 mL** of glutaraldehyde 25%.

#### ANTI-FREEZE SOLUTION:

Prepare the following solution in the same order:

For **1000 mL** solution:

| A  | B      |
|--|--------|
| Sodium phosphate monobasic, NaH <sub>2</sub> PO <sub>4</sub> ·H <sub>2</sub> O             | 13.8 g |
| Sodium phosphate dibasic heptahydrate, Na <sub>2</sub> HPO <sub>4</sub> ·7H <sub>2</sub> O | 25.8 g |
| Distilled water  | 400 ml |
| Ethylene glycol  | 300 ml |
| Glycerol   | 300 ml |

Sections can be stored in this solution at **-20 °C** freezer for many months without loss of ultrastructural preservation and tissue antigenicity.

#### Equipment and room needed:

- Necropsy room equipped with proper exhaust and ventilation to reduce exposure to aldehydes.

General-Purpose Peristaltic Pumps  
Peristaltic Pump  
Fisherbrand™ GP1000 13310656 [↗](#)

- Oxygen tank (95% O<sub>2</sub>, 5%CO<sub>2</sub>).
- Surgical instruments (scalpels, hemostats of various sizes, bone cutter, scissors, bone saw, forceps, rongeur).

Wear PPE before entering the perfusion room, i.e., gloves, coat, hair net, safety glasses, face mask and face shield or goggles.



Make sure the Hazardous Waste container is not full. Replace if full.

### Perfusion-fixation procedure for adult rhesus monkeys

1d

1



Add glutaraldehyde to the fixative solution.

2

Start oxygenation of the Ringer's solution.

3

Anesthetize the animal with Ketamine ([M]10 mg/kg , i.m.) or Telazol ([M]3 mg/kg - [M]5 mg/kg , i.m.) in its home cage and bring to the necropsy room.

4

Place an i.v. catheter.

5



Perform a tracheal intubation.

6

Inject  1 mL of heparin (i.v.).

7

Inject an overdose of pentobarbital ([M]25 mg/kg i.v.).

- 8 Ensure the absence of all reflexes, including toe pinch reflexes and brainstem reflexes, such as corneal reflexes.
- 9 Cut the thoracic cage to provide access to the heart, clamp the descending aorta, open the pericardium and expose the heart muscle, make an incision of the right atrium and insert the perfusion needle in the left ventricle.
- 10 Artificially ventilate the animal through the tracheal tube during the Ringer's solution perfusion.
- 11  Infuse transcardially ~300- **400 mL** of cold oxygenated Ringer's solution at a rate of 80-90 ml/min through a needle (14G; 1.5 inch long) inserted in the left ventricle. Let the infused solution exit the vascular system through a hole in the right atrium.
- 12  Perfuse **2.5 L** of fixative, starting at a rate of 80-90 ml/min for the first liter and then reduce the rate to 50 ml/min for the remaining solution.
- 13 After perfusion, fix the animal's head in a stereotaxic frame and open the skull using a bone saw and rongeur.
- 14 Cut the brain in **10 mm**-thick blocks in the coronal stereotaxic plane, remove the resulting blocks of tissue from the skull and post-fix them for **24:00:00** in 4% paraformaldehyde solution in PB (**0.1 Molarity (M)**, **pH 7.4**) at **40 °C**.<sup>1d</sup>
- 15 Transfer the tissue into a phosphate-buffered saline (PBS, **0.01 Molarity (M)**, **pH 7.4**) solution.
- 16 Cut the brain in **50 µm**-thick sections using a vibrating microtome or a freezing microtome

and store in a  $-20\text{ }^{\circ}\text{C}$  freezer in an anti-freeze solution ( **Please refer material section**) until further processing.