

VERSION 1
MAR 30, 2023

OPEN ACCESS

DOI:
dx.doi.org/10.17504/protocols.io.n2bvj8oowgk5/v1

Protocol Citation: Antoine Champie, Amélie De Grandmaison 2023. HTTM : Illumina libraries. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.n2bvj8oowgk5/v1> version created by Antoine Champie

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Oct 27, 2022


Last Modified: Mar 30, 2023

PROTOCOL integer ID:
71943

HTTM : Illumina libraries V.1

Antoine Champie¹, Amélie De Grandmaison¹

¹Université de Sherbrooke



Antoine Champie

ABSTRACT

Part three of the HTTM protocol. A low-cost and high-throughput Tn-seq protocol. This part cover the preparation of Illumina sequencing libraries form genomic DNA.

MATERIALS

Preparation of Nextera adapters :

Nextera (Nxt) adapters are prepared by hybridization of the following primers :

A	B
Nxt-XTv2-B-N701-T	CAAGCAGAAGACGGCATACGAGATTGCGCTTAGTCTCGTGGGCTCG GAGATGTGTATAAGAGACAGT
Nxt-XTv2-B-3R-ac3-phos5'	/5Phos/CTGTCTCTTATACACATCTCCGAGCCCACGAGAC/3InvdT/

■ Preparation of the 5X annealing buffer (5X Tris NaCl buffer : 50 mM Tris, pH 7.5-8, 250 mM NaCl) :

- 500 µl Tris-HCl 1M pH 7.5
- 500 µl NaCl 5M
- 9 ml H2O mol.-grade

■ Preparation of the adapters (40 µM 50 µL) :

- Resuspend both primers in water to obtain 100 µM stocks
- Mix 20 µl of each (Nxt-XTv2-B-N701-T and Nxt-XTv2-B-3R-ac3-phos5')
- Add 10 µl of 5X annealing buffer
- Annealing reaction in a thermocycler (decrease temperature from 98 °C to 4 °C (-0.1 °C/cycle(10s/cycle)))

Primers used for the first PCR :

A	B
---	---

A	B
Nxt_A	AATGATACGGCGACCACCGAGATCTACAC
Nxt_B	CAAGCAGAAGACGGCATACGAGAT

Primers template for barcoding PCR :

A	B
Nxt_i5_barco ding	AATGATACGGCGACCACCGAGATCTACAC [8 Nu Index] TCGTCGGCAGCGTCAGATGTGTA
Nxt_i7_barco ding	CAAGCAGAAGACGGCATACGAGAT [8 Nu Index] GTCTCGTGGGCTCGGAGATGTGTATAAG

Kit used for library preparation :

NEBNext Ultra II DNA Library Prep Kit for Illumina

NEB CAT#: E7645S

PCR mix used :

Supermix 2X

Homemade

SPRI beads used:

Ampure XP DNA beads

Beckman Coulter CAT#: A63882

BEFORE START INSTRUCTIONS

- All steps and master mixes need to be kept on ice as much as possible. Thermocyclers need to be cooled at 4 °C before inserting sample plate.

Libraries

1 h 34 m

1 Transfer  2.5 µL of DNA from the DNA extraction plate to a new PCR plate.

2 Prepare a fragmentation master mix with :





A	B
NEB Ultra II FS buffer	77 µl
NEB Ultra II FS enzyme	22 µl


A	B
Molecular grade water	11 µl

3 Add  1 µL of the fragmentation master mix to each well.

4 Incubate in a thermocycler with the following protocol :


45m

-  00:15:00 at  37 °C
-  00:30:00 at  65 °C

5 Add  1 µL of 4µM Nextera (NxT) adaptors to each well.





6 Prepare a ligation master mix with :

A	B
NEB Ultra II ligation master mix	377.4 µl
NEB Ultra II ligation enhancer	12.1 µl

7 Add  3.5 µL of ligation master mix to each well.











8 Incubate in a thermocycler with the following protocol :


40m

-  00:30:00 at  20 °C
-  00:10:00 at  65 °C

9 Prepare a PCR master mix with :


A	B
NxT_A primer 20 μ M	880 μ l
Nxt_B primer 20 μ M	880 μ l
Molecular grade water	8360 μ l
PCR Supermix 2X	11000 μ l

- 10** Add  192 μ L of PCR master mix to each well.
- 11** Split the PCR reaction into 4 different plates (50 μ l per plate).
- 12** Incubate each plate in a thermocycler with the following cycles : 3m 15s
-  00:00:30 at  98 °C
 -  00:00:15 at  98 °C
 -  00:00:30 at  72 °C
 - Repeat from step 2 for 20~25 cycles*
 -  00:02:00  72 °C
- 13** Pool the 4 PCR replicates together in a PCR plate.
- 14** Transfer  2 μ L of DNA from the pool plate to a new PCR plate.

- 15 Add  2 μL of each barcoding primer to the DNA :
- Nxt_i5_barcoding
 - Nxt_i7_barcoding









- 16 Prepare a PCR master mix with :

A	B
Molecular grade water	2090 μl
PCR supermix 2X	2750 μl


- 17 Add  44 μL of the PCR master mix to each well of the plate.

- 18 Incubate in a thermocycler using the following protocol :

3m 45s

-  00:00:30 at  98 °C
-  00:00:15 at  98 °C
-  00:01:00 at  72 °C (no anneal step)
- Repeat from step 2 for 5 cycles
-  00:02:00 at  72 °C

- 19 Pool together  2 μL of each sample.

- 20 Purify with Ampure XP SPRI beads using a 0.8 ratio. Resuspend with  50 μL of molecular grade water.



21 Proceed with QC and sequencing.