

Jul 10, 2024

Immunohistochemistry on mouse brain sections

DOI

dx.doi.org/10.17504/protocols.io.261ge3yywl47/v1

Shiyi Wang¹

¹Duke University

ASAP Collaborative Rese...



Shiyi Wang

Duke University

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.261ge3yywl47/v1

Document Citation: Shiyi Wang 2024. Immunohistochemistry on mouse brain sections. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.261ge3yywl47/v1>

License: This is an open access document distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Created: May 23, 2023

Last Modified: July 10, 2024

Document Integer ID: 82330

Keywords: ASAPCRN

Funders Acknowledgement:

Aligning Science Across

Parkinson's (ASAP) initiative

Grant ID: ASAP-020607



Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to **protocols.io** is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with **protocols.io**, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Abstract

Immunohistochemistry on mouse brain sections



1. Mice used for immunohistochemistry were anesthetized with 200 mg/kg Avertin and perfused with TBS/Heparin and 4% PFA.
2. Brains were collected and post-fixed in 4% PFA overnight, cryoprotected in 30% sucrose, frozen in a solution containing 2 parts 30% sucrose and 1-part O.C.T. (TissueTek), and stored at -80°C .
3. Floating coronal tissue sections of 30 μm , 40 μm or 100 μm thickness were collected and stored in a 1:1 mixture of TBS/glycerol at -20°C .
4. For immunostaining, sections were washed in $1\times$ TBS containing 0.2% Triton X-100 (TBST), blocked in 10% NGS diluted in TBST, and incubated in primary antibody for 2-3 nights at 4°C with gentle shaking.
5. Primary antibodies used were anti-LRRK2 (Rabbit, 1:500; ab133474, Abcam), phospho-ERM (Rabbit, 1:500; #3141, Cell Signaling), Sox9 (Rabbit, 1:500; AB5535, Millipore), GFAP (Rabbit, 1:500; Z0334, Agilent DAKO), VGluT1 (Guinea pig, 1:2000; 135304, Synaptic Systems), PSD95 (Rabbit, 1:300; 51-6900, Invitrogen), VGAT (Guinea pig, 1:1000; 131004, Synaptic Systems), and GEHPYRIN (Rabbit, 1:1000; #14304S, Cell Signaling).
6. Following the primary incubation, sections were washed in TBST, incubated in Alexa Fluor conjugated secondary antibodies diluted 1:200 (Life Technologies) for 2-3 hours at room temperature, washed with TBST, and mounted onto glass slides using a homemade mounting media (90% Glycerol, 20 mM Tris pH 8.0, 0.5% n-Propyl gallate) and sealed with nail polish.
7. For DAPI staining, DAPI (1:50,000) was added to the secondary antibody solution for the final 10 minutes of incubation.
8. Images were acquired with an Olympus FV 3000 microscope.