



Sep 11, 2020

# Isolation, Identification, and Antimicrobial Susceptibility Testing of *Streptococcus pneumoniae*

Korrie Salsabila<sup>1</sup>, Wisnu Tafoji<sup>1</sup>, Wisiva Tofriska Paramaiswari<sup>1</sup>, Miftahuddin Majid Khoeri<sup>1</sup>, Dodi Safari<sup>1</sup>

<sup>1</sup>Eijkman Institute for Molecular Biology

1 Works for me [dx.doi.org/10.17504/protocols.io.bkuakwse](https://dx.doi.org/10.17504/protocols.io.bkuakwse)

Eijkman Institute for Molecular Biology



Korrie Salsabila  
Eijkman Institute for Molecular Biology

## ABSTRACT

*Streptococcus pneumoniae* (*S. pneumoniae*) is a major cause of bacterial pneumonia and meningitis globally. This bacterium is a normal flora of nasopharynx but can move to sterile sites such as cerebrospinal fluid, blood, and pleural fluid then causes invasive pneumococcal diseases. However, antibiotics resistance causes the failure of pneumococcal disease treatment. *S. pneumoniae* penicillin-resistant is one of the priority antimicrobial resistant (AMR) pathogens in WHO list. Since pneumococcus is a fastidious bacteria, isolation and identification of *S. pneumoniae* is challenging. This bacteria need good quality of blood agar plate to grow and proper method to identify *S. pneumoniae*. Performing antimicrobial susceptibility testing of *S. pneumoniae* isolates is important to monitor and handle its resistance. Therefore, we decided to write the protocol regarding isolation, identification, and antimicrobial susceptibility testing of *S. pneumoniae*.

## DOI

[dx.doi.org/10.17504/protocols.io.bkuakwse](https://dx.doi.org/10.17504/protocols.io.bkuakwse)

## PROTOCOL CITATION

Korrie Salsabila, Wisnu Tafoji, Wisiva Tofriska Paramaiswari, Miftahuddin Majid Khoeri, Dodi Safari 2020. Isolation, Identification, and Antimicrobial Susceptibility Testing of *Streptococcus pneumoniae*.

**protocols.io**

<https://dx.doi.org/10.17504/protocols.io.bkuakwse>

## KEYWORDS

*Streptococcus pneumoniae*, isolation, identification, antimicrobial

## LICENSE

———— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Sep 04, 2020

## LAST MODIFIED

Sep 11, 2020

## PROTOCOL INTEGER ID

41570

## MATERIALS

NAME	CATALOG #	VENDOR
TSA II Trypticase Soy Agar Modified	212305	Becton-Dickinson

NAME	CATALOG #	VENDOR
DD1 Optochin	X3483A	Oxoid Microbiology Products - Thermo Fischer
Sodium deoxycholate	SRE0046-100G	Sigma
Mueller-Hinton Agar	CM0337	Oxoid Microbiology Products - Thermo Fischer
Tris-EDTA buffer solution	933283	Sigma – Aldrich
Skim Milk	232100	Becton-Dickinson
Tryptone Soya Broth	CM0129	Oxoid Microbiology Products - Thermo Fischer
Dextrose	215530	Becton-Dickinson
Glycerol	G5516	Sigma
McFarland Equivalence Turbidity Standard	R20421	
Clindamycin	CT0064B	Oxoid Microbiology Products - Thermo Fischer
Erythromycin	CT0020B	Oxoid Microbiology Products - Thermo Fischer
Oxacillin	CT0159B	Oxoid Microbiology Products - Thermo Fischer
Chloramphenicol	CT0013B	Oxoid Microbiology Products - Thermo Fischer
Tetracycline	CT0054B	Oxoid Microbiology Products - Thermo Fischer
Sulphamethoxazole/Trimethoprim	CT0052B	Oxoid Microbiology Products - Thermo Fischer
Gentamicin Solution	G1397	Sigma
FLOQ Swab	503CS01	Copan

## Nasopharyngeal Swab Collection

### 1 Nasopharyngeal swab collection follow method as described previously



Satzke C, Turner P, Virolainen-Julkunen A, Adrian PV, Antonio M, Hare KM, et al. (2013). Standard method for detecting upper respiratory carriage of *Streptococcus pneumoniae*: Updated recommendations from the World Health Organization Pneumococcal Carriage Working Group. *Vaccine*. 2013 Dec;32(1):165–79..  
<http://dx.doi.org/10.1016/j.vaccine.2013.08.062>

. A flexible nasopharyngeal flocked swab is slowly passed into nasopharynx. The swab should not be resistant before it was reached nasopharynx (NP). The nasopharynx distanced one-half to two-third from nostril to ear lobe

### 2 Once the swab reaches nasopharynx, it is rotated 3-5 times to let the mucus absorb into the swab

### 3 The swab is removed slowly and immediately placed into 1 mL of Skim milk Tryptone Glucose and Glycerol (STGG). STGG **100 mL** was made by mixing **2 g** skim milk powder, **3 g** Tryptone Soya Broth, **0.5 g** glucose/dextrose, glycerol **10 mL** and distilled water **90 mL**

### 4 The NP-STGG is vortexed for 10-20 sec to disperse the bacteria from swab to the STGG then placed it in cryobox that

filled by ice packs. The temperature should be  $\pm 4^{\circ}\text{C}$  and it is regularly monitored through the digital temperature controller

- 5 The inoculated STGG is kept at  $\pm 4^{\circ}\text{C}$  for 4-6 hours after collection until reaches targeted laboratory. After that, the swab-STGG is frozen at  $\pm -80^{\circ}\text{C}$  to secure survival of pneumococcus
- 6 The specimens are shipped to Molecular Bacteriology Unit, Eijkman Institute using dry ice.
- 7 Once it is arrived Eijkman Institute, the swabs were immediately stored at  $\pm -80^{\circ}\text{C}$  until culture.

#### Optochin Susceptibility Test

- 8 NP-STGG medium that stored in  $\pm -80^{\circ}\text{C}$  is thawed and vortexed to disperse the bacteria from the swab.
- 9  $20\ \mu\text{L}$  of swab-inoculated STGG was transferred and streaked onto blood agar plate (TSA II Trypticase Soy Agar Modified) with  $8\% \text{ (v/v)}$  sheep blood and supplemented with 5 mg/L gentamicin.



Hadinegoro SR, Prayitno A, Khoeri MM, Djelantik IGG, Dewi NE, Indriyani SAK, et al.. Nasopharyngeal Carriage of *Streptococcus pneumoniae* in Healthy Children Under Five Years Old in Central Lombok Regency, Indonesia.. Southeast Asian J Trop Med Public Health. 2016;47(3):9..

- 10 The plate is incubated at  $\pm 37^{\circ}\text{C}$  with 5%  $\text{CO}_2$  for 18-24 hours  $0\ \text{Mass Percent}$
- 11 Single colony of alpha-hemolytic and flat-depressed colony resembling pneumococci is picked and streaked onto blood agar plate (TSA II Trypticase Soy Agar Modified) with  $8\% \text{ (v/v)}$  sheep blood.
- 12 Quadrant 1 and 2 is overlapped to test optochin susceptibility of the isolate. If there are 2 different colonies in the plate, each of them should be tested for optochin susceptibility separately.
- 13 Optochin disk (ethylhydrocupreine hydrochloride) is placed on the overlap area between quadrant 1 and 2.

- 14 Susceptible isolate (zone  $\geq 14$  mm) with distinct characteristics of pneumococci (flat depressed center, shiny or wet colonies), is defined as *Streptococcus pneumoniae* (*S. pneumoniae*)
- 15 Overnight culture is harvested and put into STGG then stored in  $-80^{\circ}\text{C}$

#### Bile Solubility Test

- 16 Resistant optochin isolate (zone  $\leq 14$  mm) which has characteristics of *S. pneumoniae* or susceptible isolate which has no pneumococcus distinct characteristics (too dry colonies), is confined to bile solubility testing



<https://www.cdc.gov/meningitis/lab-manual/chpt08-id-characterization-streppneumo.pdf>

- 17 For each isolate tested, two glass tubes are prepared and labeled as 'test tube' and 'control tube'
- 18 2 mL and 1 mL of saline 0.85% are added into test tube and control tube, respectively
- 19 Bacterial suspension that equivalent to MacFarland standard 1 is made in the test tube. Bacterial density is measured by densitometer (make sure to calibrate prior to using)
- 20 1 mL of bacterial suspension from test tube is transferred into control tube then vortexed.
- 21  $1\text{ mL}$  of bile salt or sodium deoxycholate  $2\text{ Mass Percent}$  is added into test tube. The final volume of test tube and control tube are  $2\text{ mL}$  for each
- 22 ATCC 49619 of *S. pneumoniae* must be included in every bile testing
- 23 Test tube and control tube are incubated room temperature for 10 minutes. If the solution in the test tube of isolate is as clear as the test tube of ATCC 49619 (the solution should be as clear as water), the isolate is bile-soluble and defined as *S. pneumoniae*. Otherwise, if there is no turbidity difference between test tube and control tube of the isolate, the isolate was not bile-soluble
- 24 If the isolate is negative for the first 10 minutes incubation, then continue for incubation at room temperature for 2 hours
- 25 Bile soluble isolate is streaked onto blood agar plate (TSA II Trypticase Soy Agar Modified) with  $8\% \text{ (v/v)}$  sheep

blood and incubated at  $37^{\circ}\text{C}$  with 5%  $\text{CO}_2$  for 18-24 hours

- 26 Bile-soluble isolates were harvested and put into labeled STGG for freezing stock then vortexed. The bacteria should be well-homogenized in STGG
- 27 STGG-inoculated isolate is stored in  $-80^{\circ}\text{C}$

#### DNA Extraction of *S. pneumoniae*

- 28 The DNA extraction of *S. pneumoniae* is performed by boil method as described previously



Pai R, Gertz RE, Beall B. Sequential Multiplex PCR Approach for Determining Capsular Serotypes of Streptococcus pneumoniae Isolates. J Clin Microbiol. 2006;44:8.. <http://10.1128/JCM.44.1.124-131.2006>

Pure colonies of *S. pneumoniae* stored in STGG at  $-80^{\circ}\text{C}$  are streaked onto 8% sheep blood agar plate. Inoculated STGG must be kept at  $4^{\circ}\text{C}$

- 29 Incubation at  $37^{\circ}\text{C}$  for 18-24 h with 5%  $\text{CO}_2$ .
- 30 Single colony of the isolate is picked and streaked onto 8% sheep blood agar plate then incubated for 18-24 h at  $37^{\circ}\text{C}$  with 5%  $\text{CO}_2$ .
- 31 Fresh culture is harvested into TE buffer in  $1.5\text{ mL}$  microcentrifuge tube and vortexed until homogenized
- 32 Bacterial suspension is heated at  $100^{\circ}\text{C}$  for  $00:05:00$  and immediately placed at freezer  $-20^{\circ}\text{C}$   $00:05:00$
- 33 Centrifuge at  $13000 \times g$ ,  $00:10:00$  then kept at  $-20^{\circ}\text{C}$  until serotyping

#### Serotyping of *S. pneumoniae*

- 34 Serotyping of *S. pneumoniae* was performed by sequential multiplex conventional PCR targeted *wzy* gene and consisted of 8 reactions



Carvalho M d. G, Pimenta FC, Jackson D, Roundtree A, Ahmad Y, Millar EV, et al. (2010). Revisiting Pneumococcal Carriage by Use of Broth Enrichment and PCR Techniques for Enhanced Detection of Carriage and Serotypes. J Clin Microbiol. 2010 May 1;48(5):1611–8.  
<http://JOURNAL OF CLINICAL MICROBIOLOGY, May 2010, p. 1611–1618>

. Latin America scheme as described by Centers for Disease Control and Prevention, is followed for serotyping the isolates [Latin America Scheme Multiplex PCR.pdf](#). Every reaction should be detect *cpsA* gene. Each reaction has different serotypes to be identified [PCR Oligonucleotide Primers.pdf](#).



Reaction 1: serotype 14, 6A/6B/6C/6D, 23F, 19A, and 9V/9A.  
 Reaction 6C: serotype 6A/6B/6C/6D and 6C/6D.  
 Reaction 2: serotype 19F, 3, 15B/15C, 18C/18F/18B/18A, and 17F.  
 Reaction 3: serotype 1, 5, 9N/9L, 7F/7A, and 16F.  
 Reaction 4: serotype 8, 2, 4, 20, and 22F/22A.  
 Reaction 5: serotype 7C/7B/40, 12F/12A/44/46, 11A/11D, 10A, and 23A.  
 Reaction 6: serotype 21, 33F/33A/37, 15A/15F, 35F/47F, and 13.  
 Reaction 7: serotype 39, 23B, 35A/35C/42, 38/25F/25A and 35B.  
 Reaction 8: serotype 24F/24A/24B, 10F/10C/33C, 34, 19Fvar, and 31.

35 Mastermix for each reaction (25 µl) is prepared as described in Latin America scheme spreadsheet. For example, reaction 1 is performed by adding PCR H<sub>2</sub>O [0.5 µl](#) of 5× PCR buffer (Promega), [3.5 µl](#) of MgCl<sub>2</sub> 25 mM, [1 µl](#) of dNTPs (Promega) 5 mM, varied volume for 25 µM forward primer and 25 µM reverse primer, [0.2 µl](#) Taq Polymerase 2U and [2.5 µl](#) of DNA template, into 1.5 mL of microcentrifuge tube.

36 For isolates serotyping, conventional PCR is run under following condition: pre-denaturation at [94 °C](#) for [00:04:00](#), followed by 30 cycles of denaturation at [94 °C](#) for [00:00:45](#), annealing at [54 °C](#) for [00:00:45](#), and extension at [65 °C](#) for [00:02:30](#).

37 The DNA is visualized through electrophoresis. Agarose [2 Mass Percent](#) is made by adding 2 gram of agarose into 100 mL of Tris Acetat EDTA (TAE) 1× then heated in the microwave.

- 38 Agarose is poured into electrophoresis tank with the comb.
- 39 Agarose gel is placed into electrophoresis tank filled with cold TAE 1×
- 40 PCR product and DNA ladder 50 bp are added into each well of agarose gel
- 41 Voltage is set to 100V and run for 90 minutes
- 42 DNA bands are visualized using GelDoc machine
- 43 Isolate that serotype-negative is confirmed by real-time PCR targeting *lytA* to confirm that the isolate was *S. pneumoniae*. For isolate with positive *lytA* is assigned as non-typeable (NT) by PCR *S. pneumoniae*



Carvalho M d. GS, Tondella ML, McCaustland K, Weidlich L, McGee L, Mayer LW, et al (2007). Evaluation and Improvement of Real-Time PCR Assays Targeting *lytA*, *ply*, and *psaA* Genes for Detection of Pneumococcal DNA. J Clin Microbiol. 2007 Aug 1;45(8):2460–6.  
<http://10.1128/JCM.02498-06>

#### Antimicrobial Susceptibility Testing

- 44 Antimicrobial susceptibility testing is performed by following Clinical Laboratory Standard Institutes (CLSI) 2017 guidelines



Clinical and Laboratory Standards Institute. M100 Performance standards for antimicrobial susceptibility testing 27th Edition (2017).

. Frozen isolate in STGG was streaked onto **10.1128/JCM.02498-06** 8 % (v/v) sheep blood agar plate. The STGG must be kept at **4 °C** condition.

- 45 Single colony of pneumococci is picked and cultured onto 8% sheep blood agar plate then incubated for 18-24h at 37°C with 5% CO<sub>2</sub>.
- 46 Bacterial suspension is made by harvesting freshly grown *S. pneumoniae* and put into 5 mL saline

#### [M]0.85 Mass Percent

- 47 The suspension is adjusted to be equivalent as McFarland standard 0.5. Densitometer is used to measure the turbidity. Calibration using McFarland standards is done prior to measurement.
- 48 When the suspension reaches 0.5 McFarland standard, sterile cotton swab is dipped into bacterial suspension then pressed gently on the wall of glass tube to remove excess fluid.
- 49 The cotton swab is lawned confluent onto Mueller Hinton Agar (MHA) plate with [M]5 % (v/v) sheep blood and streaked 3 times by rotating in 60° direction.
- 50 Antibiotic disks are placed onto the agar and pressed gently by sterile pinset.



Clindamycin 2 µg, chloramphenicol 30 µg, oxacillin 1 µg, tetracycline 30 µg, erythromycin 15 µg, and trimethoprim/sulfamethoxazole 1,25/23.75 µg, are used. Oxacillin disk is applied to measure susceptibility of isolates to penicillin.

- 51 Inoculated MHA is incubated at 37 °C with 5% CO<sub>2</sub> for 20-24 h.
- 52 Inhibition zone is measured and recorded. Interpretation of clear zone follow Clinical Laboratory Standards Institute guideline 2017.