



DEC 27, 2022

OPEN ACCESS

DOI:

[dx.doi.org/10.17504/protocols.io.n2bvj88nwgk5/v1](https://doi.org/10.17504/protocols.io.n2bvj88nwgk5/v1)

Protocol Citation: Sarah M Prostack, Edgar M Medina, Erik Kalinka, Lillian Fritz-Laylin 2022. Protocol 1: Electroporation of *Agrobacterium tumefaciens* with a plasmid of interest.

protocols.io
<https://dx.doi.org/10.17504/protocols.io.n2bvj88nwgk5/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: Nov 02, 2022

Last Modified: Dec 27, 2022

PROTOCOL integer ID:
72187

Keywords: Electroporation, *Agrobacterium tumefaciens*, *Spizellomyces punctatus*, chytrid fungi

Protocol 1: Electroporation of *Agrobacterium tumefaciens* with a plasmid of interest

In 1 collection

Sarah M Prostack¹, Edgar M Medina¹, Erik Kalinka¹, Lillian Fritz-Laylin¹

¹University of Massachusetts at Amherst



fritzlaylinlab.umass

ABSTRACT

Electroporation is a widespread method of transforming competent *Agrobacterium tumefaciens* (*Agro*) cells with a plasmid containing a T-DNA of interest. The resulting *Agro* can be used to transform various plants and fungi, resulting in transformed cell lines. This protocol outlines the standard electroporation protocol we use to transform *Agro* in preparation for *Agrobacterium*-mediated transformation of the chytrid fungus *Spizellomyces punctatus*.

ATTACHMENTS

[Spizellomyces transformation steps.pdf](#)

GUIDELINES

Any binary plasmid that works in *Agrobacterium tumefaciens* strain EHA105 (GoldBio #CC-225-5x50) can be used for this procedure. We have had success using plasmids derived from pPZP201-BK (PMID: 7919218). It is imperative that all steps be carried out at 4°C up until electroporation. Ensure proper sterile technique throughout this protocol; perform all steps except centrifugation and electroporation charge delivery in a laminar flow hood or in the sterile area around an open flame.

Materials:

- Use a fresh streak of *Agrobacterium tumefaciens* EHA105 (GoldBio #CC-225-5x50) to spread a lawn into a LB plate, (a new streak from a -80°C stock will take ~48 hours to display visible growth at 28°C, plan accordingly).
- LB agar plates (1.5% w/v) with and without selection antibiotics, sterile (see recipe)
- Purified plasmid(s) of interest resuspended in molecular biology grade water
- Molecular biology grade water, sterile such as MilliQ or equivalent
-  Invitrogen™ S.O.C. Medium **Fisher Scientific Catalog #15-544-034**
- 10% (v/v) Glycerol, sterile (
 -  Glycerol (Certified ACS) Fisher Chemical™ **Fisher Scientific Catalog #G33-1**)
- 1.5 mL centrifuge tubes, sterile (such as
 -  Fisherbrand™ Low-Retention Microcentrifuge Tubes **Fisher Scientific Catalog #02-681-331**
)
- 0.5 mL centrifuge tubes, sterile (such as
 -  Fisherbrand™ Snap-Cap™ Flat-Top Graduated Microcentrifuge Tubes **Fisher Scientific Catalog #02-681-268**
)
- 2 mm electroporation cuvettes, sterile (
 -  Pkg of 50 individually wrapped 2mm electroporation cuvettes **Bulldog Bio Catalog #12358-346**
)
- Culture tubes, sterile (such as
 -  VWR® Culture Tubes Plastic with Dual-Position Caps **VWR Avantor Catalog #60818-703**
)
- 100-1,000 µL micropipette such as
 -  Eppendorf Research Plus Single Channel pipette 100-1000 uL blue operating button for use with 100 **pipette.com Catalog #3123000063 ES-1000**
- 20-200 µL micropipette such as
 -  Eppendorf Research Plus single channel pipette 2-20 uL yellow operating button for use with 20 uL **pipette.com Catalog #3123000039**
- 0.1-2.5 µL micropipette such as
 -  Eppendorf Research Plus single channel pipette 0.1-2.5 uL dark grey operating button for use with **pipette.com Catalog #3123000012**
- Filter tips for the micropipettes, sterile such as
 -  TIPONE® FILTER TIPS **USA Scientific Catalog #1122-1830**
- 5 mm Glass beads, sterile
- Ice bucket with ice
- Centrifuge capable of cooling to  4 °C
- Laminar flow hood and/or open flame, for maintaining sterility.
- 70% (v/v) ethanol for maintaining sterility (if using laminar flow hood)
- Exponential decay electroporator such as Gene Pulser Xcell (
 -  Gene Pulser Xcell Total System **Contributed by users Catalog ##1652660**)
- Shaking incubator at  28 °C

BEFORE START INSTRUCTIONS

Electroporation should occur at least 4 days prior to the intended *Spizellomyces* transformation time to ensure that active, single colonies are available to be selected and transferred to liquid culture the night before transformation day (see Protocol "Growing liquid cultures of *Agrobacterium* prior to transformation day").

ATTACHMENTS


[Spizellomyces transformation steps.pdf](#)



4h 25m

Steps

20m

1 Cool the following materials  On ice at least  00:20:00 prior to starting:



1. Plate with a lawn of wild-type *Agrobacterium* grown overnight at  28 °C .
2. Purified plasmid(s) of interest.
3. Molecular biology grade water, sterile.
4. 10% (v/v) glycerol, sterile.
5. 1.5 mL centrifuge tube(s), enough to hold the volume of Agro harvested.
6. 0.5 mL centrifuge tube(s), one per plasmid to be transformed, plus controls.
7. 2 mm electroporation cuvettes, one per plasmid to be transformed, plus controls.



2 Add  1 mL to  2 mL of ice-cold water to the plate of *Agrobacterium*. Hold the plate at ~45 degrees and run the water over the surface at least 3 times, gently scraping along the agar if necessary to recover the lawn of bacteria.



Note



Try to not drag too many big clumps.
The resulting harvest should have the consistency, color and density of whey.

3 Transfer the  1 mL of harvested cells to a 1.5 mL centrifuge tube, immediately place back  On ice .



4 Pellet the cells at  4000 rcf, 00:05:00 in a rotor prechilled to  4 °C .




Note

If in a pinch, a rotor for a centrifuge without cooling capabilities can be stored  Overnight at  4 °C

5m

or for  00:10:00 at  -20 °C .



- 5 Remove the supernatant and gently resuspend the cells in  1 mL of water. Do not vortex.

Note

Keep cells  On ice when not in use.

- 6 Repeat steps 4 and 5, 2 more times for a total of 3 washes.





- 7 Remove water and resuspend the cells in  800 μ L of cold 10% (v/v) glycerol, place tubes back  On ice .

- 8 Add  50 μ L of cells to new 0.5 mL tubes, place back  On ice .




Note

Use one 0.5 mL tube for each plasmid to be transformed, plus more for planned controls.

- 9 Add  1 μ L of the plasmid of interest ( 200 ng/ μ L to  300 ng/ μ L ; directly from a miniprep should work) to its appropriate tube of cells, place back  On ice and mix gently by pipetting .





Note

Use  1 μ L of water as a negative control.

- 10 Transfer the cells to cold 2 mm electroporation cuvettes, place back  On ice .

Note

Again, use one cuvette for each control or plasmid to be transformed individually.

11 Prepare the recovery media in advance by placing  150 μ L of  Room temperature SOC medium to one 15 mL culture tube for each electroporation.

12 Turn on the electroporator and create an electroporation program with the following settings:


1. Voltage= 2400V
2. Capacitance = 25 μ F
3. Resistance = 200 Ω
4. Cuvette size = 2 mm

13 Fully dry the cuvette before placing it into the electroporator chamber.

Note



Failure to fully dry the cuvette will lead to current arcing and improper electroporation.


14 Electroporate your cuvette.

15 As quickly but as gently as possible, remove a little less than  150 μ L of SOC medium from the appropriate tube for the plasmid and pipette it into the cuvette. Remove the full volume from the cuvette and return it to the original culture tube.

Note

Do this by a flame and with good sterile technique.

The best way to do this is to set a p200 to  150 μ L, but to not pull the entire volume, leaving enough space left to fit the  50 μ L that is in the cuvette.

16 Incubate culture tubes at  28 $^{\circ}$ C, shaking at  225 rpm, 04:00:00.

4h




Note

Meanwhile, place the appropriate number of LB plates with and without selection antibiotics to pre-warm


at  28 °C .

17 Add 4-6 sterile glass beads to each LB plate.


18 Add  10 µL cells to the appropriate plates.




Note

To make spreading this small volume easier, add  40 µL sterile water to the middle of the plate before adding the cells.

The electroporation efficiency for this protocol is very high: $\geq 1.6 \times 10^8$ cfu/µg pCAMBIA1391z DNA (GoldBio) or 1.25×10^5 cfu/µg plasmid (PMID: 29487777).

Do not add more than  10 µL of cells or you risk overgrowth and a lack of individual colonies.

19 Seal and invert the plates and incubate them at  28 °C for about 4 days.



Note

Colonies should appear within 4 days. If colonies of appreciable size (2-3 mm) appear earlier than that, continue on with Sp transformation.

Grow any colonies of interest in liquid media and freeze 25% glycerol stocks to avoid needing to re-electroporate Agro.