







Feb 25, 2022

## © Quick Protocol for DNA Cleanup and Concentration Using the Monarch® PCR & DNA Cleanup Kit (5 μg) (NEB #T1030) V.6

New England Biolabs<sup>1</sup>

<sup>1</sup>New England Biolabs





dx.doi.org/10.17504/protocols.io.bg9rjz56

New England Biolabs (NEB)

Tech. support phone: +1(800)632-7799 email: info@neb.com



This is the "quick" version of Monarch® PCR & DNA Cleanup Kit (5  $\mu$ g) Protocol (NEB #T1030) DNA Cleanup and Concentration. For the full protocol, please click <u>here</u>.

DOI

dx.doi.org/10.17504/protocols.io.bg9rjz56

https://www.neb.com/protocols/2015/12/08/quick-protocol-for-monarch-pcr-dnacleanup-kit-5-g-t1030

New England Biolabs 2022. Quick Protocol for DNA Cleanup and Concentration Using the Monarch® PCR & DNA Cleanup Kit (5 µg) (NEB #T1030). **protocols.io** https://dx.doi.org/10.17504/protocols.io.bg9rjz56 Isabel Gautreau

Monarch PCR Cleanup Kit, Monarch DNA cleanup kit

\_\_\_\_\_ protocol,

Jun 07, 2020

Feb 25, 2022

37905



1

Citation: New England Biolabs Quick Protocol for DNA Cleanup and Concentration Using the Monarchî PCR & DNA Cleanup Kit (5 þg) (NEB #T1030) https://dx.doi.org/10.17504/protocols.io.bq9rjz56

For detailed protocol and more information, visit www.neb.com/T1030

The full protocol is available here.

The video protocol is available <u>here</u>.

There are two protocols available for this product:

**DNA Cleanup and Concentration (below):** for the purification of **up to**  $\sqsubseteq$ **5**  $\mu$ **g** of **DNA** (ssDNA > 200 nt and dsDNA > 50 bp) from PCR and other enzymatic reactions.

## **MATERIALS**

**⊠** Monarch® PCR & DNA Cleanup Kit (5 μg) **New England** 

Biolabs Catalog #T1030

For hazard information and safety warnings, please refer to the SDS (Safety Data Sheet).

Add isopropanol to Monarch DNA Cleanup Binding Buffer prior to use\*:

- For the 50-prep kit, add **□14 mL** of isopropanol to the DNA Cleanup Binding Buffer.
- For the 250-prep kit, add **□63.5 mL** of isopropanol to the DNA Cleanup Binding Buffer.

Add ethanol to Monarch DNA Wash Buffer prior to use (4 volumes of  $\geq$  95% ethanol per volume of Monarch DNA Wash Buffer)

- For the 50-prep kit, add **20 mL** of ethanol to the Monarch DNA Wash Buffer
- For the 250-prep kit, add ■100 mL of ethanol to the Monarch DNA Wash Buffer

Always keep all buffer bottles tightly closed when not in use.

All centrifugation steps should be carried out at **316.000** x g (~ **13.000** rpm).



Dilute sample with DNA Cleanup Binding Buffer according to the table below. Mix well by pipetting up and down or flicking the tube. Do not vortex.

Α	В	С
Sample Type	Ratio of Binding Buffer:	Example
	Sample	
dsDNA > 2 kb (plasmids, gDNA)	2:1	200 μl: 100 μl
dsDNA < 2 kb (some amplicons,	5:1	500 μl: 100 μl
fragments)		
ssDNA > 200 nt*	7:1	700 μl: 100 μl

<sup>\*</sup>Please note that recovery of ssDNA < 200 nts can be increased by using the <u>Oligonucleotide Cleanup Protocol</u>, but doing so will shift the cutoff size for DNA binding to 18 nt (versus 50 nt).

A sample volume of  $20-100~\mu l$  is recommended. For smaller samples, TE can be used to adjust the volume. For diluted samples larger than  $800~\mu l$ , load a portion of the sample, proceed with step 2, and then repeat as necessary.

## 2

Insert column into collection tube and load sample onto column. Spin at **316000** x g for **00:01:00**, then discard flow-through.



Re-insert column into collection tube. Add  $\ \Box 200 \ \mu L$  DNA Wash Buffer and spin at  $\ @16000 \ x \ g$  for  $\ @00:01:00$ .

Discarding flow-through is optional.



Repeat Step 3: Re-insert column into collection tube. Add ■200 µL DNA Wash Buffer



and spin at @16000 x g for @00:01:00.

Discarding flow-through is optional.

5 Transfer column to a clean 1.5 ml microfuge tube.

Use care to ensure that the tip of the column does not come into contact with the flow-through. If in doubt, re-spin for © 00:01:00.

6

Add  $\geq -6 \,\mu$ L DNA Elution Buffer to the center of the matrix. Wait for  $\bigcirc 00:01:00$ .

7 Spin for **© 00:01:00** to elute DNA.

Typical elution volumes are  $\Box 6~\mu L - \Box 20~\mu L$ . Nuclease-free water (pH 7–8.5) can also be used to elute the DNA. Yield may slightly increase if a larger volume of DNA Elution Buffer is used, but the DNA will be less concentrated. For larger size DNA ( $\geq 10~kb$ ), heating the elution buffer to & 50~°C prior to use can improve yield. Care should be used to ensure the elution buffer is delivered onto the matrix and not the wall of the column to maximize elution efficiency.