



VERSION 2

FEB 27, 2024

OPEN ACCESS



Protocol Citation: rebecca.bennett, Nimalka Weerasuriya 2024.
General Media Recipes:
WA, PDA, LB, CMA, V8 .
protocols.io
<https://protocols.io/view/general-media-recipes-wa-pda-lb-cma-v8-c9t5z6q6>Version created
by Nimalka Weerasuriya

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Feb 27, 2024

Last Modified: Feb 27, 2024

🌐 General Media Recipes: WA, PDA, LB, CMA, V8 V.2

rebecca.bennett¹, Nimalka Weerasuriya²

¹USDA ARS; ²USDA

Peanut and Small Grain Research Unit, Stillwater, OK



Nimalka Weerasuriya
USDA, Oklahoma State University

DISCLAIMER

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

ABSTRACT

These are current general media recipes for (acidified) water agar, (acidified) potato dextrose agar, LB broth, cornmeal agar, P5ARP/P10ARP, and clarified V8 juice agar.

GUIDELINES

These are media recipes.

SAFETY WARNINGS





Use proper PPE for toxic additives in the P5ARP/P10ARP recipes inside a fume hood. Wear gloves, eye protection, and dust mask. Terraclor May cause sensitization by skin contact. Experimental carcinogen and neoplastigen.

ETHICS STATEMENT

NA

Prepare the following:

- Turn on water bath if making more than 3-4 bottles of poured media, **or** if making P₅ARP/P₁₀ARP. Set to  50 °C .
- Prepare hood by spraying surfaces with ethanol and allowing hood to equalize for  00:15:00 minutes.
- Open the necessary quantity of Petri plate bags and stack in hood (500 mL agar =~ 1.5 sleeves of 90 mm plates, 2 sleeves of 60 mm plates).
- Standardize pH meter, if applicable (aWA and aPDA).
- Thaw any antibiotics, if applicable.

Making Media

30m

1 Select media types from below:

Preparation

STEP CASE

Water Agar (1.5%)

6 steps

WA and aWA (1.5%)

2 Water Agar preparation:

A	B	C
	500 mL	250 mL
Bottle size (mL)	1000	500
Agar (g)	7.5	3.75
RO Water (mL)	500	250

500 mL ~ 1.5 sleeves of 90 mm plates or 2 sleeves of 60 mm plates.



3 Add stir bar and stir contents well.

Place a loose orange cap. Label cap with autoclave tape (1.5% WA) and date.



4 Autoclave liquid setting for 20 minutes to sterilize.

Acidified WA (aWA)

- 5 To make acidified WA (aWA):
- Adjust to  3.9 **after autoclaving** with [M] 10 % (v/v) lactic acid
-  ~0.5 mL for 500 mL of media

After Pouring Media

- 6
- Mark plates with media color codes.
 - Turn off water bath.
 - Turn off autoclave, if no longer needed.
 - Spray inside of hood with ethanol.
 - Close filling port in pH meter, if necessary.
 - Tightly cap empty media bottles so agar doesn't dry out.
- 7 When media has solidified, place in bags labeled with media name, your initials, and date.

Color Codes	
Medium	Code
water agar	I
acidulated w/lactic acid	add I after blue line
P ₅ ARP	I
P ₁₀ ARP	II
PDA	I
½ PDA	II
¼ PDA	III
acidulated w/lactic acid	add I after black lines
Cornmeal	I
V8-juice + CaCO ₃	II
5/2 agar	I
Malt agar	II
Luria-Bertani agar	III
+ streptomycin	+
+ chloramphenicol	+
+ kanamycin	+
+ rifampicin	+
+ mefoxin	+
+ ampicillin	+

Date & initial all media on sleeves!

Bennett Lab media color codes.