



Sep 08, 2020

GNRI's COVID-19 Antigen Lateral Flow Test Device Use Steps

Dr Georgios Gerardos¹¹Great North Research and Innovation*In Development*dx.doi.org/10.17504/protocols.io.bk34kyqw

Dr Georgios Gerardos

DOI

dx.doi.org/10.17504/protocols.io.bk34kyqw

PROTOCOL CITATION

Dr Georgios Gerardos 2020. GNRI's COVID-19 Antigen Lateral Flow Test Device Use Steps. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bk34kyqw>



LICENSE

— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

IMAGE ATTRIBUTION

Copyright © GNRI 2020. All rights reserved.

CREATED

Sep 08, 2020

LAST MODIFIED

Sep 08, 2020

PROTOCOL INTEGER ID

41820

GUIDELINES

These guidelines illustrate the instructions on how to use GNRI's COVID-19 Antigen Lateral Flow Test Device and the accompanying safety warnings.

MATERIALS TEXT

Foil-Packaged COVID-19 Antigen Lateral Flow Test Device

SAFETY WARNINGS

- Please maintain the social distancing guidelines during all times even if the test results are negative.
- Please wash your hands thoroughly for at least 20 seconds before and after using the test.

DISCLAIMER:

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors,

BEFORE STARTING

Make sure you are in a safe place away from traffic and not operating any machinery. The test is undertaken in room temperature and preferably indoors. If the test is to be taken outdoors please ensure that there is no rain, strong wind or direct sunlight as these may affect the test results.

Please follow all steps 10m 30s

- 1 Obtain a clean container (glass or cup) and produce a saliva sample in the container. 10s
- 2 Remove a new COVID-19 Antigen Lateral Flow Test Device from its foil packaging. 10s
- 3 Place the test strip in the container and leave for 10min. 10m
- 4 Remove the test strip from the container and read results. 10s

Results interpretation:

One Line (Control Line) = COVID-19 Negative

Two Lines = COVID-19 Positive

No Lines = Invalid Test. Please repeat test.

One Line (Test Line) = Invalid Test. Please repeat test.