



Aug 13, 2020

# Universal sandwich enzyme linked immunosorbent assay for investigating protein-LG (SpLG) interactions with immunoglobulins using a SpA-HRP conjugate.

Angel A Justiz-Vaillant<sup>1</sup>, Norma McFarlane-Anderson<sup>2</sup><sup>1</sup>University of the West Indies St. Augustine; <sup>2</sup>University of West Indies. Mona Campus**1** Works for me [dx.doi.org/10.17504/protocols.io.bjp4kmqw](https://dx.doi.org/10.17504/protocols.io.bjp4kmqw)[University of the West Indies](#) [angel.vaillant@sta.uwi.edu](mailto:angel.vaillant@sta.uwi.edu)Angel Justiz-Vaillant  
University of the West Indies St. Augustine

DOI

[dx.doi.org/10.17504/protocols.io.bjp4kmqw](https://dx.doi.org/10.17504/protocols.io.bjp4kmqw)

## PROTOCOL CITATION

Angel A Justiz-Vaillant, Norma McFarlane-Anderson 2020. Universal sandwich enzyme linked immunosorbent assay for investigating protein-LG (SpLG) interactions with immunoglobulins using a SpA-HRP conjugate.. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.bjp4kmqw>

## LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Aug 13, 2020

## LAST MODIFIED

Aug 13, 2020

## PROTOCOL INTEGER ID

40412

- 1 This ELISA is used to study the interaction of a protein-LG (SpLG) with different immunoglobulin preparations from mammalian species.
- 2 The 96 well microtitre plate is coated overnight at 4°C with 2 µg/µl per well of a mixture of protein-L and protein-G in carbonate-bicarbonate buffer pH 9.6.
- 3 Then plate is treated with bovine serum albumin solution and washed 4X with PBS-Tween.
- 4 50 µl of animal serum (1 mg/ml) is added and incubated for 1h at room temperature and the microplate is rewashed 4X with PBS-Tween.

- 5 Then 50 µl of peroxidase-labeled SpA conjugate diluted 1:5000 in PBS-non-fat milk is added to each well and incubated for 1h at RT. The plate is washed 4X with PBS-Tween.
- 6 50 µl of 4 mg/ml o-phenylenediamine solution (OPD) is added and the plate is incubated 15 minutes at RT in the dark.
- 7 The reaction is stopped with 50 µl of 3M H<sub>2</sub>SO<sub>4</sub> solution.
- 8 The plate is visually assessed for the development of colour and read in a microplate reader at 492 nm.
- 9 A cut-off point should be calculated as the mean of the optical density of negative controls x 3. The higher the OD value the higher will be the affinity of SpLG to mammalian immunoglobulins.