


Aug 26, 2024

## Pollen germination on wheat stigmas

DOI

[dx.doi.org/10.17504/protocols.io.kxygxy42ol8j/v1](https://dx.doi.org/10.17504/protocols.io.kxygxy42ol8j/v1)



Marina Millán Blánquez<sup>1</sup>

<sup>1</sup>John Innes Centre



Marina Millán Blánquez

John Innes Centre

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.kxygxy42ol8j/v1](https://dx.doi.org/10.17504/protocols.io.kxygxy42ol8j/v1)

**Protocol Citation:** Marina Millán Blánquez 2024. Pollen germination on wheat stigmas . protocols.io

<https://dx.doi.org/10.17504/protocols.io.kxygxy42ol8j/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** August 23, 2024

**Last Modified:** August 26, 2024

**Protocol Integer ID:** 106332

**Keywords:** pollen germination, stigma, wheat, aniline blue

### Abstract

Aniline blue staining of pollen tube growth on the wheat stigma and through the style 4 h and 30 min after pollination.

## Hand pollination

- 1 Please refer to wheat-training.com for detailed explanations on how to perform hand pollinations. Link to pdf: [https://www.wheat-training.com/wp-content/uploads/Wheat\\_growth/pdfs/How-to-cross-wheat-pdf.pdf](https://www.wheat-training.com/wp-content/uploads/Wheat_growth/pdfs/How-to-cross-wheat-pdf.pdf)

## Carpel dissection and fixation

- 2 Using a pair of tweezers, dissect carpels 4.5 h after pollination to allow sufficient time for pollen tube emergence.
- 3 Store samples in a fixative solution of 95% ethanol and absolute acetic acid (75% v/v) and kept at 4 °C until sample preparation for fluorescence microscopy.

## Aniline blue staining of pollinated stigmas

- 4 On the day of sampling, prepare a solution of 0.1% aniline blue in 0.1 M  $K_3PO_4$ . You can prepare a stock solution of 1% aniline blue dissolved in 1x PBS. The stock solution should be kept in the fridge at 4 °C. Use tin foil to avoid exposure to light.
- 5 Wash fixed samples three times for 5 minutes in sterile water and transferred to 0.1% aniline blue solution and kept overnight at 4 °C. Use tin foil to avoid exposure to light.
- 6 Without washing the samples, dissect out the ovary using a sharp razor blade. Try not to damage the stigma. Leave the remaining stigmatic tissue to dry at 45 °C in a hot plate for a few minutes until most of the aniline blue solution is evaporated. Cover the hot plate with an opaque lid to avoid exposure to light.
- 7 Use Vectashield (catalogue No. H-1000-1, 2BSCIENTIFIC LTD) as an antifade mounting media to preserve fluorescence.
- 8 Happy microscopy.