

Sep 13, 2024

RNAi imaging

 In 1 collection

DOI

dx.doi.org/10.17504/protocols.io.q26g7y8wkgwz/v1

Ben Jenkins¹

¹University of Oxford



Ben Jenkins

University of Oxford

OPEN  ACCESS



DOI: **dx.doi.org/10.17504/protocols.io.q26g7y8wkgwz/v1**

Protocol Citation: Ben Jenkins 2024. RNAi imaging. protocols.io **<https://dx.doi.org/10.17504/protocols.io.q26g7y8wkgwz/v1>**

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: October 11, 2022

Last Modified: September 13, 2024

Protocol Integer ID: 71179


Abstract

Imaging of *P. bursaria* cultures




Preparing plates

- 1 In a sterile hood, resuspend plate using a multichannel pipette.

#Ensure pipettes do not overlap with multiple Wells
#Pipette up and down ~20 times (10x in the middle, 10x moving clockwise around the edge)
- 2 Spin plate at  800 x g, 4°C, 00:05:00 5m


800* (*=xg)
Ensure rad is set to 16.3
Ensure that centrifuge is loaded with balance plate

Imaging plates 18m


- 3 Transfer plate to the HC screener, aligned to the left hand corner.
- 4 Run protocol "Paramecium_RNAi_A" (this takes ~  00:10:00) 10m

Name: rnai_[experiment name]_plate[*]A_t[*]_[date]
Select storage: Ben_2

[*] corresponds to the specific plate number, and day, respectively
- 5 Repeat Step 2, and transfer plate back to the HC screener. 10m

Run protocol "Paramecium_RNAi_B" (this takes ~  00:10:00)

Name: rnai_[experiment name]_plate[*]B_t[*]_[date]
Select storage: Ben_2

[*] corresponds to the specific plate number, and day, respectively
Running both protocols A and B will collectively image the entire plate, these are run separately so that cells do not have time to move from the bottom of the well for imaging
- 6 Return plate to  23 °C culture room.