


JUN 01, 2023

Abundance of fungal hyphae in seawater by epifluorescence microscopy using Calcofluor White stain method

M.H. Gutiérrez¹

¹COPAS Sur-Austral, Universidad de Concepción, Chile



Amy Gladfelter

ABSTRACT

Protocol to stain fungal hyphae from seawater with Calcofluor White.

ATTACHMENTS

[Abundance of fungi in seawater.pdf](#)

GUIDELINES

References

Damare S, Raghukumar C (2008) Fungi and macroaggregation in deep-sea sediments. *Microb Ecol* 56:168–177

Jebarag CS, Raghukumar C (2009) Anaerobic denitrification in fungi from the coastal marine sediments off Goa, India. *Mycol Res* 113:100-109

Newell, S.Y., Statzell-Tallman, A. (1982) Factors for conversion of fungal biovolume values to biomass, carbon and nitrogen: variation with mycelial ages, growth conditions, and strains of fungi from a salt marsh. *Oikos* **39**: 261-268

Rasconi S, Jobard M, Jouve L, Sime-Ngando T (2009) Use of calcofluor white for detection, identification, and quantification of phytoplanktonic fungal parasites. *Appl Environ Microbiol* 75:2545-2553

SAFETY WARNINGS



See SDS (Safety Data Sheet) for hazards and safety guidelines.

OPEN ACCESS

DOI:
dx.doi.org/10.17504/protocols.io.36wgq933klk5/v1

Protocol Citation: M.H. Gutiérrez 2023. Abundance of fungal hyphae in seawater by epifluorescence microscopy using Calcofluor White stain method. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.36wgq933klk5/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited


Protocol status: Working
We use this protocol and it's working

Created: Nov 02, 2018



Last Modified: Jun 01, 2023

PROTOCOL integer ID:
17442


Collection



- 1 Collect samples in sterile 50 mL polypropylene tubes and fix with formaldehyde or glutaraldehyde (2% final concentration). Store samples at  4 °C in the dark.

Filtration


- 2 Filter  5 mL to  30 mL (depending on the environment) of seawater by 0.22-µm mesh black 25 mm diameter polycarbonate filters (Millipore Corp.).

Stain

- 3 Stain filters with the retained material from filtration by applying  600 µL of aqueous 0.1% Calcofluor White directly to sample and filter, making sure to cover the entire area of the filter.

Allow  00:05:00 to  00:10:00 then use vacuum to remove the excess stain from the filter. **Avoid any light source.**

Abundance

- 4 Place the filter (sample side up) onto a slide and add 1 drop of non-fluorescent immersion oil to the top of the filter, then cover with a cover slip.
- 5 Count slides immediately on epifluorescent microscope or store at  -20 °C .
- 6 Use an epifluorescence microscope equipped with UV filter (used for DAPI, eg. filter set 49 Carl Zeiss Ltd., 365 nm excitation and 445-450 nm emission band pass) to examine at 1000X the entire effective area of the filters.
- 7 Count all hyphae identified and record their length and width. Use cylinder volume as a

morphological approximation to estimate the biovolume of fungal filaments.

Note

Biomass can be estimated based from biomass: biovolume ratios described from fungi (Newell and Statzell-Tallman, 1982).