



JAN 24, 2024

## Immunohistology

In 1 collection

Mai-Anh Vu<sup>1</sup>, mwhowe<sup>1</sup>

<sup>1</sup>Department of Psychological and Brain Sciences, Boston University, Boston, MA, United States

Team Cragg



Mai-Anh Vu

OPEN ACCESS



### ABSTRACT

We have developed a new micro-fiber array approach capable of chronically measuring and optogenetically manipulating local dynamics across over 100 targeted locations simultaneously in head-fixed and freely moving mice, enabling investigation of cell-type and neurotransmitter-specific signals over arbitrary 3-D volumes . This protocol includes the immunohistology steps we used to confirm viral expression, and/or stain for neurons after our micro-CT scanning and contrast agent soaking procedure.

#### DOI:

[dx.doi.org/10.17504/protocols.io.rm7vzxed5gx1/v1](https://dx.doi.org/10.17504/protocols.io.rm7vzxed5gx1/v1)

**Protocol Citation:** Mai-Anh Vu, mwhowe 2024. Immunohistology. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.rm7vzxed5gx1/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**Created:** Jan 09, 2024

**Last Modified:** Jan 24, 2024

PROTOCOL integer ID: 93202

**Keywords:** ASAPCRN

**Funders Acknowledgement:**

Aligning Science Across  
Parkinson's (ASAP)  
Grant ID: ASAP-020370

## Removing Lugol's contrast agent from tissue

- 1 In preparation for immunohistology, Lugol's solution was removed from brain tissue after CT imaging by soaking the implanted brains in a solution of 5% sodium thiosulfate (STS, Carolina) in 1% PBS for 4-6 days (Hopkins et al., 2015).
- 2 The implant was then removed, the brain returned to the STS solution for one hour.

## Slicing

- 3 The brains were then moved to a solution of 30-40% sucrose in PBS.
- 4 Once the brains sank, coronal sections (50  $\mu$ m) were sliced with a cryostat (Leica CM3050 S) and transferred to PBS for storage.

## Immunostaining

- 5 To stain for dLight1.3b or Chr2 expression:

- 5.1 Sections were then blocked and treated in a 5% normal goat serum (NGS) and 0.2% PBS triton (Sigma Aldrich) solution (PBST).
- 5.2 Treated sections were then incubated for 24-48 hours at 4°C in 5% NGS, 0.2% triton, and a GFP primary antibody (chicken polyclonal antibody; 1:1000, ThermoFisher Scientific, No. A10262).
- 5.3 Slices were then washed, and incubated for 2 hours at room temperature in 5% NGS, 0.2% PBST, and an Alexa 488 secondary antibody (goat anti-chicken; 1:200, ThermoFisher Scientific, No. A-11039).

## 6 For neuronal staining:

- 6.1 Sections were then blocked and treated in a 5% bovine serum albumin (BSA) and 0.2% PBS triton (Sigma Aldrich) solution (PBST).
- 6.2 Treated sections were then incubated for 24-48 hours at 4°C in 5% BSA, 0.2% triton, and a NeuN primary antibody (guinea pig antibody; 1:500, Synaptic Systems, 266 004).
- 6.3 Slices were then washed, and incubated for 2 hours at room temperature in 5% BSA, 0.2% PBST, and an Alexa 647 secondary antibody (goat anti-guinea pig; 1:200, ThermoFisher Scientific, No. A-21450).

## Mounting and imaging

- 7 Finally, stained sections were mounted onto glass slides using Vectashield Mounting Medium with DAPI (Vector Laboratories, H-2000).

8 Confocal images were acquired with a Zeiss LSM 800.