



Nov 02, 2022

cDNA protocol - Thermo Scientific RevertAid RT

Lauren Ponisio¹, Jocelyn Zorn¹¹University of Oregon

1 Works for me

 Sharedx.doi.org/10.17504/protocols.io.36wgq4zwxvk5/v1

FFAR Ponisiolab

 Jocelynz

ABSTRACT

cDNA protocol - Thermo Scientific RevertAid RT

DOI

dx.doi.org/10.17504/protocols.io.36wgq4zwxvk5/v1

PROTOCOL CITATION

Lauren Ponisio, Jocelyn Zorn 2022. cDNA protocol - Thermo Scientific RevertAid RT . **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.36wgq4zwxvk5/v1>

FUNDERS ACKNOWLEDGEMENT

Foundation of Food and Agricultural Research
Grant ID: CA18-SS-0000000009

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Sep 02, 2021

LAST MODIFIED

Nov 02, 2022

Prepare workspace

- 1 Thaw RNA and reagents on ice. Label four foil lids: one to replace RNA lid, one for reaction and two for aliquots
- 2 UV sterilize 3 RNase free PCR plates, 8 boxes 10uL filter tips, and DEPC treated H₂O
- 3 Once thawed, flash spin down RNA in plate centrifuge and spin down reagents in mini centrifuge

cDNA conversion

- 4 Add **1 µL** RNA (**1 ng** - **5 µg**) to empty PCR plate

- 5 Add **1 µL** [Random Hexamer Primer Thermo](#)
Fisher Catalog #S0142

[Random Hexamer Primer Thermo](#)

If not using **Fisher Catalog #S0142**
will be different. See product handbook

protocol

[RevertAid RT Reverse Transcription Kit Thermo](#)

Fisher Catalog #K1691

- 6 Add **10 µL** DEPC treated H₂O

If used amount other than **1 µL** RNA (step 1), adjust H₂O volume to bring total reaction volume to **12 µL**

7 Add  4 μL 5X reaction buffer

8  RiboLock RNase Inhibitor Thermo


Add  1 μL Scientific Catalog #E00381
U/ μL)

(20

9  dNTP Mix (10 mM each) Thermo Fisher

Add  2 μL Scientific Catalog #R0191

10 Add  1 μL

 RevertAid Reverse Transcriptase (200 U/ μL) Thermo
Fisher Catalog #EP0441

Total reaction volume should be  20 μL

11 Mix gently with pipette and centrifuge briefly on "short"

12 Place in thermocycler and run program "cDNA"

Thermocycler cDNA program: Lid 105°C. 25°C 5 min, 42°C 60 min, 70°C for 5 min, ∞ hold 4°C.

13 Aliquot product into microplates and seal with labeled foil lid

14 Store at -20°C for up to 1 week or -80°C for longer term storage