

Jun 01, 2021

Expression and purification protocol of the Human (Homo sapiens) mCherry-LC3B Ubiquitin-like modifier

Dorotea Fracchiolla¹¹Team Hurley

3

Works for me



Share

dx.doi.org/10.17504/protocols.io.btiunkew

Dorotea Fracchiolla

Sascha Martens lab, University of Vienna, Max Perutz Labs - ...

ABSTRACT

This protocol describes expression and purification procedures for obtaining mCherry-tagged human recombinant Ubiquitin-like modifier LC3B (MAP1LC3B, Microtubule-associated proteins 1A/1B light chain 3B) lacking five C-terminal

amino acids to allow in vitro protein conjugation to target PE, Phosphatidyl Ethanolamine.

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

<https://doi.org/10.1083/jcb.201912098>

DOI

dx.doi.org/10.17504/protocols.io.btiunkew

PROTOCOL CITATION

Dorotea Fracchiolla 2021. Expression and purification protocol of the Human (Homo sapiens) mCherry-LC3B Ubiquitin-like modifier. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.btiunkew>



MANUSCRIPT CITATION please remember to cite the following publication along with this protocol

<https://doi.org/10.1083/jcb.201912098>

KEYWORDS

expression, purification, recombinant, Microtubule-associated proteins 1A/1B light chain 3B, LC3B, mCherry-tag

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

IMAGE ATTRIBUTION

© Dorotea Fracchiolla

CREATED

Mar 22, 2021

LAST MODIFIED

Jun 01, 2021

OWNERSHIP HISTORY

Mar 22, 2021



zahara93.zs

May 25, 2021




Dorotea Fracchiolla


Sascha Martens lab, University of Vienna, Max Perutz Labs - Vienna (AT)

GUIDELINES

Affinity purification




Chromatogram of His-tag affinity purification for mCherry-hLC3B.




Coomassie BB stained gel of His-tag affinity purification for mCherry-hLC3B.

Size Exclusion Chromatography



Chromatogram of Size Exclusion purification for mCherry-hLC3B.



Coomassie BB stained gel of Size Exclusion Chromatography purification for mCherry-hLC3B.

MATERIALS TEXT

Materials and Reagents

- Escherichia coli Rosetta pLyss cells
 - Luria Bertani (LB) medium with antibiotics (final conc. 50µg/ml Ampicillin, 34µg/ml Chloramphenicol)
 - IPTG (isopropyl-β-D-thiogalactopyranoside)
 - 37°C shaker incubator
 - sterile flasks/sterile pipettes
 - tip sonicator
-
- **Lysis Buffer:** 50mM Hepes pH=7.5, 300mM NaCl, 10mM Imidazole, 2mM MgCl₂, 2mM β-mercaptoethanol, 1mM Pefablock, Complete Protease Inhibitors (EDTA-free CIP tablet, Roche), DNase (Sigma).
 - **Buffer A:** 50mM Hepes pH=7.5, 300mM NaCl, 10mM Imidazole (filtered and degassed) + 1mM β-mercaptoethanol;
 - **Buffer B:** 50mM Hepes pH=7.5, 300mM NaCl, 300mM Imidazole (filtered and degassed) + 1mM β-mercaptoethanol;
 - **Size Exclusion Chromatography (SEC) Buffer:** 25mM Hepes pH=7.5, 150mM NaCl, 1mM DTT (Dithiothreitol).
- Note: all purification buffers are filtered and degassed. Reducing agents (β-mercaptoethanol and Dithiothreitol) are added after degassing step.

Columns:- HT 5ml column (GE Healthcare)
- S75_16/60 (GE Healthcare)

Gels:10% SDS-PAGE

SAFETY WARNINGS

Please refer to the Safety Data Sheets (SDS) for health and environmental hazards.

BEFORE STARTING



General information

Insert: Homo sapiens LC3B, [NP_073729.1](#); Expression system: E.Coli Rosetta pLyss; plasmid origin: Sascha Martens Lab, Addgene 169168, lab internal construct database number SMC948; backbone: pET-Duet1; plasmid resistance: Ampicillin; tags & cleavage sites: N-term 6xHis, followed by Tobacco Etch Virus (TEV) cleavage site, mCherry tag, LC3B ORF. Ext coeff: 41830 M⁻¹ cm⁻¹, MW 47,89 kDa.




Protein Expression

16h 15m




1 

Transform plasmid DNA (Addgene 169168, SMC948) into E.Coli Rosetta pLyss cells and plate on Ampicillin/Chloramphenicol LB agar plate for  **Overnight** at  **37 °C**.

2 

The following day, inoculate a  **100 mL LB + Amp/Cam pre-culture** with 1-2 colonies and grow  **Overnight** at  **37 °C** shaking.

3

The following day, use  **20 mL pre-culture** to inoculate  **1 L LB medium + Amp/Cam** at  **37 °C** until an OD₆₀₀ (Optical Density) of 0.4 is reached.

4

Cool down the culture to  **18 °C** and grow to an OD₆₀₀ = 0.8.

5 Induce protein expression with **100 Micromolar (μM) IPTG** and grow for a further **16:00:00** at **18 °C** shaking. ^{16h}

6  15m

Pellet cells at **4000 rpm, 4 °C, 00:15:00** in a Sorvall RC6+ centrifuge (Thermo Scientific), discard supernatant and resuspend pellets in ice cold lysis buffer (25 ml/1 lt culture).

7 

Flash freeze resuspended pellets in liquid nitrogen and store at **-80 °C** until purification.

Protein Purification

8 Perform His-Trap affinity purification followed by Size Exclusion Chromatography.

Purification of LC3B Ubiquitin-like modifier 30m 30s

9 Cells are lysed via freeze/thaw cycles and sonication: thaw pellet corresponding to **1 L** culture by freeze/thawing in **Room temperature** water bath. All following steps are to be executed at **4 °C** or **On ice**.

10  30s

Lyse cells by sonicating them with an immersion Tip Sonicator (**2x 00:00:30**). Note: adjust times and intensity according to the available instrument.


11 Clear lysate by spinning in a Beckman centrifuge, at **40000 x g, 4 °C, 00:30:00**, **Ti45 Rotor**. 30m

12 Filter supernatant through a 0.45 μm filter and inject onto a 5ml HT column operating at **4 °C** pre-equilibrated in Buffer A at 1ml/min flow rate.

13 

Wash column for **5 column volumes (CV) with Buffer A** at 2 ml/min flow rate to remove unspecific bound proteins.

14 Elute protein through a step elution gradient in 50mM, 75mM, 100mM, 150mM, 200mM and 300mM Imidazole concentration. Perform elution at 1ml/min flow rate. Collect peak fractions at each step.




Chromatogram of His-tag affinity purification for mCherry-LC3B.

- 15 Check peak fractions of each step on a SDS-PAGE (see gel below). Pool and concentrate those containing the protein of interest (75 and 100mM Imidazole step) by spinning at **4 °C** down in an Amicon Filter 30kDa cut-off to


2 mL final volume .

Centrifugation steps are kept short (**00:05:00**) to avoid protein local concentration/aggregation on the filter.



Coomassie BB stained gel of His-tag affinity purification for mCherry-hLC3B.


- 16 Inject **2 mL protein** onto a S75_16/60 column pre-equilibrated in SEC buffer.



Chromatogram of Size Exclusion purification for mCherry-hLC3B.

- 17 Check fractions on a 10% SDS-PAGE (see gel below). Pool and concentrate those containing the proteins of interest down at **4 °C** in an Amicon Filter 30kDa cut-off.

Centrifugation steps are kept short (**00:05:00**) to avoid protein local concentration/aggregation on the filter.



Coomassie BB stained gel of Size Exclusion Chromatography purification for mCherry-hLC3B.

18 

Measure protein absorbance at A_{280} with a Spectrophotometer against Size Exclusion Chromatography purification buffer (MW = 47,89 kDa ; Extinction coefficient = $41830 \text{ M}^{-1} \text{ cm}^{-1}$).

19 Aliquot protein, flash freeze in liquid Nitrogen and store at -80°C .