




Aug 22, 2022

DAB immunostaining of thin, fixed mouse brain tissue sections using HNA or NCAM to characterize human iPSC-derived cell xenografts

Benjamin Trist¹, Ashish Mathai¹, Asheeta Prasad¹¹The University of Sydney

1 Works for me

 Sharedx.doi.org/10.17504/protocols.io.14egn7pb6v5d/v1 Benjamin Trist

ABSTRACT

This protocol describes our use of chromogenic 3,3'-diaminobenzidine (DAB) immunohistochemistry to identify human iPSC-derived cells within thin, fixed mouse brain tissue section series'. We apply this workflow for post-mortem assessment of the survival and growth of human iPSC-derived cells which have been transplanted into the living brain of athymic mice.

ATTACHMENTS

[it75bj7ap.docx](#)

DOI

dx.doi.org/10.17504/protocols.io.14egn7pb6v5d/v1

PROTOCOL CITATION

Benjamin Trist, Ashish Mathai, Asheeta Prasad 2022. DAB immunostaining of thin, fixed mouse brain tissue sections using HNA or NCAM to characterize human iPSC-derived cell xenografts. **protocols.io**
<https://protocols.io/view/dab-immunostaining-of-thin-fixed-mouse-brain-tissue-ce7fthjn>



FUNDERS ACKNOWLEDGEMENT

Michael J Fox Foundation
Grant ID: ASAP-000497

KEYWORDS

NCAM, HNA, Human-to-mouse xenograft, Human iPSC, Immunohistochemistry

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Aug 12, 2022

LAST MODIFIED

Aug 22, 2022

OWNERSHIP HISTORY

Aug 12, 2022  madhavi.d

Aug 19, 2022  Benjamin Trist

PROTOCOL INTEGER ID




68551

MATERIALS TEXT

Equipment:

- Horizontal rocker
- Vortex
- Microcentrifuge
- Glass petri dish
- Oven

Consumables:

-  **20 mL** scintillation vials
- Paint brushes
- Gelatin-Chrom Alum-coating microscope slides
 1. See related protocol - Coating superfrost microscope slides with gelatin-chromium potassium sulfate
- Microscope slide coverslips (no. 1.5,  **25 mm** x  **75 mm**)
- Glass pipettes
- Rubber teats
- Transfer pipettes

Key reagents:

- Optimal Cutting Temperature (OCT) compound
- Bovine Serum Albumin (BSA)

- Casein
- Sodium citrate
- Tween-20 and Triton X-100
- Ethanol
- Hydrogen peroxide (H₂O₂)
- DEPEX
- 3, 3'-diaminobenzidine (Sigma #D5905)
- Antibodies

 **NCAM**

1 **Abcam Catalog #ab75813**

2. HNA (Novus #NOVNB-313912)

3. Biotinylated anti-rabbit secondary antibody (Vector Labs #BA-1000)

 **Avidin/Biotin HRP Complex Vector**

4. **Laboratories Catalog #PK-6100**

Solutions:

- 1x PBS, pH 7.4

A	B
Antigen retrieval (AR) buffer	
Sodium citrate	2.94 g (10 mM)
Tween-20	500 µL (0.05%)
Up to 1L with dH ₂ O, pH 6.0	

- Quenching solution
10mL (3.3%) 33% H₂O₂, 50mL (50%) ethanol up to 100mL with 1x PBS
- 1x PBST
500 µL (0.05%) Tween-20 in 1L 1x PBS

A	B
Blocking solution	
Casein	1 g (1% w/v)
Triton X-100	250 µL (0.25% v/v)
Glycine	1.5 g (1.5% w/v)
BSA	5 g (5% w/v)
Up to 100mL with 1xPBS	

Material input (animal, cell, tissue, fraction details):

Thin, fixed athymic mouse brain tissue sections prepared from whole mouse brains grafted with human iPSC-derived neural progenitor cells.

Day 1 (~4-6 hrs)

- 1 Pre-heat oven and AR buffer to **70 °C**.
- 2 Label scintillation vials to match labels on section storage plates (mouse and section series IDs, name, date etc.).

3



Transfer sections into scintillation vials using a transfer pipette or fine paintbrush.

4



Remove anti-freeze solution and perform 3x 7 min washes in 1x PBS at

Room temperature with gentle agitation.

- Slow shaking on an orbital rocker recommended for washes/incubations to ensure even contact with solutions.
- Use a glass pipette and rubber teat to remove solution during wash changes.
- Anti-freeze solution must be rinsed off prior to immunostaining.

4.1 Remove anti-freeze solution and wash in 1x PBS at **Room temperature**^{7m} for **00:07:00** (1/3).

4.2 Remove anti-freeze solution and wash in 1x PBS at **Room temperature**^{7m} for **00:07:00** (2/3).

4.3 Remove anti-freeze solution and wash in 1x PBS at **Room temperature**^{7m} for **00:07:00** (3/3).

Antigen retrieval (AR)

5



30m

Incubate sections in AR buffer for **00:30:00** at **70 °C**.

6 Preheat AR buffer to **70 °C** prior to addition to sections.

7 After antigen retrieval, allow sections to cool for **00:30:00** before proceeding with staining. ^{30m}


8 

Perform 3x 7 min washes in 1x PBST with agitation.

8.1 

7m

Wash in 1x PBST with agitation for **00:07:00** (1/3).

8.2 

7m

Wash in 1x PBST with agitation for **00:07:00** (2/3).

8.3 

7m

Wash in 1x PBST with agitation for **00:07:00** (3/3).

Quenching step

9 

30m

Incubate sections in quenching solution for **00:30:00** at **Room temperature** with gentle agitation.

10 

Perform 3x 7 min washes in 1x PBST with gentle agitation.

10.1 


7m

Wash in 1x PBST with gentle agitation for **00:07:00** (1/3).




7m

10.2 

Wash in 1x PBST with gentle agitation for  **00:07:00** (2/3).

7m



10.3 

Wash in 1x PBST with gentle agitation for  **00:07:00** (3/3).

Blocking step



11 

1h

Incubate sections in blocking solution for  **01:00:00** at  **Room temperature** with gentle agitation.

Primary antibody step

12  

Incubate sections with NCAM (1:20,000) or HNA (1:15,000) primary antibodies diluted in blocking buffer  **Overnight** at  **4 °C** with gentle agitation.


Day 2 (~8hrs)

13 

Perform 3x 7 min washes in 1x PBST with gentle agitation.


13.1 

7m

Wash in 1x PBST with gentle agitation for  **00:07:00** (1/3).


13.2 

7m

Wash in 1x PBST with gentle agitation for  **00:07:00** (2/3).

13.3 

7m

Wash in 1x PBST with gentle agitation for  **00:07:00** (3/3).

Secondary antibody step



2h

14 

Incubate sections in anti-rabbit biotinylated secondary antibody 1:500 diluted in blocking buffer that has been diluted 2-fold for 🕒 **02:00:00** at 🌡 **Room temperature** .

15 

Perform 3x 7 min washes in 1x PBST with gentle agitation.

15.1 

7m

Wash in 1x PBST with gentle agitation for 🕒 **00:07:00** (1/3).

15.2 

7m

Wash in 1x PBST with gentle agitation for 🕒 **00:07:00** (2/3).

15.3 

7m

Wash in 1x PBST with gentle agitation for 🕒 **00:07:00** (3/3).

Tertiary complex step

16 

2h

Incubate sections in Avidin-Biotin Complex (ABC) kit solution (Vector Laboratories) for 🕒 **02:00:00** at 🌡 **Room temperature** .

17 Prepare tertiary complex 🕒 **00:30:00** prior to use according to the manufacturer's instructions. 30m

- 🧴 **100 µL** A + 🧴 **100 µL** B + 🧴 **9800 µL** 1x PBS (1:100 A + 1:100 B in 1x PBS).

18 

Perform 3x 7 min washes in 1x PBST with gentle agitation.


18.1 

7m

Wash in 1x PBST with gentle agitation for 🕒 **00:07:00** (1/3).


18.2 

7m

Wash in 1x PBST with gentle agitation for  **00:07:00** (2/3).




18.3 

7m

Wash in 1x PBST with gentle agitation for  **00:07:00** (3/3).

Chromogen step



19 Perform DAB staining as follows;

20 Prepare DAB solution by dissolving  **10 mg** DAB tablet into  **20 mL** PBS ( **0.5 mg/mL**).

21 Filter through Whatman paper #1 or 0.22 um syringe filter before use.


NB: DAB is a suspected mutagen and should be handled with care.

22 

Prepare DAB-H2O2 solution **immediately** prior to use by adding  **10 µL** 30% H2O2 per  **5 mL** DAB solution and mix thoroughly.

23 

8m

Incubate sections with DAB-H2O2 solution for  **00:08:00** .

24 

Perform 3x 7 min washes in 1x PBST with gentle agitation.

24.1 

7m

Wash in 1x PBST with gentle agitation for 🕒00:07:00 (1/3).

24.2



7m

Wash in 1x PBST with gentle agitation for 🕒00:07:00 (2/3).

24.3



7m

Wash in 1x PBST with gentle agitation for 🕒00:07:00 (3/3).

25



Mount tissue sections onto super-frost slides pre-coated with gelatin-chrome alum and allow to dry at 📍 **Room temperature** 🕒 **Overnight** .

Day 3 (2 days later)

33m

26

Process slide-mounted tissue sections through the following solutions;

27

dH₂O 🕒00:03:00 .

3m

28

50% ethanol 🕒00:03:00 .

3m

29

70% ethanol 🕒00:03:00 .

3m

30


95% ethanol 🕒00:03:00 .

3m

31

100% ethanol 🕒00:03:00 .

3m

- 32 100% ethanol ⌚ 00:03:00 . 3m
- 33 Xylene ⌚ 00:05:00 . 5m
- 34 Xylene ⌚ 00:05:00 . 5m
- 35 Xylene ⌚ 00:05:00 . 5m
- 36 
- Coverslip slides with DEPEX mounting media** and allow to dry in the fume hood
⌚ **Overnight** before proceeding with microscopy.
- 37 Image sections using bright field microscopy for subsequent xenograft characterization.