



VERSION 2

MAY 23, 2023

OPEN ACCESS

DOI:

dx.doi.org/10.17504/protocols.io.bp2l695zklqe/v2

External link:

<https://rdcu.be/c9RyP>

Protocol Citation: Miguel Chuapoco, Sripriya Ravindra Kumar, Timothy F. Shay, Xinhong Chen, David Brown, Tatyana Dobрева, Qin Huang, Xiaozhe Ding, Yicheng Luo, Pétur H. Einarsson, Alon Greenbaum, Min J. Jang, Benjamin E. Deverman, Viviana Gradinaru 2023. Library tissue handling, viral DNA extraction, and NGS sample preparation.

protocols.io

<https://dx.doi.org/10.17504/protocols.io.bp2l695zklqe/v2> Version created by Miguel Chuapoco

MANUSCRIPT CITATION:


Ravindra Kumar, S., Miles, T.F., Chen, X. *et al.* Multiplexed Cre-dependent selection yields systemic AAVs for targeting distinct brain cell types. *Nat Methods* **17**, 541–550 (2020). <https://doi.org/10.1038/s41592-020-0799-7>

Library tissue handling, viral DNA extraction, and NGS sample preparation V.2

Forked from a private protocol

Sripriya Ravindra
Miguel Chuapoco¹, Kumar¹,
David
Xinhong Chen¹, Brown¹,
Xiaozhe
Ding¹,
Yicheng Luo¹, Einarsson¹,
Min J.
Alon Greenbaum¹, Jang¹,
Viviana Gradinaru¹
Tatyana Dobрева¹, Huang¹,
Pétur H.
Benjamin E.
Deverman¹,
Timothy F. Shay¹,
Qin

¹California Institute of Technology



Miguel Chuapoco

California Institute of Technology

ABSTRACT

This protocol describes the procedure to isolate viral DNA from AAV-transfected tissue and prepare it for next-generation sequencing.

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: May 23, 2023

Last Modified: May 23, 2023




PROTOCOL integer ID:
82290

Marmoset tissue extraction (library selections)

- 1 Raise and house marmosets in compliance with your local Institutional Animal Care and Use Committee (IACUC) and ensure that all protocols and procedures have been approved by the appropriate ethical and regulatory committees.

For adult marmosets, ensure that there is no detectable neutralizing antibodies for AAV9. This can be conducted by The Penn Vector Core at the University of Pennsylvania (<https://gtp.med.upenn.edu/intranethome/core-facilities-internal/immunology-core>)
- 2 Inject marmosets with desired dose of library (e.g. 2×10^{12} vector genomes of library) intravenously (e.g. via the femoral vein).
- 3 At four weeks post-injection, euthanize marmosets and perfuse with 1X phosphate buffered saline.
- 4 Flash freeze tissue (e.g. using 2-methylbutane chilled with dry ice). Separate the brain into coronal blocks and flash freeze the blocks. Store tissue at -80 °C until ready for processing.

DNA Extraction

- 5 Add  100 mg tissue sample (brain, liver, or spinal cord) and  1 mL TRIzol reagent
-  TRIzol Reagent Thermo Fisher Scientific Catalog #15596026 to bead homogenizer tubes. Use prefilled tubes with 1.5 mm Zirconium beads or 2.8 mm stainless steel beads.

Equipment

| | |
|---|-------|
| Prefilled 2.0ml tubes, Zirconium Beads, 1.5mm Triple-Pure - High Impact, 50pk | NAME |
| Homogenizer tubes (1.5 mm Zirconium beads) | TYPE |
| Benchmark Scientific | BRAND |
| D1032-15 | SKU |
| https://www.benchmarkscientific.com/product/d1032-group/ | LINK |

Equipment

| | |
|---|-------|
| Prefilled 2.0ml tubes, Stainless Steel, 2.8mm Acid Washed, 50pk | NAME |
| Homogenizer tubes (2.8 stainless steel) | TYPE |
| Benchmark Scientific | BRAND |
| D1033-28 | SKU |
| https://www.benchmarkscientific.com/product/d1032-group/ | LINK |

- 6 Homogenize tissue in using the following settings:

- Speed: 5.0 m/s
- Time: 30 seconds
- Pause: 1 minute
- Cycles: 2

5m

Incubate for  00:05:00 .

Equipment



| | |
|---|-------|
| BEADBUG 6, SIX POSITION HOMOGENIZER, 115V | NAME |
| Tissue homogenizer (6 position) | TYPE |
| Benchmark Scientific | BRAND |
| D1036 | SKU |
| https://www.benchmarkscientific.com/product/d1036/ | LINK |

Equipment

| | |
|---|-------|
| BEADBLASTER 24 MICROTUBE HOMOGENIZER, 115V | NAME |
| Tissue Homogenizer (24 position) | TYPE |
| Benchmark Scientific | BRAND |
| D2400 | SKU |
| https://www.benchmarkscientific.com/product/d2400/ | LINK |

Note

Samples can be stored at -20 °C for up to year in TRIzol.

- 7 Centrifuge the homogenizer tubes containing the TRIzol solution and homogenized tissue using the following parameters:  12000 x g, 4°C,  00:05:00 . Transfer the supernatant to a new tube (microcentrifuge tube or similar).

5m

Equipment

Centrifuge 5425/5425 R - Microcentrifuge

Refrigerated centrifuge

Eppendorf

2231000909

[https://www.eppendorf.com/us-en/eShop-](https://www.eppendorf.com/us-en/eShop-Products/Centrifugation/Microcentrifuges/Centrifuge-5425-5425R-p-PF-934144)

[Products/Centrifugation/Microcentrifuges/Centrifuge-5425-5425R-p-PF-934144](https://www.eppendorf.com/us-en/eShop-Products/Centrifugation/Microcentrifuges/Centrifuge-5425-5425R-p-PF-934144)

NAME

TYPE

BRAND

SKU

LINK

Equipment

DNA LoBind® Tubes

Microcentrifuge tubes

Eppendorf

022431021

<https://www.eppendorf.com/us-en/eShop-Products/Laboratory-Consumables/Tubes/DNA-LoBind-Tubes-p-PF-56252>






NAME

TYPE

BRAND

SKU

LINK

- 8 Add  600 µL chloroform to each tube for every  1 mL TRIzol used for lysis, vortex briefly, and incubate for  00:03:00 . 3m
- 9 Centrifuge  12000 x g, 4°C, 00:15:00 to separate the nucleic acid phase from the protein phase. Transfer the top aqueous phase to a new tube (approximately  500 µL aqueous phase) 15m

Equipment

Centrifuge 5425/5425 R - Microcentrifuge

NAME

Refrigerated centrifuge

TYPE

Eppendorf

BRAND

2231000909

SKU

[https://www.eppendorf.com/us-en/eShop-](https://www.eppendorf.com/us-en/eShop-Products/Centrifugation/Microcentrifuges/Centrifuge-5425-5425R-p-PF-934144)

LINK

[Products/Centrifugation/Microcentrifuges/Centrifuge-5425-5425R-p-PF-934144](https://www.eppendorf.com/us-en/eShop-Products/Centrifugation/Microcentrifuges/Centrifuge-5425-5425R-p-PF-934144)

Equipment

DNA LoBind® Tubes

NAME

Microcentrifuge tubes

TYPE

Eppendorf





BRAND

022431021

SKU

<https://www.eppendorf.com/us-en/eShop-Products/Laboratory-Consumables/Tubes/DNA-LoBind-Tubes-p-PF-56252>

LINK

- 10** Add 1 equivalent volume of isopropanol, 1/10 volume of sodium acetate, and co-precipitant (e.g.  500 µL isopropanol,  50 µL sodium acetate,  2-3 µL co-precipitant) and vortex briefly. Incubate for  00:10:00.















10m



Sodium Acetate 3M, pH 5.2 Thermo Scientific Catalog
#R1181










GlycoBlue™ Coprecipitant Thermo Scientific Catalog
#AM9516

- 11** Centrifuge  12000 x g, 4°C, 00:10:00 to pellet nucleic acids. Discard supernatant and wash pellet with  1 mL 75% ethanol. Centrifuge again  7500 x g, 4°C, 00:05:00 and discard supernatant. 15m
- 12** Air dry pellet and resuspend in  84 µL PCR clean water
 UltraPure Distilled Water Invitrogen - Thermo Fisher Catalog #10977-015
- 13** Remove RNA by digestion with  3 µL RNase cocktail and digest with  3 µL SmaI.
 Supplement reaction with  10 µL CutSmart. Incubate at  Room temperature for 2-3 hours
 and  37 °C overnight.
-  RNase Cocktail™ Enzyme Mix Thermo Fisher Catalog #AM2286
-  SmaI - 2,000 units New England Biolabs Catalog #R0141S
-  CutSmart Buffer - 5.0 ml New England Biolabs Catalog #B7204S
- 14** Purify with  Zymo DNA Clean & Concentrator - 5 Zymo Research Catalog #D4014

Genome Recovery

30m 40s

- 15** Amplify the region around the diversified genome insertion by PCR with 25 cycles of  98 °C 30m 40s
 for  00:00:10,  60 °C for  00:00:30 and  72 °C for  00:30:00 using
 Q5 High-Fidelity DNA Polymerase - 500 units New England Biolabs Catalog #M0491L and
 50% of the total extracted viral DNA as a template using primers XF
 (ACTCATCGACCAATACTTGTACTATCTCTCTAGAAC) and 588-R2lib-R
 (GTATTCCTTGTTTTGAACCAACCG).

Index addition







1m 40s

- 16** Dilute PCR product 1:100 and use as template for an additional round of PCR amplification around the variable region with primers containing Read1 and Read2 sequences by 10 cycles of 98 °C for 00:00:10 , 59 °C for 00:00:30 , and 72 °C for 00:00:10 using Q5 High-Fidelity DNA Polymerase - 500 units **New England Biolabs Catalog #M0491L** and primers 588i-lib-PCR1-6bpUID-F (CACTCATCGACCAATACTTGACTATCTCTCT) and 588i-lib-PCR1-R (GTATTCCTTGTTTGAACCCAACCG).
- 17** Purify with Zymo DNA Clean & Concentrator - 5 **Zymo Research Catalog #D4014**
- 18** Append Illumina flow cell adapters and unique indices by PCR amplification with NEBNext Multiplex Oligos for Illumina (Dual Index Primers Set 1) - 96 rxns **New England Biolabs Catalog #E7600S** by 10 cycles of 98 °C for 00:00:10 , 59 °C for 00:00:30 , and 72 °C for 00:00:10 using Q5 High-Fidelity DNA Polymerase - 500 units **New England Biolabs Catalog #M0491L**

Clean up and validation

- 19** Run PCR products on a freshly-prepared 2% UltraPure™ Low Melting Point Agarose **Thermo Fisher Scientific Catalog #16520050** gel.
- 20** Verify the expected size of the band and extract from the gel.
- 21** If desired, verify the nucleotide diversity at the randomized insertion site by Sanger sequencing.

Note

If additional material is needed for Sanger sequencing, perform an additional PCR amplification using 15-20 cycles of  98 °C for  00:00:10,  60 °C for  00:00:30 and  72 °C for  00:00:10 with primers NGS-QC-F (AATGATACGGCGACCACCGAG) and NGS-QC-R (CAAGCAGAAGACGGCATACGA).