



Mar 19, 2022

Nuclei isolation from snap-frozen human placental tissue for bulk ATACseq

Scott Lindsay-Hewett¹¹University of California, San Diego

1

dx.doi.org/10.17504/protocols.io.kqdg3p2kel25/v1

Human BioMolecular Atlas Program (HuBMAP) Method Development Community
Tech. support email: Jeff.spraggins@vanderbilt.edu



Scott Lindsay-Hewett

This protocol describes the isolation of nuclei from snap-frozen human placental tissue for bulk ATACseq. It is a modified version of a protocol (SCBL Protocols - 10x Multiome Nuclei Isolation) supplied by the Jackson Laboratory (JAX), which is itself modified from Sigma Aldrich's Nuclei Isolation Kit: Nuclei EZ Prep protocol, and 10x Genomics' Demonstrated Protocol – Nuclei Isolation for Single Cell Multiome ATAC + Gene Expression Sequencing. It is recommended to work with no more than four samples at a time so that all samples can be processed in a timely manner.

DOI

dx.doi.org/10.17504/protocols.io.kqdg3p2kel25/v1

Scott Lindsay-Hewett 2022. Nuclei isolation from snap-frozen human placental tissue for bulk ATACseq. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.kqdg3p2kel25/v1>



protocol ,

Mar 17, 2022

Mar 19, 2022

59544

Equipment:



Benchtop microfuge with cooling function
Fume hood
Laminar flow hood
Cordless motor for disposable pestles (Z359971-1EA)
Metal plate
Automated cell counter
Heating block
Set of pipettes
Microscope


Supplies/reagents:

60mm Petri dishes (Fisher, cat # FB0875713A)
LoBind Eppendorf tubes (Eppendorf, cat # 022431021)
#22 sterile disposable scalpels (Fisher, cat # 22-079-714)
Trizma hydrochloride solution, 1M, pH 7.4 (Sigma, cat # T2194-100ML)
Magnesium chloride solution, 1M (Sigma, cat # M1028-100ML)
Nuclei Isolation Kit: Nuclei EZ Prep (Sigma, cat # NUC101-1KT)
Magnesium acetate tetrahydrate (Sigma, cat # M2545-250G) - use to prepare 300mM solution, filter sterilize
N,N-Dimethylformamide (Sigma, cat # 227056-100ML)
Trypan Blue solution, 0.4% (Sigma, cat # T8154-20ML)
Nuclease-free water
Bel-Art Disposable Pestles, polypropylene, With 1.5 mL Tubes, pack of 100 (Sigma, cat # BAF199230000-100EA)
Single-use tweezers, angled, blue, sterile (Thomas Scientific, cat # 19A00N163)
Dithiothreitol (DTT), 0.1M Solution (Fisher, cat # 707265ML)
Pluriselect pluriStrainer Mini 40 UM, 100/Bag (Fisher, cat # NC1469671)
Pluriselect PluriStrainer Mini 20 µm Cell Strainer, 100/Pk, Sterile (Fisher, cat # NC1423042)
Digitonin, 5% (Fisher, cat # BN2006)
Tween-20, 10%
NaCl (5 M), RNase-free (Fisher, cat # AM9760G)
Bovine Serum Albumin solution, 10% in DPBS (Sigma, cat # A1595-50ML)
200ul wide bore filter tips (Fisher, cat # 2069GPK)
1000ul wide bore filter tips (Fisher, cat # 2079GPK)
Cell Counting Slides for TC10™/TC20™ Cell Counter, Dual-Chamber (Bio-Rad, cat # 1450011)
Tris base (Fisher, cat # BP152-500) - use to prepare 1M solution; adjust to pH 7.8 using glacial acetic acid, filter sterilize
Potassium acetate (Sigma, cat # P5708-500G) - use to prepare 3M solution, filter sterilize
Bucket of ice
Filter tips

Preparation

1 Pre-chill microfuge to **4 °C**.

- 2 Prepare bucket of ice, and chill metal plate (allows for a flat surface to place Petri dish on top of - to keep tissue cool while cutting).
- 3 Dispense  **200 μ L Nuclei EZ Storage Buffer** into microfuge tube and place on ice
- 4 Thaw  **110 μ L** aliquot of **[M]100 millimolar (mM) DTT** and place on ice.

Prior to starting this protocol, prepare and freeze  **110 μ L** aliquots of **[M]100 millimolar (mM) DTT** from stock solution.

- 5 Thaw  **500 μ L** aliquot of tagmentation buffer and place on ice.

Prior to starting this protocol, prepare  **10 mL** tagmentation buffer as follows:

 **330 μ L 1M Tris-acetate, pH 7.8**

 **220 μ L 3M Potassium acetate**

 **367 μ L 300mM Magnesium acetate**

 **1.6 mL Dimethylformamide**

 **7.483 mL Nuclease-free water**

Freeze  **500 μ L** aliquots of tagmentation buffer.



Take care when handling dimethylformamide. Wear appropriate personal protective equipment and work in a fume hood. Dispose of waste appropriately.

- 6 Pull snap-frozen placenta samples (~150-200 mg each) from -80 °C and place on dry ice until ready to process.

This protocol describes nuclei isolation from 4 snap-frozen samples at a time. It is helpful to fill out a checklist (see below for example) to keep track of where each sample is in the protocol as you work through the various steps. Samples are staggered in such a way that the first sample is one step ahead of the second sample, which is one step ahead of the third sample, and so on.

 [ATACseq_sample_tracker.pdf](#)

7 Label 4 microfuge tubes (the ones that are supplied with the disposable pestles) with sample names and pre-chill on ice.

8 Get the following materials ready; you won't have time to do this in between steps:

- 4 x #22 disposable scalpels
- 4 x 60mm Petri dishes
- 4 x disposable pestles/tubes, along with cordless motor
- 4 x disposable tweezers
- 4 x 40um filters
- 4 x 20um filters
- 4 x BioRad counting slides
- 0.4% Trypan Blue stain
- 200ul wide bore filter tips
- 1000ul wide bore filter tips



Make fresh buffers

9 Prepare  **1 mL fresh lysis buffer** as detailed below, and keep on ice:

 **988 µL EZ lysis buffer**

 **2 µL 5% Digitonin**

 **10 µL 100mM DTT**

Check carefully that there is no precipitate in the 5% digitonin stock solution. If a white precipitate is present, heat at  **94 °C**  **00:05:00**. Re-dissolved digitonin is good for 1 week at room temperature.

10 Prepare  **8.5 mL wash buffer** as detailed below, and keep on ice:

- ▢ 85 μ L 1M Tris-Cl, pH 7.4
- ▢ 17 μ L 5M Sodium chloride
- ▢ 25.5 μ L 1M Magnesium chloride
- ▢ 850 μ L 10% Tween-20
- ▢ 85 μ L 10% BSA
- ▢ 85 μ L 100mM DTT
- ▢ 7352.5 μ L Nuclease-free water

Dispense BSA solution in a laminar flow hood to maintain sterility of the stock solution.

Nuclei isolation 25m

- 11 Place a sterile Petri dish onto the chilled metal plate sitting on ice. Move tube containing first placenta sample from dry ice to wet ice for ~30 seconds, then tap slightly-thawed tissue into pre-chilled Petri dish.
- 12 Cut tissue into small pieces (~ 2 to 3mm in size) using a fresh disposable #22 scalpel. Use disposable tweezers to hold tissue in place while cutting. Transfer tissue into labeled pre-chilled microfuge tube.



Take care when cutting. If tissue is too hard and/or too much pressure is applied with the scalpel, tissue can fly out of the Petri dish. Cut cautiously and slowly to begin with - it is sometimes helpful to rock the blade from side to side to gradually split the frozen tissue. As the tissue thaws it gets a little easier to cut, but if it is allowed to thaw too much then it gets more difficult.


- 13 Add ▢ 125 μ L **lysis buffer** and grind tissue with a disposable pestle attached to a cordless motor.

It should only take a few 5 second pulses to achieve a "smoothie-like" consistency. Do not over-grind.

- 14 Place ground-up tissue on ice for ⌚ **00:05:00** . During this time you can start step 11 for the ^{5m} second placenta sample.
- 15 Centrifuge sample ⌚ **500 x g, 4°C, 00:05:00** . During this time you can start step 11 for the ^{5m} third placenta sample.
- 16 Remove and discard supernatant, and resuspend pellet gently in ^{5m}
▢ **100 µL fresh lysis buffer** using 200ul wide bore filter tips. Incubate on ice for
⌚ **00:05:00** . During this time you can start step 11 for the fourth placenta sample.
- 17 Using a 1000ul wide bore filter tip, add ▢ **500 µL wash buffer** , mix gently, and centrifuge ^{5m}
⌚ **500 x g, 4°C, 00:05:00** .
- 18 Remove and discard supernatant and repeat wash twice, saving the last suspension on ice until all samples reach this step.
- 19 Filter each sample though a 40um filter into a fresh microfuge tube. To maximize recovery of nuclei, wash old microfuge tube with ▢ **200 µL fresh wash buffer** , and pass through the same 40um filter.
- 20 Repeat previous step, this time with a 20um filter. This filter will remove debris that wasn't removed in the previous step.
- 21 Centrifuge ⌚ **500 x g, 4°C, 00:05:00** . ^{5m}
- 22 Remove and discard supernatant carefully, and gently resuspend pellet in
▢ **50 µL tagmentation buffer** using 200ul wide bore filter tip. Clumps of nuclei should break up easily. Place resuspended nuclei on ice.
- 23 Take ▢ **5 µL** aliquot of resuspended nuclei; add ▢ **20 µL tagmentation buffer** and

▢ **25 μ L 0.4% Trypan blue stain** . Flick to mix, then load ▢ **10 μ L** to both sides of a BioRad counting slide. Determine nuclei concentration using an automated cell counter.

Nuclei should stain light blue; not too dark as that would indicate damage to the nuclear membrane. Nuclei should not be clumped together, and there should be little debris.



- 24 Based on nuclei concentration, prepare **10 µL** of a 5,000 nuclei/ul solution (50,000 total nuclei), using tagmentation buffer for dilution. Proceed immediately to tagmentation and library preparation.

It helps to use a pre-formulated Excel spreadsheet to quickly determine the volume of nuclei and tagmentation buffer to use. Several intermediate dilutions are often required to ensure that pipetting volumes are sensible.

Post-steps

5m

- 25 Centrifuge remaining nuclei from step 22 **500 x g, 4°C, 00:05:00**. Remove and discard ^{5m} supernatant, and gently resuspend nuclei in **50 µL Nuclei EZ Storage Buffer** using 200ul wide bore filter tips, making sure to break up clumps of nuclei. Store nuclei at **-80 °C** for up to several months.