




Aug 12, 2022

Transforming yeast


Brian Teague¹¹University of Wisconsin - Stout1 *Works for me* Share

This protocol is published without a DOI.

Yeast ORFans CURE**Brian Teague**
University of Wisconsin - Stout

ABSTRACT

In the same way that we transformed E. coli by making them take up a plasmid, we will transform yeast by inducing them to take up both a plasmid (to cut your genomic location) and some linear DNA to repair that location.



PROTOCOL CITATION

Brian Teague 2022. Transforming yeast. **protocols.io**

<https://protocols.io/view/transforming-yeast-ce78thrw>

KEYWORDS

yeast, transformation, plasmid, dna

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

IMAGE ATTRIBUTION

Masur, Public domain, via Wikimedia Commons

CREATED

Aug 12, 2022

LAST MODIFIED

Aug 12, 2022

PROTOCOL INTEGER ID

68576

PARENT PROTOCOLS

In steps of

[Transforming Yeast \(Instructor Protocol\)](#)

MATERIALS TEXT

Equipment

- Hot water bath (or dry bath or incubator) set to **42 °C**
- Vortexer
- Microcentrifuge
- Incubator set to **30 °C**

Materials

- 1 tube of frozen competent yeast cells, stored at **-80 °C**

▪

[Polyethylene Glycol \(PEG\) 3350](#) **Electron Microscopy**

Sciences Catalog #19760 Step 4

solution, **50 Mass / % volume**

[Lithium Acetate Dihydrate](#) **Sigma**

- **Aldrich Catalog #L4158** Step 4

solution,

1 Molarity (M)

▪

[Salmon Sperm DNA](#) **Research Products International**

Corp Catalog #D52150 Step 4

, **2 mg/mL**

- Plasmid DNA to transform
- URA3 PCR product
- Sterile water
- 1 yeast media plate without uracil
- 1 yeast media plate without leucine

[Glass beads 5 mm](#) **VWR**

- **Scientific Catalog #26396-596**

SAFETY WARNINGS

Lithium acetate can cause serious eye and skin irritation. Wear appropriate PPE, including a lab coat, gloves, and safety glasses.

- 1
 - Compute the volume of your Cas9 plasmid that contains **1 µg** of DNA.
 - Compute the volume of your URA3 PCR product that contains **500 ng** of DNA.
 - If you don't have enough DNA for one or both of these, check with an instructor.

2 Immediately upon retrieving a tube of yeast cells from the $-80\text{ }^{\circ}\text{C}$ freezer, thaw them rapidly by putting them in the $42\text{ }^{\circ}\text{C}$ water bath for 00:00:30

3 Centrifuge $13.000 \times g$, 00:02:00 . Remove and discard the supernatant. 2m

4 Add the following to the tube of yeast cells in order:

- 260 μL

[Polyethylene Glycol \(PEG\) 3350 Electron Microscopy Sciences Catalog #19760](#)

[M]50 Mass / % volume

[Lithium Acetate Dihydrate Sigma](#)

- 36 μL Aldrich Catalog #L4158

[M]1 Molarity (M)

- 50 μL

[Salmon Sperm DNA Research Products International Corp Catalog #D52150](#)

[M]2 mg/mL

- Enough plasmid DNA to equal 1 μg of DNA (computed above)
- Enough PCR product to equal 500 ng of DNA (computed above)






5 Vortex vigorously to resuspend the pellet in the transformation mix.

(This may take from 30 seconds to a minute -- be patient!)

6 Incubate in the $42\text{ }^{\circ}\text{C}$ water bath for 00:30:00

30m

Use a water bath float so you don't have to stand there!

- 7 Centrifuge  **13000 x g, 00:00:30** . Aspirate the supernatant. 30s
- 8 Pipette  **1 mL** of sterile water into the transformation tube. Stir the pellet with a micropipette tip to break up the cell pellet, then vortex to thoroughly resuspend the pellet.
- 9 Pour 5-10 glass beads on a **Uracil dropout** plate. Pipette  **100 µL** of cells onto it. Shake the plate to spread the cells out.
- 10 Pour 5-10 glass beads on a **Leucine** plate. Pipette  **100 µL** of cells onto it. Shake the plate to spread the cells out.
- 11 Incubate 48-72 hours **upside down** at  **30 °C** .