

Jul 21, 2021

# An Optimized Protocol for the Generation of Cortical Neurons from human iPSCs under Defined Conditions

Clare L Parish<sup>1</sup>, Lachlan L Thompson<sup>1</sup><sup>1</sup>The Florey Institute of Neuroscience and Mental Health

1 Works for me

Share

[dx.doi.org/10.17504/protocols.io.bu6znzf6](https://dx.doi.org/10.17504/protocols.io.bu6znzf6)

Benjamin Trist

## ABSTRACT

This protocol is about cortical differentiation.

## ATTACHMENTS

[Cortical Differentiation Protocol.pdf](#)

## DOI

[dx.doi.org/10.17504/protocols.io.bu6znzf6](https://dx.doi.org/10.17504/protocols.io.bu6znzf6)

## PROTOCOL CITATION

Clare L Parish, Lachlan L Thompson 2021. An Optimized Protocol for the Generation of Cortical Neurons from human iPSCs under Defined Conditions. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.bu6znzf6>

## KEYWORDS

passage, FACS, cortical differentiation, differentiation, fix, staining

## LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

May 20, 2021

## LAST MODIFIED

Jul 21, 2021

## OWNERSHIP HISTORY


May 20, 2021  Julia Rossmanith protocols.io

Jun 03, 2021  Benjamin Trist

## PROTOCOL INTEGER ID

50105

## GUIDELINES



## MATERIALS TEXT

### Fix for staining:

A	B	C
d11, d16	OTX2 (rb) 1: 500 PAX6 M) 1:50 Fluorescent reporter (GFP/RFP)	Nestin (m) 1:800 SOX2 (g) 1:300 Fluorescent reporter (GFP/RFP)
d21, d28, d35	OTX2 (rb) 1: 500 PAX6 M) 1:50 Fluorescent reporter (GFP/RFP) Tbr1 (rb): Brn2 (g): Fluorescent reporter (GFP/RFP)	Nestin (m) 1:800 SOX2 (g) 1:300 Fluorescent reporter (GFP/RFP) Tuj1 (m): Tbr1 (rb): Fluorescent reporter (GFP/RFP)

A	B	C
<b>CTX Base media</b> <b>P/S 1X (100U/ml)</b>		
DMEM/F12	48.2ml	241ml
NBM	48.2ml	241ml
B27 + VitA	1ml	5ml
N2	500ul	2.5ml
ITS-A	500ul	2.5ml
GMAX	500ul	2.5ml
P/S	1ml	5ml
B-mercap	90ul	450ul
TOTAL	100ml	500ml

## SAFETY WARNINGS

For hazard information and safety warnings, please refer to the SDS (Safety Data Sheet).

## BEFORE STARTING

### Standard setup

6x48 wells

#### Prepare CTX Base media:

A	B	C
CTX Base media P/S 1X (100 U/ml)		
DMEM/F12	48.2 ml	241 ml
NBM	48.2 ml	241 ml
B27 + VitA	1 ml	5 ml
N2	500 µl	2.5 ml
ITS-A	500 µl	2.5 ml
GMAX	500 µl	2.5 ml
P/S	1 ml	5 ml
B-mercap	90 µl	450 µl
TOTAL	100 ml	500 ml

Day -1: 6h

- 1 Coat plastic cell culture plate (laminin 521 **5 µg/ml**).

6h

Make single cells of hPSCs using accutase for **06:00:00** - **07:00:00** @ **37 °C**.

Seed cells @ **375K/cm<sup>2</sup>**.

d0:

- 2 mTESR1 media change ~1-2 hrs (up to 6 hrs) prior to differentiation initiation.

- 3 

Wash cells with PBS-/- ( **500 µl** for 48 wells, **200 µl** for 96 wells).

- 4 Add **500 µl d0 media** to initiate differentiation.

d0-11: Daily media change

- 5
  - **500 µl CTX Base supplemented with SB (1:1000) LDN (1:2500)**

d11: Passage and fix 2h

- 6
  - 4 wells will be passaged 1:2.5 (10 in total) into 8 new hLaminin-coated wells.
  - Remaining 2 wells will be fixed for staining
  - hLaminin coat 8x 48 wells for passaging and allow > **02:00:00** for coating

2h

d11: Passage 7m

7 

Aspirate media and wash wells with **500 µl/well PBS-/-**.

8 

4m

Add **500 µl EDTA** to each well and incubate for ~ **00:04:00** - **00:08:00**.

9 When clumps are ready to detach, collect into a 15 ml tube with **3 mL PBS-/-**.

*Important note: in these stages of EDTA digestion, cells will be in clumps/rosettes, do not over-triturate suspension as this will break up cells resulting in poor rosette formation.*

10 

Wash each well again with 1x **500 µl EDTA** to ensure all cells have been collected.

11 

3m

Centrifuge at **1100 rpm, 00:03:00**.

12 Aspirate supernatant and gently flick pellet.

13 

Pipette **1 mL CTX base with SB/LDN + Ri (1:1000)** gently to get clumps/cells in suspension.

DO NOT vigorously pipette to single cell. Sometimes just the addition of the 1 ml is enough to distribute cells/clumps.

14 Add remaining CTX base with SB/LDN + Ri (1:1000) volume for **600 µl/well** cell seeding.

15 Aspirate hLam from new plate and add **200 µl CTX base with SB/LDN + Ri (1:1000)**.

16 

Pipette **600 µl cell suspension** into each well.

17 Ensure even distribution by pipetting the peace sign gently.

18



Allow clumps to attach, check on microscope and place gently into incubator.

#### d12: Media change every 2nd day

2h

19

- Aspirate the CTX media and refresh with **1 mL CTX base supplemented with FGF2 (1:5000)**
- From d14, perform media change every second day of CTX + FGF2.

#### d19: Passage and fix

2h

20

- 5 wells will be passaged into 1:2.5; 10x 48 wells for maintenance
- Remaining 2.5 wells will be fixed for staining
- hLam coat 13x 48 wells for passaging and allow > **02:00:00** for coating

2h

NB: At this point in the CTX differentiation protocol, a significant amount of progenitor expansion will have occurred in response to FGF2 exposure. Typically, a 5/6-fold increase in the cell numbers is achieved from Day 0 seeding. Cell numbers from 1 well of a 48-well plate can yield between  $1.5\text{--}2.0 \times 10^6$  neural progenitors. Cells seeded at  $4.5 \times 10^5$  cells/cm<sup>2</sup>. Coat the required number of wells depending on specific application.

#### d19: Passage

7m

21



Aspirate media and wash wells with **500 µl/well PBS-/-**.

22



4m

Add **500 µl Accutase** to each well and incubate for ~ **00:04:00 - 00:08:00**.

23

When cells are ready to detach, collect into a 15 ml tube with **3 mL PBS-/-**.

24



Wash each well again with 1x **500 µl PBS-/-** to ensure all cells have been collected.

25



3m

Centrifuge at **1400 rpm, 00:03:00**.

26

Aspirate supernatant and gently flick pellet.

27



Pipette **1 mL CTX base with FGF2 + Ri (1:1000)** gently to get cells in suspension.

28 Perform cell count and transfer required volume of cell suspension for  $4.5 \times 10^5 \text{ cell/cm}^2$  into a fresh centrifuge tube containing CTX base media supplemented with ROCKi (10  $\mu\text{M}$ ) and FGF2 (20 ng/ml).

29 Aspirate hLam from new plate and add **200  $\mu\text{l}$  CTX base with FGF2 + Ri (1:1000)**.

30 

Pipette **300  $\mu\text{l}$  cell suspension** into each well.

31 Ensure even distribution by pipetting the peace sign gently.

32 

Allow cells to attach, check on microscope and place gently into incubator.

#### d20: Media change every 2nd day

- 33
- Media change with **1 mL CTX base (no factors or small molecules)**.
  - Change media every 2<sup>nd</sup> day

#### d25: Passage and fix 2h

- 34
- Cells are split at a 1:2 ratio. Coat enough plates in accordance with specific application.
  - hLam coat wells for passaging and allow > **02:00:00** for coating

2h

#### d25: Passage 7m

35 

Aspirate media and wash wells with **500  $\mu\text{l}$ /well PBS-/-**.

36 

14m

Add **500  $\mu\text{l}$  Accutase** to each well and incubate for ~ **00:06:00** - **00:08:00**.


37 Add **300  $\mu\text{l}$  CTX media** to each well and transfer cell suspension into a clean 15 ml centrifuge tube.

38 

3m

Centrifuge at  **1400 rpm, 00:03:00** .

39 Aspirate media and flick tube to dislodge pellet. Resuspend cells in (at least 5 ml) CTX base media supplemented with ROCKi.

40 Aspirate hLam from new plate and add  
 **500 µl cell suspension in CTX base media supplemented with ROCKi** .

41 Ensure even distribution by moving the plate in a 'figure 8' pattern.

42  

Allow cells to attach, check on microscope and place gently into incubator.

#### **d26: Media change every 2nd day**

43 ■ Media is replenished every other day.

#### **d29: prepare coverslips for passaging**

44 Remove sterile coverslips from Ethanol using sterile tweezers.

45 Place coverslips into 48-well plates and wash x 3 with 500 µL PBS-/-, with 5 min incubations, at Room temperature to remove all residual Ethanol.

46 Add 300 µL of Poly-L-Ornithine (PO) to each well and incubate overnight at 4°C. Each well should contain a single coverslip.

#### **d30: passage**

47 Remove 48-well plate, containing PO-coated coverslips from the fridge and wash x 2 with 5-minute incubations at RT.

48 Add 190 µl of 5 µg/ml of L521 solution diluted in 1 x PBS +/- to each well with a coverslip and incubate for 2 hours at 37°C.

49 Dissociate cells as performed at D25 timepoint.

50 Add 300 µl of Accutase to each of the wells and place in the incubator for 10 minutes at 37°C.

- 51 Following incubation, add 300 µL of CTX media to each well, harvest cells and transfer cell suspension into a clean 15 ml centrifuge tube.
- 52 Centrifuge cell suspension at 1,400 rpm for 4 minutes.
- 53 Aspirate media and flick tube to dislodge pellet. Resuspend cells in (at least 5 ml) CTX base media supplemented with ROCKi, 10 µM.

#### d37-55: maturation

- 54 Media is aspirated and replaced with maturation media (NBB27 media) supplemented with BDNF (40 ng/ml), GDNF (40 ng/ml), dcAMP (0.05 mM), AA (200 nM) and laminin (1 µg/ml; Sigma-Aldrich).
- 55 Cortical neurons are maintained in this media with growth factors until Day 55, where a final immunocytochemical analysis for cortical patterning is performed.

#### Fix: 8m

- 56 

Wash wells with 200 µL PBS-/-.

- 57 Add  200 µl 4% PFA for  00:08:00 .

8m

- 58 Aspirate PFA.

- 59 

Wash 3x with PBS-/-.

59.1 Wash with PBS-/-.

59.2 Wash with PBS-/-.

59.3 Wash with PBS-/-.



60 Fill wells with PBS/Azide or begin staining.

Staining: 20m

61 Permeabilize cells with  200 µl -  500 µl PADT for  00:20:00 -  00:30:00 .

20m

62  

Aspirate PADT and pipette on primary antibodies to incubate  Overnight at  4 °C .

63 Aspirate primaries.

64 

Wash wells 2x with PBS-/-.

64.1 Wash wells with PBS-/-.

64.2 Wash wells with PBS-/-.

65 Block with  500 µl PADT for  00:30:00 .

30m

66 

1h

Aspirate PADT and pipette on secondaries to incubate for  01:00:00 -  01:30:00 at room temperature.

67 

Aspirate secondaries and wash 2x with  500 µl PBS-/- .

67.1 Wash with  500 µl PBS-/- .

67.2 Wash with  500 µl PBS-/- .

68 DAPI stain for  00:08:00 -  00:10:00 .

69 

Wash DAPI with 1-2x PBS/- wash.

69.1 Wash DAPI with PBS/- wash.

69.2 Optional: Wash DAPI with PBS/- wash.

70 Fill well with PBS/Azide for long-term storage.

71 

Image.