



Dec 14, 2020

## 1.2 Bradford assay

Elizabeth Fozo<sup>1</sup><sup>1</sup>In-house protocol

1

Works for me

This protocol is published without a DOI.



Eadewunm

### ABSTRACT

Protein concentration estimation using Bradford assay

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

<http://www.bio-rad.com/webroot/web/pdf/lsr/literature/LIT33.pdf>

### PROTOCOL CITATION

Elizabeth Fozo 2020. 1.2 Bradford assay. **protocols.io**  
<https://protocols.io/view/1-2-bradford-assay-bqqqmvmvw>

MANUSCRIPT CITATION please remember to cite the following publication along with this protocol

<http://www.bio-rad.com/webroot/web/pdf/lsr/literature/LIT33.pdf>

### KEYWORDS

Bradford assay

### LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

### CREATED

Dec 14, 2020

### LAST MODIFIED

Dec 14, 2020

### PROTOCOL INTEGER ID

45552

### DISCLAIMER:

DISCLAIMER: THIS WORK IS IN PROGRESS. IT IS FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](#) is not peer-reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with [protocols.io](#), can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

### ABSTRACT

Protein concentration estimation using Bradford assay

#### Protein concentration estimation using Bradford assay

- 1 Prepare working reagent by diluting 1 part Dye Reagent Concentrate with 4 parts distilled, deionized water. This diluted reagent may be used for approximately 2 weeks when kept at room temperature.
- 2 Prepare dilutions of a protein standard (Get from Becker lab- Sarah, they have aliquots for the diluted standards). The linear range of the assay for BSA is 0.2 to 0.9 mg/ml
- 3 Add 50ul of each standard into the cuvettes and add 2.5 ml of Bradford reagent (optimized from 100 ul + 5 ml)
- 4 Dilute your samples with 0.1% SDS in the cuvettes. For 1:10 dilution add 45ul of 0.1% SDS to the cuvette and add 5 ul of your sample (40 ul of 0.1% SDS and 10 ul sample for 1:5), mix properly. Add 2.5 mL of Bradford reagent to each tube and mix.
- 5 Incubate at room temperature for 5 minutes. Absorbance will increase over time; samples should incubate at room temperature for no more than 1 hour.
- 6 Measure absorbance at 595 nm.
- 7 Plot the standard curve using standards and measure OD and results for the samples. Multiply the result with the dilution factor to get the concentration of your protein sample.
- 8 Use an excel sheet and calculate the amount of you sample to be mixed with the appropriate sample buffer to get a final concentration of 0.5 ug/ul.