



AUG 03, 2023

OPEN ACCESS



DOI:
dx.doi.org/10.17504/protocols.io.8epv5jz54l1b/v1

Protocol Citation: Koki Tohara 2023. Mouse genetic models to manipulate enterochromaffin cell activity - Murine Organoid ELISA.
protocols.io
<https://dx.doi.org/10.17504/protocols.io.8epv5jz54l1b/v1>


MANUSCRIPT CITATION:
 Gut Enterochromaffin Cells are Critical Drivers of Visceral Pain and Anxiety
 (https://www.biorxiv.org/content/10.1101/2022.04.04.486775v1)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

🌐 Mouse genetic models to manipulate enterochromaffin cell activity - Murine Organoid ELISA

Koki Tohara¹

¹University of California, San Francisco



Koki Tohara
 University of California, San Francisco

ABSTRACT

This protocol describes how we maintain murine intestinal organoids and perform the enzyme-linked immunosorbent assay (ELISA) for serotonin (5-HT).

MATERIALS

DPBS, no calcium, no magnesium (Thermo Fisher #14190144)
 DPBS, no calcium, no magnesium + 2 mM EDTA
 Advanced DMEM/F-12 (Thermo Fisher #12634010)
 1 M HEPES (Thermo Fisher #15630080)
 GlutaMAX™ Supplement (Thermo Fisher #35050061)
 Penicillin-Streptomycin (10,000 U/mL) (Thermo Fisher #15140163)
 Recombinant Murine Noggin (Peprotech #250-38)
 Mouse EGF Recombinant Protein (Thermo Fisher #PMG8041)
N-Acetyl-L-cysteine (Millipore-Sigma #A7250-5G)
 B-27™ Supplement (50X), serum-free (Thermo Fisher #17504044)
 R-spondin 2 (supernatant from R-spondin expressing HEK293 cells)
 Corning® Matrigel® Matrix (Corning #356231)
 Bovine Serum Albumin (Millipore-Sigma #A9418)
 Serotonin ELISA kit (LDN #BA-E-8900)
 24-well plate (Falcon® 24-Well Plate #38021)
 Corning™ Falcon™ Cell Strainers 70 µm (Corning #08-771-2)
 Ringer's solution (140 mM NaCl, 5 mM KCl, 2 mM CaCl₂, 2 mM MgCl₂, 10 mM D-glucose, and 10 mM HEPES-Na [pH 7.4]))

Following four mouse lines were used.

RC::FL-hM3Dq (Jackson Labs, Strain #026942)
 Tac1Cre (Jackson Labs, strain #021877)
 Pet1Flp (gift of Dr. Susan Dymecki)
 RC::PFTox (gift of Dr. Susan Dymecki)

Protocol status: Working
We use this protocol and it's working

Created: Oct 12, 2022

Last Modified: Aug 03, 2023

PROTOCOL integer ID:
71277

Keywords: Intestinal organoid, Serotonin, 5-HT, ELISA

Preparation or basal media

10m

- 1 Add 5 mL of GlutaMAX, 1 M HEPES, and Penicillin-Streptomycin to 500 mL of Advanced DMEM/F-12 and filter through 0.22 μ m. This basal media can last up to 1 month.

Preparation of complete organoid culture media

- 2 Add 800 μ L of B-27, Recombinant Murine Noggin (final conc. 100 ng/mL), Mouse EGF Recombinant Protein (final conc. 500 ng/mL), *N*-Acetyl-L-cysteine (final conc. 0.5 mM), and 4 mL of R-spondin 2 to 40 mL of basal media. This complete media can last up to 2 weeks.

Mouse organoid prep

10m

- 3 Organoids were generated from male Tac1-Cre(+/-);Pet1Flp(+/-);RC::FL-hM3Dq(+/-) or Tac1-Cre(+/-);Pet1Flp(+/-);RC::OFTox(+/-) animals (5-8 weeks old).
- 4 Euthanize the animal under CO₂ and spray belly with 70% EtOH.
- 5 Isolate the whole small intestine and divide the tissue into three pieces. Use the middle piece (jejunum) for organoid generation.

- 6 Filet open the intestine along the mesentery.
- 7 Scrape off the villi using a glass coverslip. Scrape well so that all the villi come off.
- 8 Cut the intestine into 2-4 mm pieces and move to a 50 mL falcon tube.
- 9 Add 10 mL ice-cold DPBS. Pipette up and down with 10 mL pipettes.
- 10 Wait for a couple of minutes until the tissue settles down and discard the supernatant.
- 11 Repeat the step 9-10 until the supernatant becomes clear (usually 4-5 times).
- 12 Add 30 mL of ice-cold DPBS + 2 mM EDTA and rock the tube at 4°C for 30 minutes.
- 13 While waiting, prepare 6x 50 mL falcon tubes and label them No.1-6.

- 14** Wait for a couple of minutes until the tissue settles down and discard the supernatant.
- 15** Add 10 mL ice-cold DPBS. Pipette up and down with 10 mL pipettes.
- 16** Wait for a couple of minutes until the tissue settles down and move the supernatant to falcon tube #1 prepared at step #13.
- 17** Repeat the step 15-16 five times. Note that the supernatant should be put into separated tubes prepared at step #13 (tube #2-6).
- 18** Place 20 μ L of each supernatant fraction on a glass slide and check under a microscope.
- 19** Identify fractions that contain crypts.
- 20** Combine crypts-containing fractions and filter through a 70 μ m strainer.
- 21** Spin down at 900 g for 5 minutes and discard the supernatant.
- 22** Resuspend pellet with 10 mL ice-cold basal media and spin down at 900 g for 5 minutes.

- 23** Carefully remove all the supernatant.
- 24** Repeat the step #22 and #23.
- 25** Resuspend the pellet into 400 μ L ice-cold matrigel.
- 26** Pipette 50 μ L of the matrigel/organoid suspension into the center of each of four wells of a prewarmed 24-well plate to form domes in the center of each well.
- 27** Place the lid on the culture plate and quickly turn the plate upside down.
- 28** Incubate at 37°C for 10 minutes to set the matrigel.
- 29** Turn the plate back upside so that you can apply media to wells. Gently add 500 μ L complete media to each well.
- 30** Add basal media to surrounding wells to prevent the wells from drying.

Maintenance of intestinal organoids

46m

- 31** Put an aliquot of matrigel in fridge at 4°C. Put a 24-well plate to a 37°C incubator. 1m
- 32** Put desired volume of complete media at 37°C 1m
- 33** Aspirate media from wells 2m
- 34** Add 500 µL cold basal media to a well and let it sit for about 15 seconds 2m
- 35** Pipette up and down to dissolve matrigel, collect organoids, and move to a 15 mL falcon tube 3m
- 36** Add 500 µL cold basal media to a well, collect residual organoids, and move to the same 15 mL falcon tube 3m
- 37** Place a 200 µL tip inside a 1000 µL tip and pipette up and down ~25 times to break up organoids 3m
- 38** Add 10 mL cold basal media to the falcon tube and mix well 1m

- | | | |
|-----------|---|-----------|
| 39 | Centrifuge at 200 g for 3 min at 4°C | 3m |
|
 | | |
| 40 | Aspirate the supernatant (be careful not to aspirate the organoid pellet) | 2m |
|
 | | |
| 41 | Add 10 mL cold basal media to the falcon tube and mix well | 1m |
|
 | | |
| 42 | Centrifuge at 200 g for 3 min at 4°C | 3m |
|
 | | |
| 43 | Aspirate the supernatant (be careful not to aspirate the organoid pellet) | 1m |
|
 | | |
| 44 | Add 200 µL of ice-cold matrigel to the falcon tube and mix the organoids well | 2m |
|
 | | |
| 45 | Pipette 50 µL of the matrigel/organoid suspension into the center of each of four wells of a prewarmed 24-well plate to form domes in the center of each well | 3m |
|
 | | |
| 46 | Place the lid on the culture plate and quickly turn the plate upside down | 1m |

- 47** Incubate at 37°C for 10 minutes to set the matrigel 10m
- 48** Turn the plate back upside so that you can apply media to wells. Gently add 500 µL complete media to each well 2m
- 49** Add basal media to surrounding wells to prevent the wells from drying 2m
- 50** Exchange the media every 3-4 days or whenever the media turns yellow.
- 51** Split again in 6-8 days.

Preparation of organoids for ELISA

38m

- 52** Organoids are grown for 4-6 days
- 53** Aspirate the supernatant (be careful not to aspirate the organoid pellet) 2m
- 54** Add 500 µL of cold DPBS to a well and let it sit for about 15 seconds 1m

- | | | |
|-----------|--|----|
| 55 | Gently pipette up and down ~10 times with P-1000 to remove matrigel (but not to break the organoids up) and move to a 15 mL falcon tube. | 2m |
| 56 | Add 500 µL of cold DPBS to the tube | 2m |
| 57 | Put the tube on ice for 5 min | 5m |
| 58 | Gently pipette up and down ~10 times with P-1000 | 2m |
| 59 | Add 10 mL cold DPBS and spin down at 200 g for 3 min at 4°C | 3m |
| 60 | Aspirate the supernatant as much as possible | 1m |
| 61 | Add 1 mL cold DPBS and gently pipette up and down ~10 times with P-1000 | 2m |
| 62 | Add 10 mL cold DPBS and spin down at 200 g for 3 min at 4°C | 3m |
| 63 | Aspirate the supernatant as much as possible | 1m |

- 64** Repeat the step #33-35 twice to completely wash the matrigel out 12m
- 65** Add 100 μ L DPBS + 0.1% BSA and move organoids to an eppendorf tube (it is important to add 0.1% BSA to prevent organoids from sticking to tubes) 1m
- 66** Remove the supernatant and add 100 μ L Ringer's (140 mM NaCl, 5 mM KCl, 2 mM CaCl_2 , 2 mM MgCl_2 , 10 mM D-glucose, and 10 mM HEPES-Na [pH 7.4])) buffer for ELISA assay (do not lyse the organoids as you are measuring released serotonin). 1m

5-HT ELISA

- 67** ELISA is performed according to the manufacturer's protocol (https://www.ldn.de/wp-content/uploads/IFU-BA-E-5900R-V15.0_wz.pdf).