



May 16, 2020

Enzymatic PCR Cleanup Protocol (NEB #M0525)

New England Biolabs¹¹New England Biolabs

1 Works for me dx.doi.org/10.17504/protocols.io.7r4hm8w

New England Biolabs (NEB)

Tech. support phone: +1(800)632-7799 email: info@neb.com

New England Biolabs
New England Biolabs

ABSTRACT

Enzymatic PCR Cleanup Protocol (NEB #M0525) uses Quick CIP; Quick CIP is a heat-labile version of calf intestinal alkaline phosphatase (CIP) purified from a recombinant source.

Features

- Rapid and irreversible heat inactivation eliminates unwanted activity
- Improved storage stability versus native enzyme
- Faster reaction setup (no supplemental additives like zinc required) and shorter incubation time
- Flexible reaction conditions (active in any restriction enzyme buffer, no clean-up required)
- Less enzyme required (high specific activity), resulting in a lower cost per reaction
- No need for multiple phosphatases (Quick CIP removes 5'- and 3'- phosphates from DNA, RNA and dNTPs)
- Active on unincorporated dNTPs in PCR products - improves DNA sequencing and SNP analysis
- Recombinant for purity, consistency and value

EXTERNAL LINK

<https://neb.com/protocols/2019/06/04/enzymatic-pcr-cleanup-protocol-neb-m0525>

MATERIALS

NAME	CATALOG #	VENDOR
Quick CIP	M0525	New England Biolabs
Exonuclease I (E. coli)	M0293	New England Biolabs
Shrimp Alkaline Phosphatase (rSAP)	M0371	New England Biolabs

SAFETY WARNINGS

Please see SDS (Safety Data Sheet) for hazards and safety warnings.



- 1 Add 1 µl of Quick CIP ([NEB #M0525](#)) to 5 µl of PCR product.



Note that 0.5 µl of Exo I ([NEB #M0293](#)) and 1 µl of rSAP ([NEB #M0371](#)) can be used in the same volume.



- 2 Incubate the mix at the 37 °C for 00:15:00 .

- 3 Inactivate both enzymes at 80 °C for 00:15:00 .

4 PCR products are ready for downstream application.