



Sep 28, 2021

跑膠

Yin-Tse Huang¹, Tina¹¹KMU

In Development



Share

This protocol is published without a DOI.

tina

ABSTRACT

electrophoresis

PROTOCOL CITATION

Yin-Tse Huang, Tina 2021. 跑膠. protocols.io
<https://protocols.io/view/protocol-bx6nprde>

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Sep 13, 2021

LAST MODIFIED

Sep 28, 2021

PROTOCOL INTEGER ID

53166

1 80 mL TBE + 320 mL DDW 配成1倍的buffer

2 大片的膠大約 30 mL

小片的膠大約 20 µl

Agarose 1% 大塊的膠 表示加 0.3 g 和 30 mL DDW

2.1 把配好的Gel 1% 放入微波爐加熱 00:01:30 , 每 00:00:30 一次

2m

3 將Gel倒入製膠模型, 若有氣泡, 可以吸管尖挑開, 之後架上齒梳, 等 00:30:00

30m

4 輕輕將齒梳向上拉開 (注意要慢慢拉開, 否則要加DNA樣品的孔洞會破掉)

5 將膠體連同膠模一起放入電泳槽中

6 加入  350 mL TBE buffer 在電泳槽中

7  3 μ l DNA  1 μ l dye 混和滴在 parafilm 上

8