



SEP 23, 2023

OPEN ACCESS



**DOI:**  
[dx.doi.org/10.17504/protocols.io.n92ldmyynl5b/v1](https://dx.doi.org/10.17504/protocols.io.n92ldmyynl5b/v1)

**Protocol Citation:** Elias Adriaenssens 2023. Rapalog-induced chemical dimerization experiments. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.n92ldmyynl5b/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working  
 We use this protocol and it's working

**Created:** Jun 28, 2023

**Last Modified:** Sep 23, 2023

**PROTOCOL integer ID:**  
 84159

**Keywords:** Rapalog-induced chemical dimerization experiments

## Rapalog-induced chemical dimerization experiments

Elias Adriaenssens<sup>1</sup>

<sup>1</sup>Sascha Martens lab, University of Vienna, Max Perutz Labs - Vienna

ASAP Collaborative Research Network

Sascha Martens lab, University of Vienna, Max Perutz Labs - Vienna



Elias Adriaenssens

Sascha Martens lab, University of Vienna, Max Perutz Labs - ...

### ABSTRACT

This protocol describes Rapalog-induced chemical dimerization experiments.

### ATTACHMENTS

[760-1928.pdf](#)

### MATERIALS

#### Reagents

iDimerize™ Inducible Heterodimer System Takara Bio Inc. Catalog #635067

LSR Fortessa Cell Analyzer (BD Biosciences)

## Rapalog-induced chemical dimerization experiments

1d

- To perform chemical-induced dimerization (CID) experiments we will first establish cell lines expressing the FRB-Fis1 and FKBP-GFP-GOI (FKBP-GFP fused to our gene of interest).

- 2 HeLa cells are consecutively transduced with lentivirus for FRB-Fis1, which also expresses mitochondrially targeted monoKeima (mt-mKeima), and then with FKBP-GFP-GOI. See the lentivirus transduction protocol for further details.
- 3 Once cell lines are established, cells are seeded in 6-well plates the day before the experiment (density: 800k per well).
- 4 Induce FKBP-FRB dimerisation by treating the cells with **100 nM** 50 nanomolar (nM) Rapalog A/C hetero-dimerizer rapalog (635057, Takara) for **24:00:00** . 1d
- 5 After 24h, collect cells by trypsinization and take the cells to the FACS facility for mt-mKeima analysis.
- 6 Analyze the cells by flow cytometry, using the mt-mKeima ratio (561/405) as a readout with an LSR Fortessa Cell Analyzer (BD Biosciences).
- 7 Gate the cells for GFP/mt-mKeima double-positive cells, and collect 10,000 events for this population.
- 8 Analyze data using FlowJo.

