



Aug 31, 2021

Plasmid construction

Chunmei Chang¹¹Team Hurley, University of California, Berkeley

1 Works for me

Share

dx.doi.org/10.17504/protocols.io.bxrjpm4n

Chunmei Chang

ABSTRACT

plasmid construction

DOI

dx.doi.org/10.17504/protocols.io.bxrjpm4n

PROTOCOL CITATION

Chunmei Chang 2021. Plasmid construction. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bxrjpm4n>



LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Aug 26, 2021

LAST MODIFIED

Aug 31, 2021

PROTOCOL INTEGER ID

52747

MATERIALS TEXT

Q5® High-Fidelity DNA Polymerases (NEB)
restriction enzymes (NEB)
quick CIP (NEB)
T4 ligase (NEB)
Gel extraction kit (Bio Basic)

- 1 Amplify the insert gene fragment by PCR with primers including 21 nt of overlapping sequence with the target gene.
- 2 Linearize the backbone with restriction enzymes (NEB).
- 3 Treat linearized backbone with quick CIP (NEB).

- 4 Run PCR products and linearized backbone in an agarose gel to confirm the size.
- 5 Purify the DNA from gel using a Gel extraction kit (Bio Basic).
- 6 Ligate the linearized backbone and the insert with the T4 ligase (NEB).
- 7 Transform the ligation products into home made competent cells.
- 8 Perform colony PCR to screen for colons that with inserted gene.
- 9 Sequence to verify that the inserted gene is correct.