



Version 2

Jan 25, 2021

# Adipose depot innervation: whole mount staining, imaging, quantification V.2

Jake W Willows<sup>1</sup>, Kristy Townsend<sup>2</sup>, Magdalena Blaszkiwicz<sup>1</sup>

<sup>1</sup>Ohio State University, Columbus; <sup>2</sup>Ohio State University

1

Works for me

[dx.doi.org/10.17504/protocols.io.brs2m6ge](https://dx.doi.org/10.17504/protocols.io.brs2m6ge)

Townsend Lab Neurobiology & Energy Balance

Tech. support phone: +12075812541 email: [kristy.townsend@maine.edu](mailto:kristy.townsend@maine.edu)



Kristy Townsend

Ohio State University, University of Maine

## ABSTRACT

### Abstract

Little is known about the diversity and function of adipose tissue nerves due, in part, to the inability to effectively visualize the various nerve subtypes residing within these tissues. The tools currently available for researchers to image and quantify adipose tissue innervation are limited and dependent on optical clearing techniques and light sheet microscopy. Here we present a method of tissue processing that uses a method of mechanically compressing tissue to decrease tissue thickness in the z-axis by expanding it in the x and y-axes whilst leaving cells intact. This has been combined with autofluorescence quenching techniques to permit imaging of intact whole tissues on both widefield and confocal microscopes and a complementary means to perform whole tissue neurite density quantification. We have included examples of how this technique can be used to further our current knowledge of adipose-nerve communication by characterizing the nerves, nerve-subtypes, and neurovascular interactions within the inguinal subcutaneous white adipose tissue in mice.

## THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Blaszkiwicz M, Willows JW, Dubois AL, et al. (2019) Neuropathy and neural plasticity in the subcutaneous white adipose depot. PLoS ONE 14(9): e0221766. <https://doi.org/10.1371/journal.pone.0221766>

Willows JW, Blaszkiwicz M, Lamore A, et al. (2019) Visualization and Analysis of Whole Depot Adipose Tissue Innervation (pre-print) bioRxiv 788885; <https://www.biorxiv.org/content/10.1101/788885v1>

## ATTACHMENTS

1 Whole Mount scWAT Isolation.pdf

2 Whole Mount Adipose Immunostaining.pdf

3 Confocal Imaging of Whole Mount Adipose Tissue.pdf

4 Whole Mount Isolectin Staining of Adipose Protocol.pdf

## DOI

[dx.doi.org/10.17504/protocols.io.brs2m6ge](https://dx.doi.org/10.17504/protocols.io.brs2m6ge)

## PROTOCOL CITATION

Jake W Willows, Kristy Townsend, Magdalena Blaszkiwicz 2021. Adipose depot innervation: whole mount staining, imaging, quantification. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.brs2m6ge>

Version created by Kristy Townsend

MANUSCRIPT CITATION please remember to cite the following publication along with this protocol

Blaszkiwicz M, Willows JW, Dubois AL, et al. (2019) Neuropathy and neural plasticity in the subcutaneous white adipose depot. PLoS ONE 14(9): e0221766. <https://doi.org/10.1371/journal.pone.0221766>

Willows JW, Blaszkiwicz M, Lamore A, et al. (2019) Visualization and Analysis of Whole Depot Adipose Tissue Innervation (pre-print) bioRxiv 788885; <https://www.biorxiv.org/content/10.1101/788885v1>

#### KEYWORDS

white adipose tissue (WAT), adipose innervation, adipose neuropathy, aging, obesity, diabetes, metabolic health, exercise, cold, peripheral neuropathy, diabetic peripheral neuropathy, DPN

#### LICENSE

\_\_\_\_\_ This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

#### CREATED

Jan 25, 2021

#### LAST MODIFIED

Jan 25, 2021

#### PROTOCOL INTEGER ID

46650

#### GUIDELINES

**See attachments under Abstract for PDF copies of protocols.**

#### MATERIALS TEXT

**See attachments under Abstract for PDF copies of protocols.**

#### DISCLAIMER:

Laboratory of Kristy Townsend at University of Maine - Neurobiology & Energy Balance Laboratory  
Public Protocols (please cite if you download and utilize our protocols)

**1 See attachments under Abstract for PDF copies of protocols.**