



Sep 01, 2021

## FindingNemo Extraction 4: Circulomics Kit

Inswasti Cahyani<sup>1</sup>, John Tyson<sup>2</sup>, Nadine Holmes<sup>1</sup>, Josh Quick<sup>3</sup>, Nicholas Loman<sup>3</sup>,  
Matthew Loose<sup>1</sup>

<sup>1</sup>University of Nottingham; <sup>2</sup>University of British Columbia; <sup>3</sup>University of Birmingham

1 Works for me

Share

[dx.doi.org/10.17504/protocols.io.bxvwpn7e](https://dx.doi.org/10.17504/protocols.io.bxvwpn7e)

Inswasti Cahyani

### ABSTRACT

This is a sub-protocol designed to extract/isolate ultra-high molecular weight (UHMW) DNA to obtain ultra-long (UL) reads on Nanopore sequencers using the **Circulomics Nanobind CBB Big DNA Kit**.

A DNA extraction protocol that yields clean and homogeneous UHMW DNA is important for a good UL sequencing output. The choice of protocol should be based on achieving these parameters.

The **Circulomics Nanobind CBB Big DNA Kit** was initially the recommended extraction protocol for UL sequencing using SQK-ULK001.

We tested this sub-protocol in **human cell line**, with standard input cells of 6 millions. As a rule of thumb, a million cells will suffice for one load on a MinION.

### DOI

[dx.doi.org/10.17504/protocols.io.bxvwpn7e](https://dx.doi.org/10.17504/protocols.io.bxvwpn7e)

### PROTOCOL CITATION

Inswasti Cahyani, John Tyson, Nadine Holmes, Josh Quick, Nicholas Loman, Matthew Loose 2021.  
FindingNemo Extraction 4: Circulomics Kit. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.bxvwpn7e>

### KEYWORDS

ultra-long sequencing, cohex, glass bead, nanopore, MinION, UHMW DNA, Monarch, Circulomics, phenol, SDS, CTAB, GM12878, Whatman, PromethION, Nanobind

### LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

### CREATED

Sep 01, 2021

### LAST MODIFIED

Sep 01, 2021

### PROTOCOL INTEGER ID

52886

## GUIDELINES

### Acknowledgements

The FindingNemo protocol series was developed by Inswasti Cahyani for the Long Read Club with significant contributions from John Tyson and Nadine Holmes, also Josh Quick, and continuous discussion and support of Matt Loose and Nick Loman. We would also like to thank Giron Koetsier (NEB) and Kelvin Liu (Circulomics) for lending their expertise and advance product trials.

Please follow on Twitter for latest updates and results:

@NininUoN

@mattloose

## MATERIALS TEXT

### Chemicals/Compounds

[☒ Tris-HCl pH 8.0 Thermo](#)

- [Scientific Catalog #J22638-AE](#)

[☒ Ethanol](#)

- [Absolute Honeywell Catalog #32221-2.5L](#)

[☒ Isopropanol Absolute Fisher](#)

- [Scientific Catalog #P/7500/15](#)

[☒ 1X Phosphate Buffer Saline Fisher](#)

- [Scientific Catalog #15453819](#)

[☒ Nuclease-free](#)

- [Water Thermofisher Catalog #AM9920](#)

### Kits

[☒ Nanobind CBB Big DNA](#)

- [Kit Circulomics Catalog #NB-900-001-01](#) (optional)

### Disposables

[☒ DNA LoBind Tubes, 1.5](#)

- [mL Eppendorf Catalog #0030108051](#)

[☒ DNA LoBind 2.0ml PCR Clean Eppendorf](#)

- [Tubes Eppendorf Catalog #0030 108.078](#)

[☒ Protein LoBind Tubes 1.5](#)

- [mL Eppendorf Catalog #0030108442](#)

- Wide-bore (or cut off) P1000 and P200 tips

#### SAFETY WARNINGS

When handling phenol always wear PPE, keep a solution of 50% (w/v) PEG-400 nearby to treat the burn in the case of accidental splashes.

#### BEFORE STARTING

##### Things to observe at all times:

- Excessive and vigorous pipetting and vortexing should be avoided as these may shear the DNA.
- Make up buffers with nuclease-free water to avoid introducing nucleases to solutions.
- Avoid unnecessary heating and freezing; isolated DNA should be stable for storage in the fridge for months.

#### UHMW DNA Extraction

- 1 This protocol is using the **Nanobind CBB Big DNA Kit** (see Materials).
- 2 Follow the manufacturer's instructions as described [here](#).
- 3 After overnight incubation at room temperature, quantify the DNA as per "**UHMW DNA QC**" and check homogeneity by calculating %CV values.
- 4 Store at 4°C or continue to **UL Library Preparation** as per Section "**Modified ULK001**".  
If only SQK-RAD004 is available, follow library preparation as in Section "**Modified RAD004**" or "**KrazyStarFish (KSF)**".  
**4 °C for storage**