



JAN 27, 2023

OPEN ACCESS

**DOI:**  
[dx.doi.org/10.17504/protocols.io.8epv59objg1b/v1](https://doi.org/10.17504/protocols.io.8epv59objg1b/v1)

**External link:**  
<https://doi.org/10.1371/journal.pone.0279681>

**Protocol Citation:** Diana A Tapia-Sidas, Brenda Vargas-Hernández, José Abrahán Ramírez-Pool, Leandro A Nuñez-Muñoz, Berenice Calderón-Pérez, Rogelio González-González, Luis Gabriel Briebe, Rosalía Lira-Carmona, Eduardo Ferat-Osorio, Constantino López-Macías, Roberto Ruiz-Medrano, Beatriz Xoconostle-Cazares 2023. Expression and purification of recombinant MM4 reverse transcriptase (RT). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.8epv59objg1b/v1>

**MANUSCRIPT CITATION:** D. Tapia-Sidas, B. Vargas-Hernández, J. Ramírez-Pool, L. Nuñez-Muñoz, B. Calderón-Pérez, R. González-González, L. Briebe, R. Lira-Carmona, E. Ferat-Osorio, R. Ruiz-Medrano, B. Xoconostle-Cázares. (2022). Starting from scratch: step-by-step development of a diagnostic test for SARS-CoV-2 detection by RT-LAMP. *PLoS One*.

## 🌐 Expression and purification of recombinant MM4 reverse transcriptase (RT)


Diana A Tapia-Sidas<sup>1</sup>, Brenda Vargas-Hernández<sup>1</sup>, José Abrahán Ramírez-Pool<sup>1</sup>, Leandro A Nuñez-Muñoz<sup>1</sup>, Berenice Calderón-Pérez<sup>1</sup>, Rogelio González-González<sup>1</sup>, Luis Gabriel Briebe<sup>2</sup>, Rosalía Lira-Carmona<sup>3</sup>, Eduardo Ferat-Osorio<sup>4</sup>, Constantino López-Macías<sup>4</sup>, Roberto Ruiz-Medrano<sup>1</sup>, Beatriz Xoconostle-Cazares<sup>1</sup>

<sup>1</sup>Departamento de Biotecnología y Bioingeniería, Centro de Investigación y de Estudios Avanzados, Mexico City, Mexico;

<sup>2</sup>Laboratorio Nacional de Genómica para la Biodiversidad, Centro de Investigación y de Estudios Avanzados, Irapuato, Guanajuato, Mexico;

<sup>3</sup>Unidad de Investigación Médica en Enfermedades Infecciosas y Parasitarias, Hospital de Pediatría, Centro Médico Nacional Siglo XXI, Instituto Mexicano del Seguro Social, Mexico City, Mexico;

<sup>4</sup>División de Investigación en Salud, UMAE Hospital de Especialidades “Dr. Bernardo Sepúlveda Gutiérrez”, Centro Médico Nacional Siglo XXI, Instituto Mexicano del Seguro Social, Mexico City, Mexico



Diana A Tapia-Sidas

### ABSTRACT

Reverse transcriptases (RTs) are RNA-dependent DNA polymerases able of synthesizing DNA (complementary DNA or cDNA) from an RNA template. RTs are especially useful in RNA-based nucleic acid detection techniques. Due to its high catalytic activity and fidelity, one of the most widely used RTs in diagnostics and molecular biology is the RT from the Moloney Murine Leukemia Virus (MMLV). However, RT-MMLV is thermally unstable, so previous studies have produced a RT variant called RT-MM4 carrying mutations of positive charges in four amino acids (E286R/E302K/L435R/D524A). This protocol describes the optimized expression process, as well as the FPLC purification of RT-MM4 for use in isothermal amplification techniques, such as end-point colorimetric or real-time fluorometric RT-LAMP.

### GUIDELINES

During the process of protein purification maintain all samples that contain the protein of interest in a cold environment to avoid protein degradation.

### MATERIALS

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**Created:** Jun 08, 2022


**Last Modified:** Jan 27, 2023

**PROTOCOL integer ID:**  
64229

**Keywords:** Moloney Murine Leukemia Virus, thermostable reverse transcriptase, protein expression, protein purification, MM4 RT

## Reagents:

-  Chaperone Plasmid Set **Takara Bio USA, Inc. Catalog #3340**

-  Chemically Competent E. coli One Shot™ BL21(DE3) **Invitrogen - Thermo Fisher Catalog #C600003**

-  100 ng/μL pET-MM4-RT plasmid Step 3.1


- Stock  30 mg/mL chloramphenicol

- Stock  100 mg/mL carbenicillin

- Stock  1 Molarity (M) IPTG

- Stock  0.5 Molarity (M) EDTA

- Stock  1 Molarity (M) DTT

- Stock  10 % (v/v) Triton X-100

-  250 μL of SOC medium Step 1.4

-  8 mL of TFBII solution Step 2.4

-  2.5 mL of TFBII solution Step 2.5

- Cryotubes with  500 μL of 30% v/v glycerol (sterilized)

- LB agar plates

- Tubes with  3 mL Luria-Bertani (LB) medium

- Tubes with  5 mL Luria-Bertani (LB) medium

- Tubes with  3 mL Terrific Broth (LB) medium

- Flasks with  20 mL LB medium

- Flasks with  50 mL LB medium

- Flasks with  50 mL TB medium

- Flask with  100 mL LB medium

- Flasks with  1 L LB medium

-  500 mL Lysis Buffer B (LB-B) Step 5.5/8.4

-  500 mL Elution Buffer-BI (EB-BI) Step 10.5

-  1 L Desalting Buffer-B (DB-B) Step 11.2


-  500 mL Elution Buffer-BII (EB-BII) Step 12.5

-  500 mL Storage Buffer-B (SB-B) Step 13.1

-  HisTrap HP 5mL **Cytiva Catalog #17524801**

-  HiPrep 26/10 Desalting Column **Cytiva Catalog #17508701**

-  HiTrap SP XL 5mL **Cytiva Catalog #17516101**

-  Quick Start™ Bradford 1x Dye Reagent **BioRad Sciences Catalog #5000205**

-



## Precision Plus Protein™ Unstained Protein Standards **Bio-rad** **Laboratories Catalog #1610363**

- Tricine-SDS-PAGE electrophoresis solutions (Step 6.2)
- 8% polyacrylamide gels for Tricine-SDS-PAGE (Step 6.2)

### **Equipments:**

- Thermomixer

Equipment	
<b>Thermomixer® R</b>	NAME
Dry block heating and cooling shaker, 120 V, 60 Hz, 1/cs	TYPE
Eppendorf	BRAND
T3317	SKU

- Orbital shaker

Equipment	
<b>MaxQ™ HP Incubated Tabletop Orbital Shaker</b>	NAME
MaxQ™ HP, 120 V 60 Hz, 6,5 A o 230 V 50/60 Hz, 3,2 A	TYPE
Thermo Scientific	BRAND
SHKE420HP	SKU
<a href="https://www.thermofisher.com/order/catalog/product/SHKE420HP">https://www.thermofisher.com/order/catalog/product/SHKE420HP</a>	LINK

- Floor model orbital shaker

Equipment	
<b>MaxQ™ HP Incubated and Refrigerated Console Shakers</b>	NAME
MaxQ™ 481 HP, 230 V, 50 Hz	TYPE
Thermo Scientific	BRAND
SHKE481HP	SKU
<a href="https://www.thermofisher.com/order/catalog/product/SHKE481HP">https://www.thermofisher.com/order/catalog/product/SHKE481HP</a>	LINK

- Centrifuge

Equipment	
<b>Sorvall™ Legend™ XT/XF Centrifuge Series</b>	NAME
Thermo Scientific	BRAND
75004541	SKU
<a href="https://www.thermofisher.com/order/catalog/product/75004541">https://www.thermofisher.com/order/catalog/product/75004541</a>	LINK

- Ultrasonic Processor 130W

Equipment	
<b>Ultrasonic Processor</b>	NAME
130-Watt Ultrasonic Processor	TYPE
Cole-Parmer	BRAND
ML-04714-52	SKU
<a href="https://www.coleparmer.com/p/cole-parmer-130-watt-ultrasonic-processors/44347">https://www.coleparmer.com/p/cole-parmer-130-watt-ultrasonic-processors/44347</a>	LINK

- Ultrasonic Processor 750W

Equipment	
<b>750-Watt Ultrasonic Processor</b>	NAME
CPX750	TYPE
Cole-Parmer	BRAND
ML-04711-60	SKU
<a href="https://www.coleparmer.com/p/cole-parmer-500-and-750-watt-ultrasonic-processors/16401">https://www.coleparmer.com/p/cole-parmer-500-and-750-watt-ultrasonic-processors/16401</a>	LINK

-Nanodrop

Equipment	
<b>NanoDrop™ One UV-Vis Spectrophotometer</b>	NAME
spectrophotometer	TYPE
Thermo Scientific	BRAND
ND-ONE-W	SKU
<a href="https://www.thermofisher.com/order/catalog/product/ND-ONE-W">https://www.thermofisher.com/order/catalog/product/ND-ONE-W</a>	LINK
Sample Volume (Metric): Minimum 1µL; Spectral Bandwidth: ≤1.8 nm (FWHM at Hg 254 nm); System Requirements: Windows™ 8.1 and 10, 64 bit; Voltage: 12 V (DC); Wavelength Range: 190–850 nm	SPECIFICATIONS

- FPLC system

Equipment	
ÄKTA pure	NAME
Protein purification system	TYPE
Cytiva	BRAND
29046665	SKU
<a href="https://www.cytivalifesciences.com/en/us/support/products/akta-pure-150-l-29046665">https://www.cytivalifesciences.com/en/us/support/products/akta-pure-150-l-29046665</a>	LINK

- Spectrophotometer UV/Vis
- Incubator (37°C)
- Water bath (60°C)
- Ultra Low–Temperature Freezer (-80°C)
- Freezer -20°C
- Refrigerator (4°C)
- Analytical balance

**Other:**

- Ice bath
- Microcentrifuge tubes
- Sterile 0.45 µm membrane filter
- 150 mL Superloop (Cytiva)
- Dialysis membrane
- Ultrafiltration tube (Amicon Ultra-15)

Equipment	
Amicon Ultra-15	NAME
PLTK Ultracel-PL membrane, 15 ML - 30 kDa cutoff	TYPE
Millipore	BRAND
UFC903024	SKU
<a href="https://www.merckmillipore.com/MX/es/product/Amicon-Ultra-15-Centrifugal-Filter-Unit,MM_NF-UFC903024?ReferrerURL=https%3A%2F%2Fwww.bing.com%2F&amp;bd=1">https://www.merckmillipore.com/MX/es/product/Amicon-Ultra-15-Centrifugal-Filter-Unit,MM_NF-UFC903024?ReferrerURL=https%3A%2F%2Fwww.bing.com%2F&amp;bd=1</a>	LINK

- Image Lab 6.1 Software (Bio-Rad)







#### BEFORE START INSTRUCTIONS

Ensure to have all the necessary materials and reagents already cleaned, sterilized and filter (in case of the purification solutions).

## Preparation of RT expression cells

2d 12h 20m 53s



### 1 Transformation of chemically competent BL21 (DE3) cells with pKJE7 plasmid.



- 1.1 Add  1  $\mu$ L of plasmidic DNA consisting of the pKJE7 plasmid from the  Chaperone Plasmid Set **Takara Bio USA, Inc. Catalog #3340** to  50  $\mu$ L of competent cells  Chemically Competent E. coli One Shot™ BL21(DE3) **Invitrogen - Thermo Fisher Catalog #C600003** . Mix the cells gently and incubate  On ice for  00:20:00 .

20m

- 1.2 Transfer the cells to a heat block at  42 °C and incubate for  00:00:53 .

53s

**1.3** Transfer the cells immediately to an ice bath and incubate  On ice for  00:05:00 . 5m




**1.4** Add  250  $\mu$ L of SOC medium at room temperature to the transformed cells and incubate at  225 rpm, 37°C, 01:00:00 . 1h

#### Note

SOC medium composition


A	B
Tryptone	2%
Yeast extract	0.5%
NaCl	10 mM
KCl	2.5 mM
MgCl <sub>2</sub>	10 mM
MgSO <sub>4</sub>	10 mM
Glucose	20 mM

Adjust to pH 7 and sterilize by filtration.

**1.5** Plate  25  $\mu$ L of transformed cells culture onto LB agar with the corresponding selective agent. Incubate the plates  Overnight at  37 °C . 18h



#### Note

The pKJE7 plasmid requires  30  $\mu$ g/mL chloramphenicol as selective agent.



**1.6** Select a single colony of transformed cells and inoculate in  3 mL Luria-Bertani (LB) medium supplemented with the selective antibiotic. Incubate 18h



 Overnight at  180 rpm, 37°C .





1.7 Centrifugate the cell culture at  10000 x g, 4°C, 00:05:00 . Remove the supernatant and resuspend the pellet in  500 µL LB medium .

5m





1.8 Add  500 µL of 30% v/v glycerol , mix by pipetting up and down and store at  -80 °C .

5m

## 2 Preparation of chemically competent BL21 (DE3) cells harboring pKJE7 plasmid.

2.1 Take BL21(DE3) cells harboring pKJE7 plasmid from a frozen glycerol stock using a bacterial inoculating loop and inoculate  3 mL LB liquid medium with  30 µg/mL of chloramphenicol . Incubate  Overnight at  180 rpm, 37°C .

18h

2.2 Inoculate  1 mL overnight culture in  100 mL LB medium with  30 µg/mL chloramphenicol and incubate at  180 rpm, 37°C, 03:00:00 .





3h

### Note

Monitor the cell growth by measuring the optical density (OD) at 600 nm and remove the cells from incubation when the OD reaches 0.3 to 0.4.

2.3 Chill the cell culture  On ice for  00:10:00 and centrifugate the cells at  4000 x g, 4°C, 00:10:00 .

20m

2.4 Gently resuspend the cell pellet in  8 mL of TFBI solution pre-cooled and incubate  On ice for  00:45:00 . Centrifugate the cells at  4000 x g, 4°C, 00:05:00 .




50m

### Note

TFBI medium composition

A	B
Potassium acetate	30 mM
Rubidium chloride	100 mM
Calcium chloride	10 mM
Manganese chloride	50 mM
Glycerol	15% v/v

Adjust to pH 5.8 with 1M acetic acid and sterilize by filtration.

- 2.5 Gently resuspend the cell pellet in  2.5 mL of TFBII solution pre-cooled and incubate  On ice for  00:05:00 .


5m

### Note


TFBII medium composition

A	B
MOPS	10 mM
Rubidium chloride	10 mM
Calcium chloride	75 mM
Glycerol	15% v/v

Adjust to pH 6.5 with 1M sodium hydroxide and sterilize by filtration.

- 2.6 Prepare aliquots of  50 µL of competent cells using microcentrifuge tubes previously

30m

chilled on an ice bath. Place the aliquots on dry ice until frozen and store at  -80 °C .





### 3 Transformation of chemically competent BL21 (DE3)/pKJE7 cells with the pET-MM4-RT plasmid, the expression vector for the quadruple mutant (E286R/E302K/L435R/D524A, designated as MM4) of the reverse transcriptase (RT) from Moloney Murine Leukemia Virus (MMLV).

#### CITATION



Yasukawa K, Mizuno M, Konishi A, Inouye K (2010). Increase in thermal stability of Moloney murine leukaemia virus reverse transcriptase by site-directed mutagenesis.. Journal of biotechnology.

LINK

<https://doi.org/10.1016/j.jbiotec.2010.09.961>

- 3.1 Add  1 µL of plasmidic DNA of  100 ng/µL pET-MM4-RT expression vector to  50 µL of competent cells BL21 (DE3)/pKJE7. For the transformation procedure  go to step #1 .





#### Note

The pET-MM4-RT plasmid requires  100 µg/mL carbenicillin as selective agent and the pKJE7 plasmid requires  30 µg/mL chloramphenicol . Use LB medium supplemented with both antibiotics as selective media.





## Small-scale screening cultures

1d 14h 16m 45s

### 4 Preparation of bacterial cultures for RT expression.



- 4.1 Inoculate  5 µL glycerol stock of BL21(DE3)/pKJE7/pET-MM4-RT or BL21(DE3)/pKJE7 cells in  5 mL culture medium (LB or TB) supplemented with selection agent(s). Incubate  Overnight at  200 rpm, 37°C .

18h

- 4.2** For each treatment, inoculate  500 µL overnight culture in  50 mL culture medium with  100 µg/mL carbenicillin . Use LB or TB according to the medium use for the overnight culture. Incubate  200 rpm, 37°C, 03:00:00 .

#### Note

Monitor the cell growth by measuring the optical density **(OD) at 600 nm** and remove the cells from incubation when the OD reaches **0.6**.

- 4.3** Once the culture reaches an OD<sub>600</sub> of 0.6, incubate the cell cultures  On ice for  00:20:00 before adding the inducer (IPTG).

## 5 Small-scale RT expression under different induction conditions.

- 5.1** Induce the expression of the RT under different conditions. Each treatment should be evaluated in triplicate. For example:

Strain	[IPTG]	Temperature	Medium
BL21(DE3)/pKJE7*	0.5 mM	16°C	LB
BL21(DE3)/pKJE7*	0.5 mM	37°C	LB
BL21(DE3)/pKJE7*	0.5 mM	16°C	TB
BL21(DE3)/pKJE7/pET-MM4-RT	0 mM	16°C	LB
BL21(DE3)/pKJE7/pET-MM4-RT	0.1 mM	16°C	LB
BL21(DE3)/pKJE7/pET-MM4-RT	0.5 mM	16°C	LB
BL21(DE3)/pKJE7/pET-MM4-RT	1.0 mM	16°C	LB
BL21(DE3)/pKJE7/pET-MM4-RT	0 mM	37°C	LB
BL21(DE3)/pKJE7/pET-MM4-RT	0.5 mM	37°C	LB

Strain	[IPTG]	Temperature	Medium
BL21(DE3)/pKJE7/pET-MM4-RT	0.5 mM	16°C	TB

\* BL21(DE3)pKJE7 strain is used as negative expression control.

**5.2** Incubate at 16 °C or 37 °C according to each treatment at 180 rpm for 16:00:00 . 16h

**5.3** Centrifugate the cell cultures at 6000 x g, 4°C, 00:10:00 . Discard the supernatant, remove all the liquid and leave the cell pellet as dry as possible. 10m

**5.4** Weigh the centrifugation tube with the cell pellet (total weight). 10m

#### Note

Weigh the empty tube prior centrifugation and subtract it to the total weight to calculate the weight of the cell pellet and hence the biomass produced.




**5.5** Resuspend the cell pellet in 5 mL lysis buffer B (LB-B) (pre-cooled). 5m

### Note


#### Lysis buffer B composition (LB-B)

A	B
NaH <sub>2</sub> PO <sub>4</sub> /Na <sub>2</sub> HPO <sub>4</sub> pH7.8	50 mM
NaCl	300 mM
DTT	2.5 mM
Imidazole	10 mM
Glycerol	5% v/v
PMSF	3 mM


Prepare the buffer with Milli-Q water and adjust to pH 7.8. Store at 4°C.

- 5.6** Disrupt cells by ultrasonication at an amplitude of 40%. Apply five cycles of  00:00:15 on  3m 45s and  00:00:30 off.

### Note

Place the tubes  On ice while processing.

Equipment	
<b>Ultrasonic Processor</b>	NAME
130-Watt Ultrasonic Processor	TYPE
Cole-Parmer	BRAND
ML-04714-52	SKU
<a href="https://www.coleparmer.com/p/cole-parmer-130-watt-ultrasonic-processors/44347">https://www.coleparmer.com/p/cole-parmer-130-watt-ultrasonic-processors/44347</a>	LINK

**5.7** Centrifugate at  6000 x g, 4°C, 00:15:00 . Recover the supernatant (soluble protein fraction) and discard the pellet.

15m

## 6 Analysis of RT expression.

**6.1** Measure total protein concentration by measuring absorbance at 280 nm in a NanoDrop spectrophotometer.

3m

Equipment	
<b>NanoDrop™ One UV-Vis Spectrophotometer</b>	NAME
spectrophotometer	TYPE
Thermo Scientific	BRAND
ND-ONE-W	SKU
<a href="https://www.thermofisher.com/order/catalog/product/ND-ONE-W">https://www.thermofisher.com/order/catalog/product/ND-ONE-W</a>	LINK
Sample Volume (Metric): Minimum 1 µL; Spectral Bandwidth: ≤1.8 nm (FWHM at Hg 254 nm); System Requirements: Windows™ 8.1 and 10, 64 bit; Voltage: 12 V (DC); Wavelength Range: 190–850 nm	SPECIFICATIONS

- 6.2 Analyze all supernatant samples by Tricine-SDS-PAGE electrophoresis through a 8% polyacrylamide gel. Load  100 µg protein sample per well. 4h

#### CITATION

Hermann Schägger (2006). Tricine-SDS-PAGE. Nature Protocols.

LINK







[10.1038/nprot.2006.4](https://doi.org/10.1038/nprot.2006.4)






- 6.3 Select the best conditions for protein expression according to the results analysis (biomass, total protein production and electrophoretic profile).

## Large-scale production of RT

1d 14h 24m 20s

### 7 Expression of recombinant RT.

- 7.1 Inoculate  10 µL glycerol stock of **BL21(DE3)/pKJE7/pET-MM4-RT** expression cells in  20 mL LB medium supplemented with  100 µg/mL carbenicillin and  30 µg/mL chloramphenicol. Incubate  Overnight at  180 rpm, 37°C. 18h


- 7.2 Inoculate  10 mL overnight culture in  1 L LB medium with  100 µg/mL carbenicillin and  30 µg/mL chloramphenicol. Incubate  200-220 rpm, 37°C, 03:00:00. 3h

- 7.3 Place the inoculum  On ice for  00:30:00 and then add  0.5 millimolar (mM) IPTG for induction. 30m




### Note

Do not add any additional inducers. For expression of the chaperones contained in pKJE7 plasmid, the basal expression is enough to promote correct RT enzyme folding.

**7.4** Incubate at  180 rpm, 16°C, 16:00:00 for recombinant protein expression.

16h

## 8 Soluble protein fraction recovery

**8.1** Centrifugate at  6000 x g, 4°C, 00:12:00 to harvest cells. Discard the supernatant ensuring to remove all the liquid and leave the cell pellet as dry as possible.


12m


**8.2** Weigh the centrifugation tube with the cell pellet (total weight).

3m

### Note

Weigh the empty tube prior centrifugation and subtract it to the total weight to calculate the weight of the cell pellet and hence the biomass produced.

**8.3** Store the cell pellet at  -80 °C until use (just in case that the purification step is not performed immediately after expression).

**8.4** Resuspend the cell pellet in  50 mL lysis buffer B (LB-B) (pre-cooled). If necessary, defroze the cell pellet in an ice bath before adding the lysis buffer.

5m




### Note

Lysis buffer B composition (LB-B)

A	B
NaH <sub>2</sub> PO <sub>4</sub> /Na <sub>2</sub> HPO <sub>4</sub> pH7.8	50 mM
NaCl	300 mM
DTT	2.5 mM
Imidazole	10 mM
Glycerol	5% v/v
PMSF	3 mM


Prepare the buffer with Milli-Q water and adjust to pH 7.8. Store at 4°C.

8.5

Disrupt cells by ultrasonication with an ultrasonic processor at an amplitude of 40% applying pulses of  00:00:10 of ultrasonication and  00:00:10 of pause during  00:04:00 .

4m 20s

### Note

Place the sample  On ice and keep it cold while processing.

## Equipment

750-Watt Ultrasonic Processor

NAME

CPX750

TYPE

Cole-Parmer


BRAND

ML-04711-60

SKU

<https://www.coleparmer.com/p/cole-parmer-500-and-750-watt-ultrasonic-processors/16401>

LINK

- 8.6** Centrifugate at  11000 x g, 4°C, 00:30:00 . Recover the supernatant (soluble protein fraction) and discard the pellet.

30m


### Note

Place the supernatant in an ice bath or store at 4°C until use.

## Purification of recombinant RT by FPLC

### 9 Sample preparation.

#### Note

Keep all protein samples  On ice during the purification process to avoid protein degradation.

- 9.1** Filter the supernatant (soluble protein fraction) through a 0.45 µm membrane.

9.2 Load the soluble protein fraction onto a 150 mL Superloop (Cytiva). Store at 4 °C until use.

## 10 Immobilized metal affinity chromatography (Ni<sup>2+</sup>-IMAC).

10.1 Connect a  HisTrap HP 5mL **Cytiva Catalog #17524801** to a FPLC system.

Equipment	
ÄKTA pure	NAME
Protein purification system	TYPE
Cytiva	BRAND
29046665	SKU
<a href="https://www.cytivalifesciences.com/en/us/support/products/akta-pure-150-l-29046665">https://www.cytivalifesciences.com/en/us/support/products/akta-pure-150-l-29046665</a>	LINK

10.2 Equilibrate the column with 8 column volumes (CV) of lysis buffer B (LB-B) at a flow of 2.5 mL/min.

10.3 Connect the Superloop charged with the protein fraction and load the sample onto the column at a flow of 2.5 mL/min.

10.4 Wash the column with 10 CV of LB-B at a flow of 2.5 mL/min

**10.5** Wash the column with 10 CV of 2% elution buffer-BI (EB-BI) at a flow of 2.5 mL/min.

**Note**


Elution buffer-BI composition (EB-BI).


A	B
NaH <sub>2</sub> PO <sub>4</sub> /Na <sub>2</sub> HPO <sub>4</sub> pH7.8	50 mM
NaCl	300 mM
DTT	2.5 mM
Imidazole	500 mM
Glycerol	5% v/v
PMSF	1 mM

Prepare the buffer with Milli-Q water and adjust to pH 7.8. Store at 4°C.


**10.6** Elute the proteins by passing 5 CV of 100% EB-BI through the column using a flow of 2.5 mL/min.

**10.7** Immediately after elution add [M] 2 millimolar (mM) EDTA and [M] 2.5 millimolar (mM) DTT to the eluted fractions.

**10.8** Analyze all collected fractions by Tricine-SDS-PAGE electrophoresis through a 8% polyacrylamide gel. Load  10 µL protein sample per well.

**10.9** Pool all elution fractions carrying the recombinant RT protein. Store at  4 °C until use.

## 11 Desalting step.

**11.1** Connect a  HiPrep 26/10 Desalting Column **Cytiva Catalog #17508701** to the FPLC system.

**11.2** Wash the column with 2.5 CV of Mili-Q water. Then, equilibrate the column with 2 CV of desalting buffer-B (DB-B). For both steps use a flow of 10 mL/min.

### Note




Desalting buffer-B composition (DB-B).

A	B
HEPES pH 7.5	50 mM
NaCl	40 mM
EDTA	2 mM
DTT	5 mM
Glycerol	5% v/v
PMSF	1 mM

Prepare the buffer with Milli-Q water and adjust to pH 7.5. Store at 4°C.

**11.3** Load the sample onto the column at a flow of 5 mL/min.

**11.4** Wash the column with 2 CV of DB-B for protein elution at a flow of 10 mL/min.

- 11.5** Analyze all collected fractions by qualitative Bradford assay using the  Quick Start™ Bradford 1x Dye Reagent **BioRad Sciences Catalog #5000205**. Pool the fractions with higher protein concentration.
- 11.6** Load the pool of desalted fractions onto a 150 mL Superloop (Cytiva). Store at  4 °C until use.
- 12 Cation exchange chromatography (CEC).**
- 12.1** Connect a  HiTrap SP XL 5mL **Cytiva Catalog #17516101** to the FPLC system.
- 12.2** Equilibrate the column with 10 CV of DB-B at a flow of 2 mL/min.
- 12.3** Connect the Superloop charged with the protein fraction and load the sample onto the column at a flow of 2 mL/min.
- 12.4** Wash the column with 5 CV of DB-B at a flow of 2 mL/min.
- 12.5** Elute proteins by washing the column with a linear gradient of 10 CV of elution buffer-BII (EB-BII). Use a flow of 2 mL/min.


## Note

Elution buffer-BII composition (EB-BII).

A	B
HEPES pH 7.5	50 mM
NaCl	1 M
EDTA	2 mM
DTT	5 mM
Glycerol	5% v/v
PMSF	1 mM

Prepare the buffer with Milli-Q water and adjust to pH 7.5. Store at 4°C.

**12.6** Analyze all collected fractions by Tricine-SDS-PAGE electrophoresis through a 8% polyacrylamide gel. Load  10 µL protein sample per well.

**12.7** Pool all elution fractions carrying the recombinant RT protein. Store at  4 °C until use.

## 13 Purified RT enzyme concentration and formulation.

**13.1** Load the purified RT enzyme pool onto a dialysis membrane (pre-hydrated). Place the membrane into a beaker with precooled storage buffer-B (SB-B) at a ratio 1:50 (v/v).




## Note

Storage buffer-B composition (SB-B).

A	B
Tris-HCl pH 7.5	50 mM
NaCl	150 mM
EDTA	0.1 mM
DTT	1 mM
Glycerol	50% v/v

Prepare the buffer with Milli-Q water and adjust to pH 7.5. Store at 4°C.

**13.2** Dialyze  Overnight at  4 °C with slow agitation.

**13.3** Recover the dialized protein, load it onto an **Amicon Ultra-15ML - 30 kDa cutoff centrifugal filter**. Concentrate until a concentration equal or higher than  1 mg/mL .

## Equipment

### Amicon Ultra-15

PLTK Ultracel-PL membrane, 15 ML - 30 kDa cutoff

Millipore

UFC903024

[https://www.merckmillipore.com/MX/es/product/Amicon-Ultra-15-Centrifugal-Filter-Unit,MM\\_NF-UFC903024?ReferrerURL=https%3A%2F%2Fwww.bing.com%2F&bd=1](https://www.merckmillipore.com/MX/es/product/Amicon-Ultra-15-Centrifugal-Filter-Unit,MM_NF-UFC903024?ReferrerURL=https%3A%2F%2Fwww.bing.com%2F&bd=1)

NAME

TYPE

BRAND

SKU

LINK

## Note

Monitor protein concentration measuring absorbance at 280 nm using a NanoDrop spectrophotometer.

## Equipment

### NanoDrop™ One UV-Vis Spectrophotometer

spectrophotometer

Thermo Scientific

ND-ONE-W

<https://www.thermofisher.com/order/catalog/product/ND-ONE-W>

Sample Volume (Metric): Minimum 1 µL; Spectral Bandwidth: ≤1.8 nm (FWHM at Hg 254 nm); System Requirements: Windows™ 8.1 and 10, 64 bit; Voltage: 12 V (DC); Wavelength Range: 190–850 nm

NAME

TYPE

BRAND

SKU

LINK

SPECIFICATIONS

- 13.4** Prepare aliquots of 50  $\mu\text{L}$  of concentrated RT enzyme .
- 13.5** Add 0.05 % (v/v) tergitol NP-40 to the enzyme aliquots and store at  $-20\text{ }^{\circ}\text{C}$  .
- 13.6** Determine final protein concentration by measuring absorbance at 280 nm in a NanoDrop spectrophotometer.
- 13.7** Analyze the final RT enzyme formulation by Tricine-SDS-PAGE electrophoresis through a 8% polyacrylamide gel. Load 3  $\mu\text{L}$  protein sample per well. Load 3  $\mu\text{L}$  protein ladder Precision Plus Protein™ Unstained Protein Standards **Bio-rad Laboratories Catalog #1610363** .
- 13.8** Analyze the electrophoresis gel by densitometry using the **Image Lab 6.1 Software (Bio-Rad)**. Determine protein concentration for each RT enzyme aliquot analyzed using the protein ladder as weight standard.

#### Note

The protein ladder

Precision Plus Protein™ Unstained Protein Standards **Bio-rad Laboratories Catalog #1610363**

includes three reference bands: The 50 KDa with 750 ng, the 20 KDa and 100 KDa bands with 150 ng each per each 10  $\mu\text{L}$  of the protein ladder mix.