



3 ▼

Oct 07, 2021

Colorimetric determination of urea V.3

PLOS One ✓

noah.langenfeld¹, Laurenpayne¹, Bruce Bugbee¹¹Crop Physiology Laboratory, Utah State University

2

dx.doi.org/10.17504/protocols.io.byvipw4e

USU Crop Physiology Laboratory

PLOS ONE Lab Protocols

noah.langenfeld

This protocol measures the absorbance of urea in solution in complexation with diacetyl monoxime at 520 nm and is linearly proportional to concentration up to 5 mM urea.

DOI

dx.doi.org/10.17504/protocols.io.byvipw4e<https://journals.plos.org/plosone/article?id=10.1371/journal.pone.0259760>

noah.langenfeld , Laurenpayne , Bruce Bugbee 2021. Colorimetric determination of urea. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.byvipw4e>

noah.langenfeld



protocol

Langenfeld NJ, Payne LE, Bugbee B (2021) Colorimetric determination of urea using diacetyl monoxime with strong acids. PLOS ONE 16(11): e0259760.

<https://doi.org/10.1371/journal.pone.0259760>

Added purity of reagents.

urea, urea assay, colorimetric urea, urea test, thiosemicarbazide, diacetyl monoxime

protocol ,

Oct 07, 2021

Oct 07, 2021

53898

thiosemicarbazide, diacetyl monoxime, water, sulfuric acid, phosphoric acid, ferric chloride

This protocol utilizes strong acids. Ensure proper precautions and safety equipment are used when dealing with these chemicals.

Ensure a lab coat, lab goggles, and plastic gloves are worn throughout this assay.

Mixed Acid Reagent Preparation

1

Dissolve **2.5 mg** ferric chloride (97% purity) in **45 mL** deionized water in a 250 mL volumetric flask.

2 Add **80 µL** of **14.8 Molarity (M)** concentrated phosphoric acid (85 wt. % phosphoric acid in water).

3






Caution: Preparation of **4.7 Molarity (M)** sulfuric acid is highly exothermic. Slowly add acid to **100 mL** water and stir frequently to avoid spattering.

Prepare **4.7 Molarity (M)** sulfuric acid by diluting **65.25 mL** concentrated sulfuric acid (**18 Molarity (M)**, 98% purity) up to **250 mL** with deionized water.

4 Dilute ferric chloride, deionized water, and phosphoric acid mixture to **250 mL** with **18 Molarity (M)** sulfuric acid.






5 Mix until dissolved.

Mixed Color Reagent Preparation

- 6 Add  **20.9 mg** diacetyl monoxime (99% purity) and  **52.4 mg** thiosemicarbazide (99% purity) to a 250 mL volumetric flask.
- 7 Dilute to volume ( **250 mL**) with deionized water.
- 8 Mix until dissolved.

Urea Assay

10m 5s

- 9 Prepare known concentrations of urea at 0, 1, 2, 3, 4, and 5 mM for a calibration curve.
- 10 Fill a 600 mL beaker with  **200 mL** deionized water.
- 11 Place the 600 mL beaker on hot plate and bring to a boil.
- 12 To test a sample, aliquot  **1 mL** into a 20 mL glass test tube.
- 13 Add  **2 mL** Mixed Acid Reagent.
- 14 Add  **2 mL** Mixed Color Reagent.
- 15 Ensure contents of test tube are vigorously mixed by vortexing for at least  **00:00:05** . 5s
- 16 Repeat steps 12-15 for each sample and standard to be analyzed.

- 17 Ensure test tubes are labeled and place in boiling water bath for 🕒00:20:00 .
- 18 Remove tubes from water bath and let cool for 🕒00:10:00 . 10m
- 19 Fill plastic cuvette to line with cooled solution.
- 20 Place sample cuvette in spectrophotometer and record absorbance at 520 nm. Blank against the standard solution with both reagents, but no urea.
- 21 Use the calibration curve to determine urea concentrations of unknown solutions.