

6



Aug 14, 2020

• Protein isolation from the extract of coelomic fluid of the Lumbricus terrestris.

Angel A Justiz-Vaillant¹

¹University of the West Indies St. Augustine

1	Works for me	dx.doi.org/10.17504/protocols.io.bjsaknae					
University of the West Indies angel.vaillant@sta.uwi.edu							
Angel Justiz-Vaillant University of the West Indies St. Augustine							

DOI

dx.doi.org/10.17504/protocols.io.bjsaknae

PROTOCOL CITATION

Angel A Justiz-Vaillant 2020. Protein isolation from the extract of coelomic fluid of the Lumbricus terrestris.. **protocols.io**

https://dx.doi.org/10.17504/protocols.io.bjsaknae

1.	0	N I	0	г

This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Aug 14, 2020

LAST MODIFIED

Aug 14, 2020

PROTOCOL INTEGER ID

40482

- 1 Six earthworms, Lumbricus terrestris, after removal from the earth are placed in a 15 ml glass tube and washed thoroughly with PBS, pH 7.4.
- 2 The body of the annelid is pierced with an 18-gauge needle.
- 3 The coelomic fluid content of the annelid is added to 2 ml of PBS, pH 7.4.
- 4 After incubation at 4°C for 1h, 0.9 ml chloroform is added.
- 5 The mixture is centrifuged at $1000 \times g$ for 5 min at $4^{\circ}C$.

protocols.io
1
08/14/2020

Citation: Angel A Justiz-Vaillant (08/14/2020). Protein isolation from the extract of coelomic fluid of the Lumbricus terrestris.. https://dx.doi.org/10.17504/protocols.io.bjsaknae

6	The supernatant (1.1 ml) is collected and equal volume of cold ethanol is added drop wise to the preparation.
7	The mixture is then incubated at 4°C for 12 h.
8	Then, it is centrifuged at 4°C for 5 min.
9	The pellet (protein extract) is resuspended in 0.3 ml of PBS, pH 7.4.
10	The preparation is dialyzed against 1L of PBS, pH 7.4 for 12 h at 4° C.
11	The protein concentration is assessed by the Bradford method or any available method.