



MAR 01, 2024

🌐 LENTIVIRAL TRANSDUCTION OF HUMAN PLURIPOTENT STEM CELLS

Renuka Ravi Gupta^{1,2,3,4}, Nona Farbehi^{1,2,3,5}, henderson^{3,6}, Vikram Khurana^{3,7}, Gist Croft^{3,8}, Lorenz Studer^{3,6}, Joseph Powell^{1,2,3,4}

¹Garvan Institute of Medical Research, Sydney, NSW 2010, Australia;

²Garvan Weizmann Center for Cellular Genomics, Garvan Institute of Medical Research, Sydney, NSW 2010, Australia;

³Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD 20815, USA;

⁴School of Medical Science, University of New South Wales, Sydney, NSW, 2052, Australia;

⁵Graduate School of Biomedical Engineering, University of New South Wales, Sydney, NSW, 2052, Australia;

⁶The Centre for Stem Cell Biology, Developmental Biology Program, Sloan Kettering Institute for Cancer Research, New York, NY, USA;

⁷Ann Romney Center for Neurologic Diseases, Brigham and Women's Hospital and Harvard Medical School, Boston, MA 02115, USA;

⁸The New York Stem Cell Foundation Research Institute, New York, NY, USA

OPEN ACCESS



Protocol Citation: Renuka Ravi Gupta, Nona Farbehi, henderson, Vikram Khurana, Gist Croft, Lorenz Studer, Joseph Powell 2024. LENTIVIRAL TRANSDUCTION OF HUMAN PLURIPOTENT STEM CELLS. [protocols.io](https://protocols.io/view/lentiviral-transduction-of-human-pluripotent-stem-c9wtz7en) <https://protocols.io/view/lentiviral-transduction-of-human-pluripotent-stem-c9wtz7en>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development
We are still developing and optimizing this protocol

Created: Feb 29, 2024

Last Modified: Mar 01, 2024



Renuka Ravi Gupta

ABSTRACT

We have developed a protocol for lentiviral transduction of human pluripotent stem cells (hPSCs), including induced pluripotent stem cells (iPSCs) or human embryonic stem cells (hESCs). In this protocol, concentrated lentiviral supernatant with a Multiplicity Of Infection (MOI) ranging from 0.1 to 0.3 (equivalent to 10% to 30% Blue Fluorescent Protein (BFP) positive cells) is combined with E8 Flex media and added to adherent H9 CRISPRi dCAS9 cells in 48-well plates. Subsequent centrifugation (spinfection) is employed to ensure efficient transduction. Transduction efficiency is assessed by determining the percentage of cells expressing Blue Fluorescent Protein (BFP) using Fluorescence Activated Cell Sorting (FACS).

ATTACHMENTS

[LENTIVIRAL
TRANSDUCTION OF
HUMAN PLURIPOTENT
STEM CELLS_.docx](#)

Keywords: ASAPCRN

Funders Acknowledgement:

ASAP

Grant ID: ASAP-000472

| A | B | C |
|-------------------------------------|--------------------------|-----------|
| MATERIAL | COMPANY | CATALOG |
| 48 well TC treated plate | Falcon | 353078 |
| 15ml polypropylene centrifuge tubes | Falcon | 352096 |
| 5ml serological pipettes | Corning | 4487 |
| 10ml serological pipettes | Corning | 4488 |
| DNA Low-bind tubes 1.5ml | Eppendorf | 022431021 |
| P1000 tip | Neptune | BT1250 |
| FBS | Bovogen | 2008A |
| DPBS | Thermo Fisher Scientific | 14040133 |
| E8 Flex media Kit | Thermo Fisher Scientific | A2858501 |
| RevitaCell Supplement(100x) | Thermo Fisher Scientific | A2644501 |
| Accutase | StemCell Technologies | 7922 |
| Vitronectin-N (VTN-N) | Thermo Fisher Scientific | A14700 |

REAGENT COMPOSITION

| A | B |
|---------------------------|--------------|
| FACS Buffer (PBS +2% FBS) | |
| REAGENT | VOLUME IN mL |
| PBS | 49 |
| FBS | 1 |

Day 0: Coating wells with VTN-N and seeding hPSCs

- 1 Add 60 ul of VTN-N to 6 ml of DPBS.

- 2 Coat 100 ul per well in a 48 - well plate.
- 3 Incubate the plate at room temperature for an hour and the plate is ready to be used.
- 4 Seed 3×10^4 cells/cm² in a 48 well plate with E8 flex media and RevitaCell after dissociating the cells with accutase.
- 5 Incubate the cells overnight at 37°C with 5% CO₂ and 20.9% O₂.

Day 1: Transduction of hPSCs with Lentiviral CRISPRi library supernatant

- 6 The H9 dCAS9 CRISPRi cells were transduced with the pooled CRISPRi library to harvest genomic DNA (gDNA), after calculating the exact volume of viral supernatant needed to achieve the desired MOI (multiplicity of infection).
- 7 Prepare 15 ml tubes with E8 flex media and concentrated lentiviral supernatants for a MOI of 0.1-0.3.
- 8 Add 200 ul/well of the virus media cocktail with the cells from step i.
- 9 To increase the transduction efficiency, centrifuge the plate at 300g for 20 minutes at 25°C.

- 10 Incubate the cells at 37°C with 5% CO₂ and 20.9% O₂ for 16-18 hours.

Day 2: Replace media

- 11 Aspirate the viral supernatant media gently and immediately add maturation media.
- 12 Return the plate back to the incubator.

Day 3: FACs Sort

- 13 Aspirate the spent media.
- 14 Wash the cells 10 times with DPBS to remove the viral particles from the lentivirus transduced hPSCs.

Note

Be very gentle while doing the washes as the cells tend to lift off during the wash step

- 15 Add 100 ul accutase and incubate the cells for 10 mins in the incubator.
- 16 Note: Ideally the hPSCs should dissociate as single cells.

17 Neutralize the accutase with E8 flex media and collect the cells into 1.5ml eppendorf tubes.

Note

Use a P1000 tip to pipette the cells up and down to break them into single cell suspension.

18 Centrifuge the cells at 300 g for 4 minutes.

19 Aspirate the spent media gently without disturbing the pellet.

20 Resuspend the cells in 300 ul of FACs buffer.

21 Transfer the cells with the FACs buffer into FACs tubes.

22 Sort the H9 CRISPRi cells to obtain 10- 30% BFP positive cells and collect the cells in E8 flex media.

23 Centrifuge the cells at 300 g for 4 mins.

24 Aspirate the spent media to obtain a cell pellet to be frozen or freshly used for DNA extraction.

