




Oct 10, 2022

DIVA Design instructions SOP (User's section)

Henrique De Paoli¹, Christopher J Petzold¹¹Lawrence Berkeley National Laboratory

1 Works for me



Share

This protocol is published without a DOI.

LBNL (ABF)

Agile BioFoundry



Christopher J Petzold

Lawrence Berkeley National Laboratory

ABSTRACT

This is a "how to" step-by-step instruction for creating *Projects* and *Designs* in DIVA for automated design of oligos, plasmids and build instructions. User can either build plasmids on their own or submit the final Design for review/PI approval and third party (e.g., DIVA team) construction pipeline. DIVA also allows for tracking Design progress and is linked to ICE (a repository for parts, plasmids and strains used as Design input and/or construction output in DIVA).

EXTERNAL LINK

<http://diva.agilebiofoundry.org>

PROTOCOL CITATION

Henrique De Paoli, Christopher J Petzold 2022. DIVA Design instructions SOP (User's section). **protocols.io**

<https://protocols.io/view/diva-design-instructions-sop-user-39-s-section-byprpvm6>



FUNDERS ACKNOWLEDGEMENT

U. S. Department of Energy, Energy Efficiency and Renewable Energy, Bioenergy Technologies Office

Grant ID: DE-AC02-05CH11231

KEYWORDS

DIVA, Design, DNA construction

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Sep 30, 2021

LAST MODIFIED

Oct 10, 2022

PROTOCOL INTEGER ID

53713

BEFORE STARTING

1. Define backbone(s), insert(s), and assembly strategy.

a.Guided by Target/Host expert (e.g., expression of a gene cluster) and/or Learn (e.g., design of experiments), compile all ICE entries (i.e., DNA parts linked to Strains **with samples**) and/or annotated files of **sequence-verified** DNA/protein sequences (e.g., genbank) to be used.

i.Synthetic fragments will not require physical template pre-available.

ii.Synthetic fragments must be $\leq 5\text{kb}$ and verified for synthesis constraints (e.g., <https://boost.jgi.doe.gov/index.html>)

iii.Optional: Codon juggle, polish and/or partition synthetic fragments

2. Resource *Crop data* and *Relevant milestone* (contact PI for reference)

3. You will be asked to select a Project from a self-complete list. If you don't have a Project, create one first. Although Projects can be added later to each Design from the Design page, you must have a Project linked to your Design to submit it for review.

DIVA access

- 1 Go to DIVA server at <https://diva.agilebiofoundry.org/> and login by using your LBNL LDAP/affiliate username and password.

You may need access to the ABF ICE parts registry. To create an account/login go to:
<https://registry.agilebiofoundry.org/>

Define Project

- 2 Click the New Project button.




Fig.1

You can skip to step 3 if you already have a Project created that you want to use for your designs.

- 2.1 Input a *Project Name, Description* and define the DIVA Team who will review your design.

Create new project [X]

Name

Example Project 1

Description 426 characters left

Comprise all designs associated with the NOG pathway for the production of acetyl-CoA.

Principal Investigators +add

Henrique De Paoli
hcdepaoli@lbl.gov [X]

DIVA Team

System Default +

Save **Cancel**

Fig.2

The System Default is the LBNL ESE DIVA Team. If you are not using the ESE DIVA team for your designs, create a Group (with desired people) within your Profile (top right corner) so the group appears in the self-complete field.

- 2.2 Add your Principal Investigator/Supervisor's Name, who is responsible for approving your expenses, then save the project.

Create Designs

- 3 Click the New Design button.




Fig.3

- 4 Input *Design Name*, *Description*(with *Crop data*, e.g., 11-2-TRY_M1, and *Relevant milestone*, e.g., FY21Q4_DBTL_AS2) and *Project*. Hit Create (it will take you directly to *Device Editor*).




Fig.4

- 5 Input desired sequences in each cell using 'right click' and collecting from the appropriate source (e.g., Search for Part, Fig.5B; insert from ICE, Fig.5C-D, e.g., [ABF_009895](#)). Use left-most column for plasmid backbone.

First-time Users: follow links of Mock Parts to win a hands-on experience!!!
 - You will design 2 variants, 1 insert each, and Run j5 as SLIC/Gibson/CPEC.




Fig.5A.

Pro tip: Change the SynBol icon to help you with the design structure.




Fig.5D

- 6 Define *Part Source*, *Name* (with '_marker' and '_enzyme', e.g., bb_tDNA_kan and ABF_009842_EcoRI_amp), and *Start BP* and *End BP*. If using gDNA, use '_gDNA' instead of '_marker'.




Fig.6

Pro tip: Part source parameters can also be defined by going under *Part* tab, change part definition and/or 'double clicking' on desired part.

- 7 Add columns and rows as needed for the remaining parts (e.g., Insert 1 [ABF_009896](#), Insert 2 [ABF_009569](#)).




Fig.7

- 8 Under the Part tab, define Forced Assembly Strategy for each part individually (if left blank, PCR will be automatically adopted by j5). Avoid PCRs larger than 5.1 kb.
- Optional: define Eugene rules.




Fig.8A




Fig.8B

- 9 Under the Collection tab, define Columns names (double click on e.g., 'No_name_0'). Also, define combinatorial (set by default) vs 'one per row' (must have all rows filled to be toggled).
 - a. Optional: select cells and insert empty part as needed
 - b. Optional: edit DNA assembly directives and parameters
 - c. Optional: define Column Direction (must toggle in the Device Editor area, not in the table).




Fig.9

Run j5

- 10 Press the J5 (green) button to open the 'j5 controls' and select Preview assembly (this is simply putting sequences together). Inspect output visually for desired output (e.g., number of variants) and adjust as needed.




Fig.10

- 11 While in the pop-up section, select the Advanced j5 tab for more options.
 - a. Optional: Edit j5 parameters
 - b. Define master lists




Fig.11
Also serves Step 12.

- 12 Select desired assembly from the drop-down menu.

You must select SLIC/Gibson/CPEC or Golden Gate (i.e., non-mock assembly; checks for assembly feasibility, designs oligos and output warnings and Build instructions) to be able to submit design for review.

Select the 'Run j5' button.

Note: If j5 returns an error message, check the j5 User's manual from Help tab in Device Editor.

- 13 Download results and open the combinatorial file (e.g., *pmas00001_combinatorial.csv*, or simply *pmas00001.csv* if making just one construct). Must be free of critical warnings (e.g., mispriming, incompatibilities between assembly pieces, homologous sequence repeat, etc).

Examples of “non-critical” warnings:

- *Primer3 could not find a reasonable set of primers.*
- The master plasmids/oligos list file is empty.




Fig.13A

Downloads > j5_20210930162141YOTF

Search j5_20210930162141YOTF

Name	Type	Compressed size
eugeneruleslist.eug	EUG File	0 KB
j5_parameters	Microsoft Excel Comma Separated Values File	1 KB
masterdirectsyntheselist	Microsoft Excel Comma Separated Values File	1 KB
masteroligolist	Microsoft Excel Comma Separated Values File	1 KB
masterplasmidlist	Microsoft Excel Comma Separated Values File	3 KB
mastersequences	Microsoft Excel Comma Separated Values File	1 KB
masterzippedsequences	Compressed (zipped) Folder	14 KB
partlist	Microsoft Excel Comma Separated Values File	1 KB
pmas00001	GenBank DNA	4 KB
pmas00001	Microsoft Excel Comma Separated Values File	5 KB
pmas00001_combinatorial	Microsoft Excel Comma Separated Values File	6 KB
pmas00002	Microsoft Excel Comma Separated Values File	5 KB
pmas00002	GenBank DNA	4 KB
targetpartsorder	Microsoft Excel Comma Separated Values File	1 KB

Fig.13B

13	Non-degenerate Part IDs and Sources																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																			</
----	-------------------------------------	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	----

Fig.13C

- 14 If no changes are needed, you can close the Advanced j5 Controls pop up and select the link Design Details (top right corner) to go to the main Design page.





Fig.14

Submit for DIVA Review

- 15 You can now Submit for DIVA Review. Pay attention to prompted messages and (un)check boxes accordingly.






Fig.15

- 16 Wait for Reviewer feedback. Address Reviewer requests (if any) and re-submit for Review, or solicit PI budgetary approval (Submit to PI button).
- PI was defined while creating the Project you linked to the Design.

Next Steps: PI will receive an email to open Design in DIVA for reviewing expenses, input Project Charging code and activity ID, and submit design for DNA Build construction.