

JAN 30, 2024

🌐 Chemogenetic modulation of catecholaminergic neurons

Cristian González-Cabrera¹, Csilla Novák¹, Maria P. Contreras¹, Ernesto Durán¹, Matthias Prigge¹

¹Neuromodulatory Network Group, Leibniz Institute for Neurobiology, Magdeburg

ASAP Collaborative Research Network

TeamPrigge



priggelab

OPEN ACCESS



DOI:

dx.doi.org/10.17504/protocols.io.rm7vzxdz5gx1/v1

Protocol Citation: Cristian González-Cabrera, Csilla Novák, Maria P. Contreras, Ernesto Durán, Matthias Prigge 2024. Chemogenetic modulation of catecholaminergic neurons.

protocols.io

<https://dx.doi.org/10.17504/protocols.io.rm7vzxdz5gx1/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

ABSTRACT

This protocol outlines a method for modulating noradrenergic and dopaminergic neuronal activity in the LC and SNc-VTA, respectively, using Cre-dependent hM3D and hM4D expression. It details stereotactic injection techniques for these chemogenetic tools and the preparation of the J60 agonist for intraperitoneal delivery. Developed for longitudinal studies spanning up to four weeks, the protocol aims to assess neuromelanin's activity-dependent accumulation through daily agonist administration.

MATERIALS

JHU37160 aka J60

Created: Jan 30, 2024

Last Modified: Jan 30, 2024

PROTOCOL integer ID: 94402

Keywords: Chemogenetics,
hM4Di, hM3D, DREADDs, long-
term modulation

Funders Acknowledgement:

ASAP
Grant ID: 020505

Stereotactic injection

1 Stereotactic injection of viruses into DbH-Cre (JAX# 033953) or DAT-Cre (JAX# 006660) animals

Viral Construct:

pAAV-hSyn-DIO-hM4D(Gi)-mCherry ([Addgene link](#))

Volume: 400ul

pAAV-hSyn-DIO-hM3D(Gq)-mCherry ([Addgene link](#))

Volume: 400ul

(in our experiment we mix our DIO-hTyr Virus 1:1 with DREADDs)

Stereotactic coordinates:

- Substantia nigra

ML: 1.4mm AP: 3.25mm DV: 4.0mm (blunt needle / NF34BL-2)

- Locus coeruleus

ML: 1.2mm AP: 5.2mm DV: 3.5mm (blunt needle / WPI NF34BL-2)

ML: 1.5mm AP: 5.2mm DV: 4.2mm (beveled needle / WPI NF34BV-2)

Injection set-up:

- Hamilton syringe coupled to a WPI injector (UMP3T-1) ([link](#))
- Injection speed 100nl/min


Expression

- 2 Assess accepted pain management protocol for injected animals, Start experiment as early as 3 weeks post-injections

Preparation agonist solution

3 Preparation of agonist solution

Several agonists are available for chemogenetic modulation with advantages and disadvantages, i.e., clozapine N-oxid (CNO), clozapine, or JP60. We are using J60 (JHU37160)
Please look at the reference for more detail (Bonaventura et al., 2019)



JHU37160

chemical structure JHU37160 aka J60 (Bonaventura et al., 2019)

- Solve entire vial (1mg) in 10ml saline (0.9% NaCl)
(Compound can be ordered from Torcris ([link](#)))


The solubility is very good, so by vortexing for a couple of seconds it should be already dissolved.

- Aliquot in a volume of 0.3ml and freeze at -20degree

DO NOT USE IT FOR MORE THAN 1 MONTH!

Modulation of Neuronal Activity

- 4
 - Animals are injected at a dosage of 1mg/kg body weight. Therefore, each aliquot is the maximum for one animal. The injections are intraperitoneal, as shown in the picture. When administering repeated injections, use different locations in the belly to avoid the same injection site.



Intraperitoneal injection of mice. Please consult a local vet for more details.