


Expression and purification protocol of PI3KC3-C1 complex

Jun 01, 2021

Expression and purification protocol of PI3KC3-C1 complex

Chunmei Chang¹¹Laboratory of James H. Hurley, UC Berkeley - CA

1 Works for me



Share

dx.doi.org/10.17504/protocols.io.bseenbbe

Dorotea Fracchiolla

Sascha Martens lab, University of Vienna, Max Perutz Labs - ...

ABSTRACT

This protocol outlines the procedure for expression and purification of the PI3KC3-C1 complex.

ATTACHMENTS

[_ASAP_Team_Hurley_Hurley_lab_Chunmei_Chang_Exp_purif_protocol_PI3KC3.pdf](#)

DOI

dx.doi.org/10.17504/protocols.io.bseenbbe

PROTOCOL CITATION

Chunmei Chang 2021. Expression and purification protocol of PI3KC3-C1 complex. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bseenbbe>

KEYWORDS

expression, purification, PI3KC3-C1 complex, plasmid amplification, transfection, HEK GnTI cells

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Feb 14, 2021

LAST MODIFIED

Jun 01, 2021

OWNERSHIP HISTORY

Feb 14, 2021



Emily Hasser

University of Washington

Mar 01, 2021



Dorotea Fracchiolla

Sascha Martens lab, University of Vienna, Max Perutz Labs - Vienna (AT)

PROTOCOL INTEGER ID

47270

MATERIALS TEXT

Transfection in HEK GnTI cells

Materials:

- polyethylenimine (PEI, Polysciences), make 1 mg/ml stock
- hybridoma media (Gibco)
- HEK293 expression freestyle media (Life Science Technology)
- all the DNA constructs

Protein purification

Materials:

- Glutathione Sepharose 4B (GE Healthcare)
- Strep-Tactin Sepharose (IBA)
- reduced glutathione (GSH, GE Healthcare)
- desthiobiotin (Sigma)
- protease inhibitor cocktail (Roche)
- TEV protease (1 mg/ml, homemade)

Buffers:

- lysis buffer: 50 mM HEPES pH 7.4, 1% Triton X-100, 200 mM NaCl, 1 mM MgCl₂, 10% glycerol, and 1mM TCEP
- wash buffer: 50 mM HEPES pH 8.0, 200 mM NaCl, 1 mM MgCl₂ and 1 mM TCEP
- GSH elution buffer: 50 mM HEPES pH 8.0, 200 mM NaCl, 50 mM GSH, 1 mM MgCl₂ and 1 mM TCEP
- Strep elution buffer: 50 mM HEPES pH 8.0, 200 mM NaCl, 10 mM desthiobiotin, 1 mM MgCl₂ and 1 mM TCEP
- Gel filtration buffer: 20 mM HEPES pH 8.0, 200 mM NaCl, 1 mM MgCl₂ and 1 mM TCEP

SAFETY WARNINGS

Please refer to the Safety Data Sheets (SDS) for health and environmental hazards.

Plasmid amplification

1 Amplify the plasmid.

constructs: pCAG-VPS15, pCAG-TSF-VPS34, pCAG-TSF-BECN1, pCAG-GST-TEVs- A TG14

All of the constructs are Amp resistant and high copy, could be amplified in normal XL10 strain.

Transfection in HEK GnTI cells

2 Grow HEK293 GnTI suspension cells to 1.5-2 million cells/ml in the desired volume, **1 L cells** would be good for PI3KC3-C1 complex.

3 In warm hybridoma media, make 20 X solution of DNA (**20 µg/ml**) and PEI (**80 µg/ml**).

i.e. For a 1 L transfection add **4 mL PEI** to **50 mL Hybridoma media** and **1 mg DNA (in total)** in another **50 mL hybridoma media** . The ratio of the DNA is: 1.5 ATG14: 1.0 VPS34: 1.0 VPS15: 1.0 BECN1.



30m

Add PEI to the DNA, mix together and incubate for a further **00:30:00**.

5 Add the DNA and PEI mixture to cells.

Harvest cells 2d 12h

6 Harvest cells after expression for **60:00:00** - **72:00:00**.

3d



Spin down cells at **2000 rpm, 4°C, 00:20:00**.



Wash the pellet with 1 X PBS once, freeze and store in **-80 °C**.

Protein purification 8h



Lyse cells with the lysis buffer (**50 mL lysis buffer** for **1 L cell pellet**) supplemented with protease inhibitors. Clarify the lysate by centrifugation (**16000 rpm, 4°C, 01:00:00**).



4h

Incubate the supernatant with **2 mL Glutathione Sepharose 4B** at **4 °C** (or **Room temperature**) for **04:00:00**.

11 

Apply the mixture to a gravity flow column and wash the resin extensively with wash buffer.

12 Elute the protein complex with  **10 mL GSH elution buffer** .

13 

4h


Add  **1 mg TEV protease** to eluted proteins and incubate at  **4 °C**  **Overnight** .

14 

Apply the TEV-treated complexes to a 1 ml Strep-Tactin Sepharose gravity flow column, wash the resin extensively with wash buffer.

15 Elute the protein complex with  **5 mL Strep elution buffer** .

 **1 mL** per fraction.

16 Apply the eluted fraction with highest concentration (or gently concentrate all the fractions to  **1 mL**) to Superdex 6 column (10/300 Increase).

17 Collect the peak fractions and use freshly.