



Jun 10, 2025

Amplification of cytb gene of Plasmodium spp. in blood spots from wild animals

DOI

dx.doi.org/10.17504/protocols.io.bp2l619k5vqe/v1

Gabriela Ulloa Urizar¹

¹Universidade Federal Rural da Amazônia (UFRA), Belém, Pará, Brasil



Gabriela Ulloa Urizar

Universidade Federal Rural da Amazônia (UFRA), Belém, Pará, ...

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.bp2l619k5vqe/v1

Collection Citation: Gabriela Ulloa Urizar 2025. Amplification of cytb gene of Plasmodium spp. in blood spots from wild animals. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bp2l619k5vqe/v1>

License: This is an open access collection distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this collection and it's working

Created: February 14, 2022

Last Modified: June 10, 2025

Collection Integer ID: 58196



Funders Acknowledgements:

ERANet-LAC

Grant ID: ERANet17/HLH-0271

Instituto de Salud Carlos III

Grant ID: AC18/00054

CNPq

Grant ID: 140312/2020-0

CNPq

Grant ID: 201546/2020-5

CNPq

Grant ID: 400800/2019-5

CONCYTEC

Grant ID: N° 136-2018-FONDECYT

Abstract

Amplification of the Plasmodium mitochondrial cytochrome oxidase b gene (cytb) for detection of the Plasmodium genus in wildlife samples.

Guidelines

A) To detect *Plasmodium* spp. in humans and non human primates use the primers:

Procedure:


ROUND 1 PCR


1. Hydrate the primers with water without DNAases and dilute the primers for *Plasmodium* spp. (DW2 and DW4) to a concentration of 5uM.


Note


Final effective concentration in PCR is 0.5 μ M for 20uL volume.

2. Volumes for the pcr mix for 1 reaction (1 rx) were:

Primer DW2 [5uM]:  2 μ L

Primer DW4 [5uM]:  2 μ L

Master Mix (2X):  10 μ L

Free water:  1 μ L

Sample DNA:  5 μ L

3. PCR conditions:

94°C ---  00:03:00

94°C ---  00:00:20

62°C ---  00:00:30


72°C ---  00:01:00

The last 3 steps are repeated for 40 cycles.

72°C ---  00:10:00

4°C --- store at this temperature until electrophoresis

Electrophoresis

4. Add 15 ul of PCR product to each well of the agarose gel. Then run on a 1.5% agarose gel prepared with TAE buffer at 95 volts for  00:40:00



Note

Before nested PCRs are performed, dilute the PCR products for Plasmodium spp. (DW2 and DW4) by 20-fold, i.e. 2ul of PCR product in 38ul of DNAase-free water. Then use this dilution as template DNA for the following nested PCR.


B) nested PCR cytb in humans and non human primates use the primers:**ROUND 2 PCR**


1. Hydrate the primers with water without DNAases and dilute the primers for *Plasmodium* spp. (FP3 and RP3) to a concentration of 5uM.


Note


Final effective concentration in PCR is 0.5 μ M for 20uL volume.

2. Volumes for the pcr mix for 1 reaction (1 rx) were:

FP3 primer [5uM]:  2 μ L

RP3 Primer [5uM]:  2 μ L

Master Mix (2X):  10 μ L

Free water:  4 μ L

Diluted PCR product (1/20):  2 μ L

3. PCR conditions:

94°C ---  00:03:00

94°C ---  00:00:20

52°C ---  00:00:30

72°C ---  00:01:00


The last 3 steps are repeated for 40 cycles



72°C ---  00:10:00

4°C --- store at this temperature until electrophoresis

Electrophoresis

4. Add 15 ul of PCR product to each well of the agarose gel. Then run on a 2% agarose gel prepared with TAE buffer at 95 volts for  00:45:00 . Visualize bands under UV or blue light. Expected size: ~776 bp.

Materials

- DW2 and DW4 primers
- Platinum II Hot Start PCR Master Mix (2X)
- Nuclease-free water
- Template DNA (from blood spot samples)
- PCR tubes
- Thermal cycler
- Agarose, TAE buffer
- Gel electrophoresis system
- DNA stain
- DNA ladder (100 bp)



Before start

Materials

- DW2 and DW4 primers
- Platinum II Hot Start PCR Master Mix (2X)
- Nuclease-free water
- Template DNA (from blood spot samples)
- PCR tubes
- Thermal cycler
- Agarose, TAE buffer
- Gel electrophoresis system
- DNA stain
- DNA ladder (100 bp)

To detect *Plasmodium* spp. in humans and non human primates use the primers:

ROUND 1

DW2: 5'-TAA TGC CTA GAC GTA TTC CTG ATT ATC CAG-3'

DW4: 5'-TGT TTG CTT GGG AGC TGT AAT CAT AAT GTG-3'

Nested PCR cytb in humans and non human primates use the primers:

ROUND 2

FP3: 5'-TAT ATA ACT TAT TTT TTG ATA TG-3'

RP3: 5'-GTR ATW GCA TTA TCT GGA TGT GA-3'

For the PCR Mix (conventional and nested) we used the **Platinum™ II Hot-Start PCR Master Mix (2X)** from Invitrogen.



Files

 SEARCH

Protocol references

MARTINSEN, E. S., PAPERNA, I. & SCHALL, J. J. Morphological versus molecular identification of avian Haemosporidia: an exploration of three species concepts. *Parasitology* **133**, 279–288 (2006).

Acknowledgements

This protocol has been cited in our peer-reviewed publication:

<https://doi.org/10.1016/j.meegid.2024.105554>

It has now been updated to include step-by-step instructions, in compliance with protocols.io formatting requirements.