




Jul 23, 2020

Lentiviral_vector_production_with_PEI

Kenneth Matreyek¹, Anna Bruchez¹¹Case Western Reserve University School of Medicine**1** Works for me This protocol is published without a DOI.

MatreyekLab

 Kenneth Matreyek

ATTACHMENTS

[PEI_Transfections_
lentivirus protocols.pdf](#)

PROTOCOL CITATION

Kenneth Matreyek, Anna Bruchez 2020. Lentiviral_vector_production_with_PEI. **protocols.io**
<https://protocols.io/view/lentiviral-vector-production-with-pei-bhmyj47w>

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jun 17, 2020

LAST MODIFIED

Jul 23, 2020

PROTOCOL INTEGER ID

38296

MATERIALS

NAME	CATALOG #	VENDOR
PEI MAX - Transfection Grade Linear Polyethylenimine Hydrochloride (MW 40000)	NC1038561	Fisher Scientific

DISCLAIMER:

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](#) is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with [protocols.io](#), can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Preparation

1d

- 1 This protocol will require transfecting HEK 293T cells. Make sure you have some 293T cells thawed and growing for at least a couple of days before doing this protocol. Furthermore, prepare stocks of PEI and Diluent if you have not done so already.

1.1 PEI solution: 1 µg/µl PEI in H2O. I have used "Polysciences, Inc. Supplier Diversity Partner PEI

MAX - Transfection Grade Linear Polyethylenimine Hydrochloride (MW 40,000)"

1. Dissolve PEI in endotoxin-free dH₂O that has been heated to approximately **80 °C**.
2. Let cool to room temperature.
3. Adjust pH until it reaches **pH7**, filter sterilize (0.22 µm filter), aliquot, and store at **-20 °C**. A working stock can be kept at **4 °C**.

1.2 Diluent solution: **10 Milimolar (mM) HEPES**, **150 Milimolar (mM) NaCl** **pH7.05**

For **500 mL**:

1. Get **445 mL** distilled water
2. Add **5 mL HEPES**
3. Add **15 mL 5 Molarity (M) NaCl**
4. Filter sterilize (0.22 µm filter). Surround in foil to protect from light. Store at room temperature.

Transfection 20m

- 2 Prepare the mixture in a sterile microcentrifuge tube.

2.1 If transfecting a 10cm plate (eg. when trying to make larger stocks of lentiviral vectors):

1. Aim for a total volume of Diluent + DNA of **500 µl**
2. Add **10.5 µg total DNA**. When using VSV-G, I do a 5:5:1 ratio packaging plasmid, transfer vector, and pMD-VSVG.
3. Add **42 µl PEI** to the diluted DNA mixture. Immediately pipet up and down and / or vortex.
4. Incubate between **00:10:00 (minutes)** and **00:15:00 (minutes)**. Don't go over **00:20:00 (minutes)**.
5. Trypsinize, count, and plate 10 million cells per 10cm plate, in a total volume of **10 mL**.
6. After the incubation is over, add the DNA + PEI mixture dropwise to the cells, and return to the incubator.

2.2 If transfecting an individual 6-well (eg. when trying to test out some lentiviral vectors):

1. Aim for a total volume of Diluent + DNA of **80 µl**
2. Add **1.8 µg total DNA**. When using VSV-G, I do a 5:5:1 ratio packaging plasmid, transfer vector, and pMD-VSVG.
3. Add **7 µl PEI** to the diluted DNA mixture. Immediately pipet up and down and / or vortex.
4. Incubate between **00:10:00 (minutes)** and **00:15:00 (minutes)**. Don't go over **00:20:00 (minutes)**.
5. Trypsinize, count, and plate 1.5 million cells per 6-well, in a total volume of **2 mL**.
6. After the incubation is over, add the DNA + PEI mixture dropwise to the cells, and return to the incubator.

Change media 1d

- 3 The next day (between 12 and 24 hours), replace with fresh media. To increase lentivector concentration of the supernatants, use **6 mL** for 10cm plates and **1.5 mL** for 6-wells.

Collect the supernatants

3d

- 4 Collect media over the next four days and store pooled supernatant at 4°C . If using a viral envelope protein that is able to enter HEK 293T cells (eg. VSV-G or Ebola virus GP), it may be wise to collect multiple times in one day (Once when entering the lab, and one when going home), as the viral particles will be able to re-enter the producer cells essentially becoming a dead-end product. This may not be as important for viral envelopes with receptors that are not normally expressed in 293T cells (eg. SARS CoV), although its possible that virus stored at 4°C will still be more stable than virus left at 37°C .

Prepare the supernatants for use

30m

- 5 The supernatants may still have detached 293T cells floating which may "contaminate" downstream cultures.

- 5.1
 1. Centrifuge the supernatant at greater than $300 \times g$, 4°C 00:03:00.
 2. Filter the supernatant through a 0.45um filter.
 3. Aliquot and store the lentivector supernatant at -80°C .