

Jan 12, 2021

# S.O.F protocol for nuclei isolation from fresh and frozen tissues using OptiPrep® discontinuous gradient

Luciano G Martelotto<sup>1</sup><sup>1</sup>Harvard Medical School**1** Works for me [dx.doi.org/10.17504/protocols.io.brdbcm22w](https://dx.doi.org/10.17504/protocols.io.brdbcm22w)

Human Cell Atlas Method Development Community

Single Cell Core, Harvard Medical School

Luciano Martelotto  
Harvard Medical School

## ABSTRACT

This protocol is the result of the combination of various nuclei isolation protocols for single cell RNA-seq experiments using droplet-based methods, and is an extension of the Frankenstein (S.O.F means Son Of Frankenstein).

Developed to prepare nuclei isolates from fresh and frozen material of small-to-large sizes. The good thing is that it does not use FACS but OptiPrep® discontinuous gradient to remove debris. It is the alternative protocol when FACS is not available.

## EXTERNAL LINK

<https://research.unimelb.edu.au/centre-for-cancer-research/our-research/single-cell-innovation-lab>

## ATTACHMENTS

\_\_\_Frankenstein\_\_\_proto  
col\_for\_nuclei\_isolation\_fro  
m\_fresh\_and\_frozen\_tissu  
e\_FINAL.pdf

## DOI

[dx.doi.org/10.17504/protocols.io.brdbcm22w](https://dx.doi.org/10.17504/protocols.io.brdbcm22w)

## EXTERNAL LINK

<https://research.unimelb.edu.au/centre-for-cancer-research/our-research/single-cell-innovation-lab>

## PROTOCOL CITATION

Luciano G Martelotto 2021. S.O.F protocol for nuclei isolation from fresh and frozen tissues using OptiPrep® discontinuous gradient . **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.brdbcm22w>



## EXTERNAL LINK

<https://research.unimelb.edu.au/centre-for-cancer-research/our-research/single-cell-innovation-lab>

## KEYWORDS

snRNAseq, 10x, nuclei, isolation, FACS

## LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Jan 12, 2021

LAST MODIFIED

Jan 12, 2021

PROTOCOL INTEGER ID

46212

ATTACHMENTS

[\\_\\_\\_Frankenstein\\_\\_\\_protocol\\_for\\_nuclei\\_isolation\\_from\\_fresh\\_and\\_frozen\\_tissue\\_FINAL.pdf](#)

GUIDELINES

Use RNA techniques.

The use of RNase Inhibitor is HIGHLY recommended.

## Required Buffers and Reagents

 Nuclei EZ lysis buffer

1. **Sigma Catalog #EZ PREP NUC-101** (keep 4°C)
2. Nuclei wash and resuspension buffer (prepare fresh, keep 4°C)

1x PBS

1.0% BSA (MACS® BSA Stock Solution, Miltenyi)

0.2-0.5 U/uL RNase Inhibitor (Protector RNA Inhibitor, Millipore Sigma)

### 3. Gradient Solutions

G1: OptiPrep®: 60% (w/v) solution of iodixanol in water.

G2: 150 mM KCl, 30 mM MgCl<sub>2</sub>, 120 mM, Tricine-KOH pH 7.8

G3: solution containing 50% (w/v) of iodixanol --> mix 5 volumes of G1 with 1 volume of solution G2 (final: 50% iodixanol, 25 mM KCl, 5 mM MgCl<sub>2</sub>, 20 mM Tricine-KOH pH 7.8).

G4: 25 mM KCl, 5 mM MgCl<sub>2</sub>, 20 mM Tricine-KOH pH 7.8.

G5: solution containing 25% (w/v) of iodixanol --> mix equal volumes of G3 and G4 (final: 25% iodixanol, 25 mM KCl, 5 mM MgCl<sub>2</sub>, 20 mM Tricine-KOH pH 7.8).

G6: solutions of 30% --> mix 6 volumes of G3 with 4 volumes of G4 (final: 30% iodixanol, 25 mM KCl, 5 mM MgCl<sub>2</sub>, 20 mM Tricine-KOH pH 7.8).

G7: solutions of 35% --> mix 7 volumes of G3 with 3 volumes of G4 (final: 35% iodixanol, 25 mM KCl, 5 mM MgCl<sub>2</sub>, 20 mM Tricine-KOH pH 7.8).

### 4. Stock solutions and G2/G4 Set-Up

Keep the following Stock Solutions at 4°C:

-500 mM Tricine: 8.96 g in 100 ml water

-1 M KCl: 7.45 g in 100 ml water

-1 M MgCl<sub>2</sub>·6H<sub>2</sub>O: 20.3 g in 100 ml water

G2: To 50 ml water add 24 ml, 15 ml and 3 ml respectively of the Tricine, KCl and MgCl<sub>2</sub>·6H<sub>2</sub>O stock solutions (above); adjust to pH 7.8 with 1 M KOH and make up to 100 ml.

G4: Dissolve 8.5 g of sucrose in 50 ml of water; add 4 ml, 2.5 ml and 0.5 ml respectively of the Tricine, KCl and MgCl<sub>2</sub>·6H<sub>2</sub>O stock solutions (above); adjust to pH 7.8 with 1 M KOH and make up to 100 ml.

### SAFETY WARNINGS

See SDS (Safety Data Sheet) for safety warnings and hazards.

### DISCLAIMER:

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](https://protocols.io) is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with [protocols.io](https://protocols.io), can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

#### ABSTRACT

This protocol is the result of the combination of various nuclei isolation protocols for single cell RNA-seq experiments using droplet-based methods, and is an extension of the Frankenstein (S.O.F means Son Of Frankenstein).

Developed to prepare nuclei isolates from fresh and frozen material of small-to-large sizes. The good thing is that it does not use FACS but OptiPrep® discontinuous gradient to remove debris. It is the alternative protocol when FACS is not available.

#### BEFORE STARTING


All samples and reagents are kept on ice or at 4 °C (wet ice).

Prepare all buffers and reagents as described in the "Materials" section.

### Tissue Homogenization



- 1 Mince/chop tissue with a razor blade to small pieces. The tissue may be as small as a grain of rice.

For mincing the tissue, you may take the tube out of ice, however, be quick and return to ice.

- 2 Add  **300 µl** of chilled Nuclei EZ Lysis Buffer (supplemented with RNase Inhibitor 0.2-0.5 U/µL) to the tissue in 1.5 mL tube.
- 3 **Gently** homogenize the sample using a douncer by stroking 10-20 times. Keep tube on ice at all times.

### Nuclei Isolation

10m

- 4 Add an extra  **900 µl** of chilled Nuclei EZ Lysis Buffer (supplemented with RNase Inhibitor 0.2-0.5 U/µL), mix gently by pipetting using **wide-bore tips** and incubate on ice for  **00:05:00**. Repeat mixing 2-3 times during the incubation. RNase inhibitor will now be in the range of 0.2-0.5 U/µL.
- 5 Filter homogenate using a 70 µm-strainer mesh to fit a 15 ml Falcon tube (e.g. pluriStrainer Mini 70 µm Cell Strainer).

pluriStrainer Mini 70 µm  
Cell Strainer


pluriSelect 43-10070-40






- 6 Transfer flow through into a 1.5 mL LoBind tube and centrifuge the nuclei at **500 x g** for **00:05:00** at **4 °C** and remove supernatant leaving behind ~ **50 µl**.
- 7 **Optional:** add **1 mL** of EZ Lysis buffer (supplemented with RNase Inhibitor 0.2-0.5 U/µL), gently resuspend <sup>10m</sup> pellet and incubate for **00:05:00** on ice. Then centrifuge the nuclei at **500 x g** for **00:05:00** at **4 °C**.
- 8 Remove supernatant as much as possible without disturbing pellet (if pellet looks loose leave ~ **50 µl** behind) and add **1 mL** of Nuclei Wash and Resuspension Buffer (supplemented with RNase Inhibitor 0.2-0.5 U/µL) **without resuspending**.
- 9 Centrifuge the nuclei at **500 x g** for **00:05:00** at **4 °C**.
- 10 **Optional:** Remove supernatant as much as possible without disturbing pellet (if pellet looks loose leave ~ **50 µl** <sup>5m</sup> behind) and add **500 µl** of Nuclei Wash and Resuspension Buffer (supplemented with RNase Inhibitor 0.2-0.5 U/µL) **without resuspending**. Centrifuge the nuclei at **500 x g** for **00:05:00** at **4 °C**.
- 11 Remove supernatant as much as possible without disturbing pellet (if pellet looks loose leave ~ **50 µl** behind) and resuspend the pellet in **1 mL** of G5 buffer and transfer onto a 10 mL LoBind centrifuge tube. Then add another mL of G5 (final volume 2 mL).

#### Debris Removal 10m

- 12 Carefully, underlayer the sample with **2 mL** of the 30% iodixanol and then **2 mL** of the 35% iodixanol, as shown in the picture below.



- 13 Centrifuge the nuclei at **5000 x g** for **00:20:00** at **4 °C** using swinging bucket centrifuge with break set to <sup>20m</sup> OFF. (Note: using **10000 x g** is also possible).

- 14 After centrifugation, a white-ish band of nuclei between the 30%-35% iodixanol interface should be visible. Aspirate the top layers down until the white nuclei band at the interphase of 29%-35%. Using 200 uL tip, collect the nuclei band and transfer to a fresh 1.5 mL LoBind tube.
- 15 Then top up to ~1.3 mL with Wash and Resuspension Buffer (supplemented with RNase Inhibitor 0.2-0.5 U/uL), mix<sup>5m</sup> well but gently, and centrifuge the nuclei at  **500 x g** for  **00:05:00** at  **4 °C**.
- 16 Remove supernatant and repeat once more for a total of 2 washes.
- 17 Check integrity and purity under microscope and count manually or using an automatic counter. For automatic counter I recommend LUNA-FL™ Dual Fluorescence Cell Counter and Acridine Orange/Propidium Iodide (AO/PI) Cell Viability Kit (F23001).
- 18 Proceed to your amazing snRNA-Seq experiment!