



Jul 23, 2020

# Extraction of Oxycodone from Rat Brain for Mass Spec Analysis

Sierra Simpson<sup>1</sup>, Olivier George<sup>1</sup><sup>1</sup>University of California, San Diego**1** Works for me This protocol is published without a DOI.

George Lab

Tech. support email: [olgeorge@ucsd.edu](mailto:olgeorge@ucsd.edu)

Sierra Simpson

University of California, San Diego

## ABSTRACT

### Extraction of Oxycodone from Rat Brain for Mass Spec Analysis

Adapted from Momper Lab protocol

## PROTOCOL CITATION

Sierra Simpson, Olivier George 2020. Extraction of Oxycodone from Rat Brain for Mass Spec Analysis.

**protocols.io**<https://protocols.io/view/extraction-of-oxycodone-from-rat-brain-for-mass-s-825hyg6>

## KEYWORDS

Mass spec, oxycodone

## LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Nov 06, 2019

## LAST MODIFIED

Jul 23, 2020

## PROTOCOL INTEGER ID

29501

## MATERIALS

NAME	CATALOG #	VENDOR
Acetonitrile	AC1400.SIZE.1L	Bio Basic Inc.

## MATERIALS TEXT

2ml polypropylene tube, Sarsted, REF# 72.693.005

BioSpec Products, 2.3mm diameter zirconia / silica beads, catalog # 11079125z

d<sub>6</sub>-Oxycodone Cayman chemicals # 19249 - Internal Standard

Oxycodone - Standard

Syringe Filter Unit - Millex - GV, catalog # SLGVR04NL

Mass Spec Vials- [CTV-1206 - 500ul PP SR/Crimp Vial, 1000/pk](#)Mass Spec Vial Lids- [CTC-1370B - Snap TOP Cap, PTFE/Sil w/slit, Blue, 100/pk](#)

## SAFETY WARNINGS

Use ACN in the hood

Make sure you know how to use the bead beater machine

## BEFORE STARTING

- Prepare a homogenizing tube (2ml polypropylene tube, Sarsted, REF# 72.693.005) for each brain.
- Add ~100mg homogenizing beads to each tube (BioSpec Products, 2.3mm diameter zirconia / silica beads, catalog # 11079125z)
- Prepare a fresh 35mm petri dish and fresh razor blade for each brain, and boxes of wet ice to place dishes once brains are homogenized to keep tissue cool when weighing.
- Place 50ml of Acetonitrile on wet ice (4C) for 30min before homogenization step to cool before process.

#### Brain Homogenization

- 1 Brains are snap-frozen upon harvest and are stored in -80C until ready for processing. Thaw out the tissue slowly on ice for 5-10 minutes, once thawed place the brain on a sterile dish and mince the entire brain using a fresh razor blade.
- 2 Add ~20- 40mg minced tissue to each homogenizing tube.
- 3 Add ~375ul 4C acetonitrile (ACN) to each homogenizing tube containing the tissue and beads. 0.053mg tissue/ $\mu$ l ACN (example: 20mg tissue / 375  $\mu$ l ACN)
- 4 Add internal standard d6-Oxycodone to a concentration of 0.4ng/ $\mu$ l ( Cayman chemicals # 19249)
- 5 Store these tubes on ice until ready to homogenize on a bead beater machine.
- 6 Put tubes in a bead beater machine to homogenize and break up the tissue(s).

#### Bead Beating

- 7 The bead beater machine is located in the SKAGGS building in the Palmer Taylor lab / Dorrestein Labs. Make sure to contact them before use to reserve time on the machine. Make sure the sample tubes are secured – this machine shakes the samples quite quickly and if not properly secured can cause severe damage to the machine and the lab.
- 8 Set the machine at 25,000 cycles per min.
- 9 Set the vibration time for 30 seconds
- 10 After 30 seconds of vibration wait for 30 seconds with no vibration (prevents heat buildup of the tissue from the vibrating beads)
- 11 Repeat the vibration step two more times (total of 3 cycles)
- 12 Transfer all the samples back on ice and allow to cool for 5 – 10 mins.

- 13 After 5 – 10 mins centrifuge the homogenizing tubes at 5000 rpm for 5 mins.
- 14 Transfer the clear supernatant to the pre-labeled 1.5ml epi tube and place on ice (~500µl):
- 15 Place the homogenizing tube back on ice until ready to add H<sub>2</sub>O to wash along the wall of the homogenizing tube to collect any remaining sample stuck along the wall of the tube.
- 16 Add ~125µl H<sub>2</sub>O (=0.160mg/µl H<sub>2</sub>O, example 20mg tissue/125µl H<sub>2</sub>O) to the homogenizing tube.
- 17 Vortex vigorously for 3-5 seconds.
- 18 Transfer the resulting H<sub>2</sub>O mix (~125µl) to the pre-labeled 1.5ml polypropylene bullet tube containing ~500µl of homogenized supernatant – at this point the clear solution will become cloudy and place on ice (~625µl)
- 19 Centrifuge samples at 5000 rpm for 5 mins.

#### Filtration and Prep for Mass Spec

- 20 Transfer the clear supernatant (~625µl) to a 1cc syringe with an attached filter unit (Millex – GV, catalog # SLGVR04NL) - this step further cleans the solution to filter out any additional contaminants that can clog the mass spec.
- 21 Discard the tube containing the pellet
- 22 After transferring the clear supernatant to the 1cc syringe, insert the plunger of the syringe to the 1cc syringe to elute the solution through the filter unit into a pre-labeled mass spec vial.
- 23 Place a cap on the mass spec vial containing the filtered extracted tissue. The sample is ready for mass spec analysis (store the samples in -20C until ready for mass spec analysis, but we prefer to extract and run the sample in the same day if possible)