



Mar 26, 2025

## 0.01M Phosphate buffered saline (1L)

This protocol is a draft, published without a DOI.

Thomas R Clarke<sup>1</sup>

<sup>1</sup>University of Edinburgh

Duguid Lab



Thomas R Clarke

University of Edinburgh

---

OPEN  ACCESS



**Protocol Citation:** Thomas R Clarke 2025. 0.01M Phosphate buffered saline (1L). **protocols.io** <https://protocols.io/view/0-01m-phosphate-buffered-saline-1l-cimcuc2w>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** In development

**We are still developing and optimizing this protocol**

**Created:** October 31, 2022

**Last Modified:** March 26, 2025

**Protocol Integer ID:** 72068

## Abstract

Guide for making 0.01M PBS.

Used as the solution for slicing brains in the vibratome, and for transcardial perfusionsTHIS IS A DRAFT  
PROTOCOL - PLEASE EXERCISE CAUTION



- 1 Measure ~900mL distilled water 2m
- 2 Add 5 phosphate buffered saline tablets (Sigma) and stir with magnetic flea 5m
- 3 Fill up to 1000ml with distilled water 2m
- 4 Measure pH and adjust to 7.4 5m