



Oct 14, 2020

Restriction enzyme & gel electrophoresis

Hung Liang Pai¹¹Chung Shan Medical University

1

Works for me

This protocol is published without a DOI.

Hung Liang Pai

Chung Shan Medical University

PROTOCOL CITATION

Hung Liang Pai 2020. Restriction enzyme & gel electrophoresis. **protocols.io**
<https://protocols.io/view/restriction-enzyme-amp-gel-electrophoresis-bh67j9hn>

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jul 02, 2020

LAST MODIFIED

Oct 14, 2020

PROTOCOL INTEGER ID

38847

Preparation

- 1 Set the dry bath incubator, **37 °C**
- 2 Take 10x buffer from the **-20 °C fridge** to thaw.

Protocol-Restriction enzyme

- 3 To add **500 ng** DNA for each sample into each eppendorf, calculate the required DNA volume for each sample. The volume is decided by the formula below:
Volume(sample) = **500 ng** / concentration measured by spectrophotometer (ng/uL)
- 4 Calculate the required ddH₂O volume for each sample. The volume is decided by the formula below:
Volume(ddH₂O) = 20uL - 2uL - 0.5uL - 0.5uL - Volume(sample) uL
- 5 Take sufficient eppendorfs and mark them for each sample. First, add ddH₂O based on previous calculation in each eppendorf. Next, add DNA samples based on previous calculation in each eppendorf.
- 6 Take an eppendorf and mark as CT(cocktail). Add 10x buffer to the eppendorf. The volume is decided by the formula below:
Volume(10x buffer) = 2uL x (the number of samples + 1)

- 7 Take the determined restriction enzyme from the **-20 °C fridge** . Add RE to the previous eppendorf respectively.
Remember to put restriction enzyme **On ice** during the whole process and take RE back to the **-20 °C fridge** right after the process ASAP. The volume is decided by the formula below:

$$\text{Volume(RE 1)} = 0.5\mu\text{L} \times (\text{the number of samples} + 1)$$

$$\text{Volume(RE 2)} = 0.5\mu\text{L} \times (\text{the number of samples} + 1)$$
- 8 Add **3 µl** of CT in each eppendorf. Vortex and spin down each eppendorf to make the liquid stay at the bottom.
- 9 Put each eppendorf in the **37 °C** dry bath incubator for **02:00:00** or even more.

Protocol-Gel electrophoresis

- 10 Put a gel into the electrophoresis tank.
- 11 Take a duran youtility. Add **3.5 mL** TAE and **350 mL** ddH2O into it and then mix it. Put it into the electrophoresis tank. Make sure the gel is soaked into the liquid completely.
- 12 Take **5 µl** of marker (1kb) and load it in the first well.
- 13 Take loading dye and drop on the parafilm separately. The volume is decided by the formula below:

$$\text{Volume(loading dye)} = 1\mu\text{L} \times (\text{the number of samples})$$
- 14 Take **5 µl** of sample reacted with RE for **02:00:00** . After pipetting with **1 µl** loading dye, add it into each well.
Repeat this step until all the samples are filled into each well.
- 15 Connect wires with the power supply. Set time for **00:30:00** and click "start" button. Make sure there are bubbles appearing next to anode and cathode.
- 16 After gel electrophoresis, put the gel into gel reading machine and report the result.