

Aug 08, 2024

Immunoblotting using precast gels

DOI

dx.doi.org/10.17504/protocols.io.14egn6j8pl5d/v1

Francisco Bustos^{1,2}

¹Pediatrics and Rare Diseases Group, Sanford Research, Sioux Falls, SD, USA;

²Department of Pediatrics, Sanford School of Medicine, University of South Dakota, Sioux Falls, SD, USA



Francisco Bustos

Sanford Research

OPEN  ACCESS



DOI: **dx.doi.org/10.17504/protocols.io.14egn6j8pl5d/v1**

Protocol Citation: Francisco Bustos 2024. Immunoblotting using precast gels. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.14egn6j8pl5d/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 07, 2024

Last Modified: August 08, 2024

Protocol Integer ID: 104903

Abstract

Immunoblotting is a key technique to visualize changes in protein levels upon treatments. This technique can be challenging and established procedures are required to ensure reproducibility. Here we present our optimized protocol for immunoblotting of protein samples using precast gels and semiwet transfer. This protocol can be used to analyze samples from cell extracts and from in vitro reactions.



Protocol materials





- ☒ Polysorbate 20 **Thermo Fisher Scientific Catalog #233360010** In [2 steps](#)
- ☒ MES-SDS Buffer (20X) **Boston Bioproducts Catalog #BP-177** Step 7
- ☒ Tris Base **Merck MilliporeSigma (Sigma-Aldrich) Catalog #648310** Step 14
- ☒ Gel Knife **Thermo Fisher Catalog #EI9010** Step 19
- ☒ Amersham™ Protran® Western blotting membranes nitrocellulose **Merck Catalog #GE10600041** Step 16
- ☒ Ponceau S solution **Merck MilliporeSigma (Sigma-Aldrich) Catalog #Ponceau S solution** Step 26
- ☒ Glycine **Fisher Scientific Catalog #BP381-500** In [2 steps](#)
- ☒ PowerPac™ HC Power Supply **Bio-Rad Laboratories Catalog #1645052** In [2 steps](#)
- ☒ Grade 3MM Chr Blotting Paper, sheet, 46 × 57 cm **Cytiva Catalog #3030-917** Step 17
- ☒ Carnation Instant Non Fat Dry Milk **Amazon.com Catalog #12428935** Step 29
- ☒ Sodium dodecyl sulfate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #62862** Step 39
- ☒ Fisher BioReagents™ Bovine Serum Albumin Heat Shock Treated **Fisher Scientific Catalog #BP1600-100** Step 31
- ☒ Western blot boxes, 3 1/2 x 2 9/16 x 1in. 8.9 x 6.5 x 2.5cm **Fisher Scientific Catalog #NC1126730** Step 30
- ☒ PageRuler® Plus Prestained Protein Ladder, 10 to 250 kDa **Thermo Fisher Catalog #26620** Step 10
- ☒ NUPAGE LDS sample buffer (4x) **Thermo Fisher Scientific Catalog #NP0007** Step 1
- ☒ NuPAGE Sample Reducing Agent (10X) **Thermo Fisher Scientific Catalog #NP0009** Step 1
- ☒ XCell SureLock® Mini-Cell and XCell IITM Blot Module Kit **Invitrogen - Thermo Fisher Catalog #EI0002** Step 9
- ☒ NuPAGE® 4-12% Bis-Tris Protein Gels, 1.0 mm, 12-well **Thermo Fisher Catalog #NP0322BOX** Step 8
- ☒ SNAP id® 2.0 Blot Roller **Merck MilliporeSigma (Sigma-Aldrich) Catalog #SNAP2RL** Step 21
- ☒ Immobilon-P PVDF Membrane, 0.45um, roll **Merck MilliporeSigma (Sigma-Aldrich) Catalog #IPVH00010** Step 16
- ☒ Sponge Pad for Blotting **Thermo Fisher Catalog #EI9052** Step 18
- ☒ Methanol, 99.9%, for analysis **Thermo Fisher Scientific Catalog #176840025** Step 15


Safety warnings

! All steps must be performed using personal protective equipment including gloves and eye protection.


Sample preparation

- 1

Mix  and  in a 7:3 ratio (for example:  LDS sample buffer and  of Reducing Agent).


- 2



In a 1.5 mL microcentrifuge tube add the following reagents in order: LDS sample buffer/Reducing Agent mix, IPMS lysis buffer, protein sample. For example.




Sample volume	LDS sample buffer/Reducing Agent
10 or 15 microliters	5 microliters
20	7
30	10

- 3

Poke hole in top of 1.5 mL microcentrifuge tube with a hypodermic needle.
- 4

Boil samples at  for  using a ThermoMixer F1.5.



Equipment

Eppendorf ThermoMixer F1.5	NAME
ThermoMixer	TYPE
Eppendorf	BRAND
5384000020	SKU
https://www.eppendorf.com/us-en/Products/Temperature-Control-and-Mixing/Instruments/Eppendorf-ThermoMixerF-p-PF-133437	LINK

5 Centrifuge for 00:00:10 in a Mini Centrifuge at Room temperature

10s

Equipment

Fisherbrand Mini-Centrifuge

NAME

Centrifuge

TYPE

Fisherbrand

BRAND

12-006-901

SKU

<https://www.fishersci.com/shop/products/fisherbrand-standard-mini-centrifuge/12006901#?keyword=>

LINK

6 Samples can be stored at -20 °C until the day of running SDS-PAGE electrophoresis.



SDS-PAGE Electrophoresis

7 Prepare 500 mL 1X MES-SDS Buffer (50 millimolar (mM) MES, 50 millimolar (mM) Tris Base, 0.1 Mass / % volume SDS, 1 Mass / % volume EDTA, pH 7.3) by diluting MES-SDS Buffer (20X) **Boston Bioproducts Catalog #BP-177** stock 1:20 in Milli-Q water.





8 Unpack a NuPAGE®; 4-12% Bis-Tris Protein Gels, 1.0 mm, 12-well **Thermo Fisher Catalog #NP0322BOX** and remove protective tape and comb.

9 Assemble gel running tank XCell SureLock® Mini-Cell and XCell IITM Blot Module Kit **Invitrogen - Thermo Fisher Catalog #EI0002** by placing gel in the tank in front of the buffer core and a second gel or a buffer dam behind the




buffer core. Lock the gel tension wedge in place. Add 1X MES-SDS buffer between the gels and outside the buffer core.


- 10 Load  8 μL protein molecular weight marker



 PageRuler[®]; Plus Prestained Protein Ladder, 10 to 250 kDa **Thermo Fisher** **Catalog #26620**

. If needed load  2 μL at the end of the gel.

- 11 Load samples to a maximum volume of  20 μL .

- 12 Place gel tank lid and connect to power supply



 PowerPac[™] HC Power Supply **Bio-Rad Laboratories** **Catalog #1645052** matching colors (for positive and negative).



- 13 Run gel at 150 V for  02:00:00 at  Room temperature or as needed depending on the protein molecular weight and desired separation.



2h



Transfer

- 14 Prepare 12X transfer buffer:

- 58 g Tris base  Tris Base **Merck MilliporeSigma (Sigma-Aldrich)** **Catalog #648310**
- 190 glycine  Glycine **Fisher Scientific** **Catalog #BP381-500**
- Milli-Q water up to  2 L


- 15 Prepare 2.5 L of 1X transfer buffer ( 48 millimolar (mM) Tris,  39 millimolar (mM)

Glycine,  20 % (v/v) Methanol) by mixing  212 mL 12X Transfer buffer with

 1788 mL milli-Q water and  500 mL Methanol

 Methanol, 99.9%, for analysis **Thermo Fisher Scientific** **Catalog #176840025** .

- 16 Cut

 Immobilon-P PVDF Membrane, 0.45um, roll **Merck MilliporeSigma (Sigma-Aldrich)** **Catalog #IPVH00010**

or

1m




Amersham™ Protran® Western blotting membranes
nitrocellulose **Merck Catalog #GE10600041**


membrane to the required gel size (usually 9x7 cm) and activate PVDF membrane by incubating in methanol for at least 00:01:00 .

Note



Do not put Nitrocellulose membrane in methanol

- 17 Cut Grade 3MM Chr Blotting Paper, sheet, 46 × 57 cm **Cytiva Catalog #3030-917** to the required gel size (usually 9x7 cm).
- 18 Prepare Sponge Pad for Blotting **Thermo Fisher Catalog #EI9052** . If dirty, boil in warm tap water in the microwave for 00:05:00 . Soak sponge pads in 1X transfer buffer in a plastic tray. 5m
- 19 Remove the gel from the cassette by separating the cassette plates with a Gel Knife **Thermo Fisher Catalog #EI9010** . Using the gel knife, cut off the bottom part and wells of the gel as needed.
- 20 Place gel in a plastic tray with 1X transfer buffer.
- 21 Make transfer sandwich inside the cathode core of the blot module in the following order:
- 2 sponge Pad
 - 1 filter paper + gel (grab this from transfer buffer tray), Remove air bubbles by carefully rolling a
- SNAP id® 2.0 Blot Roller **Merck MilliporeSigma (Sigma-Aldrich) Catalog #SNAP2RL**
- on top of the gel.
- 1 presoaked PVDF or Nitrocellulose membrane. Remove air bubbles using the Blot Roller.
 - 1 filter paper. Remove air bubbles using the Blot Roller.
 - Sponges up to fill cavity (6-7 total)
- 22 Close blot module with the anode core and place this in the buffer chamber. Lock the blot module using the gel tension wedge.

23 Place lid and connect to power supply
 PowerPac™ HC Power Supply **Bio-Rad Laboratories Catalog #1645052** matching colors (for positive and negative).

24 Run gel at 35 V for  01:30:00 at  Room temperature . 1h 30m

25 Unlock blot module by releasing the tension wedge and remove from buffer chamber. Remove membrane from the sandwich and cut excess membrane using scissors.

26 Incubate membrane in 5m
 Ponceau S solution **Merck MilliporeSigma (Sigma-Aldrich) Catalog #Ponceau S solution**
 in a plastic tray for  00:05:00 with shaking in a rocking shaker.

Equipment

Rocking Shaker

Shaker

Ohaus

30391966

<https://us.ohaus.com/en-us/products/equipment/open-air-shakers/rocking-waving-shakers/shaker-rocking-2-tier-shrk07al2-us>

NAME

TYPE

BRAND

SKU

LINK

27 Wash with milli-Q water 3 times, quick washes. Image in Chemidoc MP if needed.





Equipment

ChemiDoc™ MP Imaging System

NAME

Imaging System

TYPE

Bio-rad

BRAND

12003154

SKU



Blocking

1h

- 28 Make 1X TBS-T (Tris buffer saline-Tween 20: [M] 20 millimolar (mM) Tris-HCl pH=7.4, [M] 100 millimolar (mM) NaCl, [M] 0.1 % (v/v) Tween-20).

For 1L:

- 20 mL [M] 1 Molarity (M) Tris pH=7.4
 - 20 mL [M] 5 Molarity (M) NaCl
 - 1 mL Tween-20
- Polysorbate 20 **Thermo Fisher Scientific Catalog #233360010**
- Milli-Q water up to 1 L

- 29 Make Blocking solution: 1X TBS-T 5% skimmed milk:

- 12.5 g powder skimmed milk
- Carnation Instant Non Fat Dry Milk **Amazon.com Catalog #12428935** .
- Fill up to 250 mL with 1X TBS-T.

- 30 Add Blocking solution to a Western Blot box

Western blot boxes, 3 1/2 x 2 9/16 x 1in. 8.9 x 6.5 x 2.5cm **Fisher Scientific Catalog #NC1126730**

and place the membrane in this solution. Incubate in a rocking shaker for at least

01:00:00 at Room temperature

1h



Primary antibody

1h

- 31 Prepare antibody in 1X blocking solution (for HRP detection) or in TBST

[M] 5 Mass / % volume BSA



Fisher BioReagents™ Bovine Serum Albumin Heat Shock Treated **Fisher Scientific Catalog #BP1600-100**

(for infrared detection) at the desired concentration.

- 32 Discard blocking solution from Western blot box and replace by diluted antibody. Incubate at

🌡️ 4 °C



Overnight in rocking shaker.

20h



- 33 Next day recover antibody in tube and wash membrane 3 times for



00:05:00

with 1X TBS-T

5m

Secondary antibody

1h

- 34 Prepare HRP or fluorescently labelled secondary antibody in blocking solution or TBST respectively.

- 35 Add diluted secondary antibody to Western blot box and incubate at



Room temperature

for

1h



01:00:00 in rocking shaker.

- 36 Wash membrane 4 times for



00:05:00

with 1X TBS-T

5m

- 37 Image in Chemidoc MP for HRP detection or LI-COR Odyssey FC imager for infrared detection.



Equipment

ChemiDoc™ MP Imaging System

NAME

Imaging System

TYPE

Bio-rad

BRAND

12003154

SKU



Equipment

Odyssey CLx

NAME

Imaging System

TYPE

LI-COR

BRAND

Odyssey CLx

SKU

<https://www.licor.com/bio/odyssey-clx/>

LINK



38 Membrane can be stored up to a week at 4 °C until performing stripping.



Stripping

10m

39 Prepare Stripping buffer:


- weight 15 g glycine Glycine **Fisher Scientific Catalog #BP381-500**





- Dissolve in  800 mL milli-Q water

- Adjust pH to 2.2

- add  1 g SDS


 Sodium dodecyl sulfate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #62862**

- 10 mL Tween 20  Polysorbate 20 **Thermo Fisher Scientific Catalog #233360010**

- Bring volume up to 1 L with milli-Q water


40 Incubate membrane  00:10:00 with stripping buffer in rocking shaker at

10m

 Room temperature .

41 Repeat step 37

42 Wash twice with TBST

43  [go to step #28](#) and probe with a different antibody.