



Aug 12, 2021

Single-cell RNA sequencing

Klaus H. Kaestner Lab¹, Suzanne Shapira¹¹University of Pennsylvania Perelman School of Medicine

Suzanne Shapira: HIRN-HPAP correspondence: sshapira@pennmedicine.upenn.edu

1 Works for me

 Sharedx.doi.org/10.17504/protocols.io.bwdepa3e

Human Islet Research Network



Lili Liang

ABSTRACT

Single-cell RNA sequencing (scRNA-seq) allows for transcriptional profiling of individual cells within a heterogenous sample. This protocol describes a method for performing scRNA-seq using handpicked pancreatic islets from organ donors.

Date submitted: August 13, 2020

DOI

dx.doi.org/10.17504/protocols.io.bwdepa3e

EXTERNAL LINK

<https://hpap.pmacs.upenn.edu/explore/workflow/islet-molecular-phenotyping-studies?protocol=4>

PROTOCOL CITATION

Klaus H. Kaestner Lab, Suzanne Shapira 2021. Single-cell RNA sequencing. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bwdepa3e>

KEYWORDS

null, Single-cell RNA sequencing, HPAP, HIRN

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jul 06, 2021


LAST MODIFIED

Aug 12, 2021

PROTOCOL INTEGER ID

51334

Steps in pre-processing

1. Transfer handpicked islets (approximately 5,000 IEQs) into  15 mL conical tube.

2. Add **10 mL** of 1xPBS w/o Ca²⁺, Mg²⁺ ([Rockland, MB-008](#)). Centrifuge for 2 min at RT, 180 xg. Aspirate the supernatant.
3. Add **1 mL** of warm (**37 °C**) 0.05% Trypsin ([Invitrogen, 25300054](#)) to the islets. Pipette up and down with p1000.
4. Incubate at **37 °C** for 9 min, or until cells are in single cells. Pipette up and down at t=7 min, 4 min, 2 min, 0 min.
5. Stop the trypsin reaction by adding **1 mL** of 100% FBS ([Hyclone, SH3091003](#)) to the dissociated islets and pass cells through BD FACs tube with strainer top ([Corning 352235](#))
6. Use **1 mL** of 100% FBS to rinse the tube and pass through the strainer.
7. Transfer cells to **15 mL** conical. Centrifuge 4 min, 400 xg.
8. Remove the supernatant and wash cells with PBS with 10% FBS. Centrifuge for 4 min, 400 xg.
9. Wash the cells with PBS with 10% FBS and centrifuge for 4 min, 400 xg. Remove the supernatant.
10. Count cells using a countess chamber.
11. For the scRNAseq, do the final resuspension in 10% FBS in PBS and adjust the volume to make the final suspension 1000 cells/microliter. Filter the cells one more time prior to loading them onto the 10X Genomics chip.

Links to Kits used in post-processing

- 2
 1. Older HPAP samples (specifically donors HPAP-001 to HPAP-019) were processed using [C1 Single-Cell mRNA Seq HT IFC and Reagent Kit v2](#) (product ID: 101-4964) which has been discontinued.
 2. All current samples (HPAP-019 onwards) were processed using the [Chromium Single Cell 3' Reagent Kit](#).
 3. Unfortunately, since 10x Genomics is going to stop manufacturing the above kit, we will be using the new kit [Chromium Next GEM Single Cell 3' Reagent Kits v3.1](#). For this protocol we target a 5000 cell recovery.