



DNA extraction from fecal samples V.2

Yoshiyuki Matsuo¹

¹Kansai Medical University

Version 2 ▾

Aug 03, 2021

1 Works for me

Share

dx.doi.org/10.17504/protocols.io.bs9inh4e



Yoshiyuki Matsuo

ABSTRACT

DNA extraction from fecal samples

DOI

dx.doi.org/10.17504/protocols.io.bs9inh4e

PROTOCOL CITATION

Yoshiyuki Matsuo 2021. DNA extraction from fecal samples. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.bs9inh4e>

Version created by Yoshiyuki Matsuo

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Mar 12, 2021

LAST MODIFIED

Aug 03, 2021

PROTOCOL INTEGER ID

48138

Materials

- 1
 - Phosphate-buffered saline (PBS)
 - EZ-beads (Promega, AMR76813M): 2 mL tube containing 0.2 mm ceramic (zirconium dioxide) spheres and one 5 mm zirconia bead
 - Maxwell RSC Blood DNA Kit (Promega, AS1400)
 - MN Bead Tube Holder (MACHEREY-NAGEL, 740469): Rubber-foam adapter for processing bead tubes with Vortex-Genie instrument

Instrument

- 2
 - Maxwell RSC (Promega, AS4500): Automated nucleic acid purification platform
 - Micro Smash (TOMY, MS-100): Beads cell disrupter
 - Vortex-Genie Mixer
 - Heating block
 - Microcentrifuge

Preparation of fecal samples

- 3 Place 50-100 mg of fecal sample into tube.

- 4 Add 1 ml of PBS per 100 mg of feces.
- 5 Mix thoroughly by vortexing or pipetting.
- 6 Allow the sample to stand for 2 min to sediment large debris.
- 7 Transfer 300 µl of the suspension to 1.5 mL tube.
- 8 Centrifuge at 18,000 x g for 3 min.
- 9 Discard the supernatant.
- 10 Resuspend the pellet (~30 mg of feces) in 300 µl of PBS.
- 11 Incubate at 70°C for 10 min and cool to room temperature.

Mechanical cell disruption by beat beating

- 12 Transfer 300 µl of the suspension to EZ-beads tube.
- 13 Lyse cells either by using disruption device (13.1) or vortex mixer (13.2).
 - 13.1 Place the EZ-beads tube in Micro Smash instrument and disrupt cells at 2,500 rpm for 2 min.
 - 13.2 Place the EZ-beads tube on MN Bead Tube Holder attached to Vortex-Genie mixer and vortex for 5 min at maximum speed.
- 14 Briefly spin the tube.

Automated DNA extraction using Maxwell RSC Blood DNA Kit

- 15 Add 300 µl of Lysis Buffer and 30 µl of Proteinase K Solution to the sample in EZ-beads tube.
- 16 Mix thoroughly by vortexing for 10 sec.
- 17 Incubate at 56°C for 20 min.
- 18 Briefly spin the tube.
- 19 Transfer the supernatant (~600 µl) to 1.5 mL tube.
- 20 Centrifuge at 18,000 x g for 3 min.
- 21 Transfer the cleared lysate to Maxwell RSC Cartridge.
- 22 Add 50 µl of Elution Buffer to elution tube.
- 23 Start the extraction run following the manufacturer's instructions.