



Jul 12, 2022

# Striatal dopamine measurement through HPLC

Pranay Srivastava<sup>1</sup>, Waijiao Cai<sup>2</sup>, Xiqun Chen<sup>1</sup>

<sup>1</sup>MassGeneral Institute for Neurodegenerative Disease, Department of Neurology, Massachusetts General Hospital, Harvard Medical School, Boston, USA Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, Towson, MD, USA;

<sup>2</sup>MassGeneral Institute for Neurodegenerative Disease, Department of Neurology, Massachusetts General Hospital, Harvard Medical School, Boston, USA Department of Integrative Medicine, HuaShan Hospital, Institutes of Integrative Medicine, Fudan University, Shanghai, China

1 Works for me

Share

[dx.doi.org/10.17504/protocols.io.dm6gpbjdplzp/v1](https://dx.doi.org/10.17504/protocols.io.dm6gpbjdplzp/v1)

Pranay Srivastava

## ABSTRACT

Protocol for striatal dopamine measurement in mouse brain

## DOI

[dx.doi.org/10.17504/protocols.io.dm6gpbjdplzp/v1](https://dx.doi.org/10.17504/protocols.io.dm6gpbjdplzp/v1)

## PROTOCOL CITATION

Pranay Srivastava, Waijiao Cai, Xiqun Chen 2022. Striatal dopamine measurement through HPLC. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.dm6gpbjdplzp/v1>



## KEYWORDS

ASAPCRN

## LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Feb 18, 2022

## LAST MODIFIED

Jul 12, 2022

MATERIALS TEXT

**PE Solution Recipe: (250 mL):**

Polished ddH <sub>2</sub> O-	249.06 mL
0.1mM EDTA (From 0.5 mM stock)-	50 uL
1uM DHBA (From 10mM stock)-	25 uL
50mM Phosphoric Acid (From 85% Fisher #A260 HPLC grade)-	855.3 uL

**Catecholamine Standard Prep:**

Prepare 10 mM Stock in water or PE solution. (May be frozen as aliquots at -80C)

-Dopamine = 18.96 mg to 10 mL

(Dopamine hydrochloride, Sigma #H8502, MW 189.64)

-DOPAC = 16.82 mg to 10 mL

(3,4 Dihydroxyphenylacetic Acid, Acros #22560050, MW 168.15)

-HVA = 18.22 mg to 10 mL

(Homovanillic acid, MilliporeSigma #69673-25MG, MW 182.17)

-DHBA= 22.01 mg to 10 mL

(3,4 Dihydroxybenzylamine Hydrobromide, Sigma #858781, MW 220.06)

\*\*\*(*INTERNAL STANDARD IN PE SOLUTION ONLY*)

Costar® Spin-X® Centrifuge Tube Filters, 0.22 µm Pore CA Membrane, Nonsterile - Cat#8161

- 1 Weigh the striatal tissue  
Note - make sure to keep it frozen – can keep it on dry ice for the process
- 2 Add 20x PE buffer supplemented with an internal standard solution to tissue.  
(eg. 10 mg sample tissue would receive 200 uL of PE solution)
- 3 Use a Teflon pestle to homogenize the tissue into the solution by crushing it against the walls of the tube.  
Note - Do not sonicate as it may oxidize our analyte)



Note - Keep homogenized samples on ice until centrifugation

Note - If using a hand homogenizer, note that the machine produces heat during prolonged use which can degrade the analytes. Allow time to cool between samples.

Alternatively, can use Dounce homogenizers


4 Centrifuge at  **14000 rpm, 4°C, 00:15:00**

15m

5 Remove supernatant and transfer to Costar Spinx (#8161) 0.22 - MCA filter tube  
Centrifuge for  **00:05:00** at  **14.000 rpm, 4°C, 00:05:00** or until all liquid has passed through the filter.

10m

Note - Filters have a max volume of 500uL, and may have to perform multiple spins/samples if the supernatant volume is >500ul.

6 The resulting solution can be run on HPLC or stored at  **-80 °C**.

7 Running a standard - HPLC

Perform serial dilution of the stocks as follows to create the standard curve.

Level	[ng/mL]	uL stock	uL diluent
1	1000	1 from each 10 mM stock (DA/DOPAC/HVA)	997
2	600	429.8 from level 1	286.59
3	181.28	216.47 from level 2	286.59
4	54.77	216.47 from level 3	286.59
5	16.55	216.47 from level 4	286.59
6	5	216.47 from level 5	286.59

8 Load at least 50 uL into each vial to run the standard curve in triplicate (on HPLC).

Run the samples