

Sep 16, 2020

mRNA Synthesis and Transfection

In 1 collection

Yingchao Xue^{1,2}, Xiping Zhan³, Shisheng Sun⁴, Senthilkumar S. Karuppagounder^{5,6,7}, Shuli Xia^{2,5}, Valina L Dawson^{5,6,7,8,9}, Ted M Dawson^{5,6,7,8,10}, John Laterra^{2,5,8,11}, Jianmin Zhang¹, Mingyao Ying^{2,5}

¹Department of Immunology, Research Center on Pediatric Development and Diseases, Institute of Basic Medical Sciences, Chinese Academy of Medical Sciences and School of Basic Medicine, Peking Union Medical College, State Key Laboratory of Medical Molecular Biology;

²Hugo W. Moser Research Institute at Kennedy Krieger; ³Department of Physiology and Biophysics, Howard University;

⁴College of Life Sciences, Northwest University; ⁵Department of Neurology, Johns Hopkins University School of Medicine;

⁶Neuroregeneration and Stem Cell Programs, Institute for Cell Engineering, Johns Hopkins University School of Medicine;

⁷Adrienne Helis Malvin Medical Research Foundation; ⁸Department of Neuroscience, Johns Hopkins University School of Medicine;

⁹Department of Physiology, Johns Hopkins University School of Medicine;

¹⁰Department of Pharmacology and Molecular Sciences, Johns Hopkins University School of Medicine;

¹¹Department of Oncology, Johns Hopkins University School of Medicine

1 Works for me This protocol is published without a DOI.

Neurodegeneration Method Development Community
Tech. support email: ndcn-help@chanzuckerberg.com



Anita Broellochs
protocols.io

ABSTRACT

This protocol explains the mRNA synthesis and transfection of lines ND1014, N1, and ND27760 from *Synthetic mRNAs Drive Highly Efficient iPS Cell Differentiation to Dopaminergic Neurons*.

EXTERNAL LINK

<https://www.ncbi.nlm.nih.gov/pmc/articles/PMC6344911/>

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Synthetic mRNAs Drive Highly Efficient iPS Cell Differentiation to Dopaminergic Neurons. Xue Y, Zhan X, Sun S, Karuppagounder SS, Xia S, Dawson VL, Dawson TM, Laterra J, Zhang J, Ying M. Stem Cells Transl Med. 2019 Feb;8(2):112-123. doi: 10.1002/sctm.18-0036. Epub 2018 Nov 1. PMID: 30387318

PROTOCOL CITATION

Yingchao Xue, Xiping Zhan, Shisheng Sun, Senthilkumar S. Karuppagounder, Shuli Xia, Valina L Dawson, Ted M Dawson, John Laterra, Jianmin Zhang, Mingyao Ying 2020. mRNA Synthesis and Transfection.

protocols.io

<https://protocols.io/view/mrna-synthesis-and-transfection-9u5h6y6>

MANUSCRIPT CITATION please remember to cite the following publication along with this protocol

Synthetic mRNAs Drive Highly Efficient iPS Cell Differentiation to Dopaminergic Neurons. Xue Y, Zhan X, Sun S, Karuppagounder SS, Xia S, Dawson VL, Dawson TM, Laterra J, Zhang J, Ying M. Stem Cells Transl Med. 2019 Feb;8(2):112-123. doi: 10.1002/sctm.18-0036. Epub 2018 Nov 1. PMID: 30387318

EXTERNAL LINK

<https://www.ncbi.nlm.nih.gov/pmc/articles/PMC6344911/>

COLLECTIONS ①



Protocols for Synthetic mRNAs Drive Highly Efficient iPS Cell Differentiation to Dopaminergic Neurons

KEYWORDS

ND1014, N1, ND27760, ipsc, SNCA, Atoh2, Ngn2, HiScribe T7 ARCA

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Nov 27, 2019

LAST MODIFIED

Sep 16, 2020

OWNERSHIP HISTORY

Nov 27, 2019  Liz Brydon [Protocols.io](#)

Sep 16, 2020  Anita Broellochs [protocols.io](#)

PROTOCOL INTEGER ID

30333

PARENT PROTOCOLS

Part of collection

[Protocols for Synthetic mRNAs Drive Highly Efficient iPS Cell Differentiation to Dopaminergic Neurons](#)

GUIDELINES

Please follow the National Institutes of Health guidelines.

MATERIALS

NAME	CATALOG #	VENDOR
HiScribe T7 ARCA mRNA Kit - 20 rxns	E2065S	New England Biolabs
MEGAclean[®] Transcription Clean-Up Kit	AM1908	Thermo Fisher
Lipofectamine[®] MessengerMAX[®] Transfection Reagent	LMRNA001	Thermo Fisher


SAFETY WARNINGS

Please refer to the Safety Data Sheets (SDS) for safety and environmental hazards.

BEFORE STARTING

Obtain approval to work with human stem cells from an appropriate Institutional Review Board.

mRNA Synthesis

- 1 Clone coding sequence of Atoh1 into a vector containing the T7 promoter and poly(A) tail for in vitro transcription.
- 2 Clone coding sequence of Ngn2 into a vector containing the T7 promoter and poly(A) tail for in vitro transcription.
- 3 
Linearize Atoh1 and Ngn2 vectors.
- 4 Subject linearized Atoh1 and Ngn2 vectors to mRNA synthesis using the HiScribe T7 ARCA mRNA Kit.

5 Purify mRNAs using the MEGAclean Transcription Clean-Up Kit.

Transfect mRNA

6 Transfect using a cationic lipid (e.g. Stem-In).

6.1

For each well of a 12-well plate, incubate **0.25 µg mRNA** with **1.5 µl lipid** in **50 µl PBS** for **00:15:00**.

6.2 Add DNA:lipid complexes to cells.