

Aug 27, 2024

## Bioluminescence Imaging

DOI

**dx.doi.org/10.17504/protocols.io.ewov191r2lr2/v1**

Santiago Unda<sup>1,2</sup>, Michael G. Kaplitt<sup>1,2</sup>

<sup>1</sup>Department of Neurological Surgery, Weill Cornell Medical College, New York, NY 10065;

<sup>2</sup>Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD 20815, USA



**Eileen Ruth Torres**

Weill Cornell Medicine

---

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.ewov191r2lr2/v1](https://dx.doi.org/10.17504/protocols.io.ewov191r2lr2/v1)

**Protocol Citation:** Santiago Unda, Michael G. Kaplitt 2024. Bioluminescence Imaging. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.ewov191r2lr2/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** August 27, 2024

**Last Modified:** August 27, 2024

**Protocol Integer ID:** 106521

**Keywords:** ASAPCRN, imaging, bioluminescence

**Funders Acknowledgement:**

**Aligning Science Across**

**Parkinson's**

**Grant ID:** 020608



## Abstract

This protocol describes a bioluminescence stepwise imaging protocol in mice with AAV-luciferase expression in the vagal afferents using the Ivis Lumina series III imaging system.

## Materials


- D-luciferin potassium salt (Xenogen)
- Isoflurane, USP (Covetrus)



## Preparation and In vivo imaging

- 1 Induce anesthesia with isoflurane 2-3%.
- 2 Shave the neck, thoracic and abdominal area.
- 3 Prepare D-luciferin potassium salt (Xenogen) stock prepared in phosphate buffered saline (PBS) at 15mg/mL.
- 4 Inject D-luciferin solution (150mg/kg) intraperitoneal.
- 5 Image mice in supine position at 5-, 10- and 15-minutes post-luciferin injection.

## Software protocol and analysis

- 6 Living imaging software (Ivis imaging) was configured to luminescent imaging mode with an exposure time of  00:00:00 and 5-minute delay.
- 7 Post-acquisition, the total flux of photons/radiance per second (s) in each pixel in the ROI area, in a square centimeter (cm<sup>2</sup>) of the tissue. This data is then multiplied by one steradian (sr), which result in average radiance displayed as (photons/s/cm<sup>2</sup>/sr).