



JUL 30, 2023

🌐 Measurement of GCase activity in lysosomes in live cells

rachel.bates¹

¹UCL



rachel.bates

ABSTRACT

OPEN  ACCESS



Protocol Citation: rachel.bates 2023. Measurement of GCase activity in lysosomes in live cells. **protocols.io** <https://protocols.io/view/measurement-of-gcase-activity-in-lysosomes-in-live-cxwzxp6>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Jul 30, 2023

Last Modified: Jul 30, 2023

PROTOCOL integer ID:
85689

Keywords: ASAPCRN

Assay uses the substrate PFB-FDGluc (5-(Pentafluorobenzoylamino)Fluorescein Di-β-D-Glucopyranoside) from ThermoScientific (Cat# P11947). Substrate is supposed to be taken up by late endosomes and lysosomes only and fluoresces green when cleaved by lysosomal GCase.

Described by Mazzulli et al., J Neurosci. 2016 Jul 20;36(29):7693-706.

MATERIALS

PFB-FDGluc (5-(Pentafluorobenzoylamino)Fluorescein Di-β-D-Glucopyranoside)

Catalog number: P11947

Opti-MEM™ I Reduced Serum Medium

Catalog number: 31985062

- 1 Resuspend 5mg of substrate in 100 ul methanol to get 50mg/mL stock. Do not expose to light. Store at -20C
- 2 Culture cells in 48 well plate. Should be over 70% confluent to ensure a robust signal. If too confluent activity will plateau rapidly. 1d
- 3 For control wells add 10uM CBE or 100uM Bafilomycin overnight.
- 4 Wash cells once carefully with prewarmed PBS. 5m
- 5 Working substrate solution is 400 mg/ml (I don't get signal with lower concentrations). Prepare in OPTIMEM prewarmed to 37 °C: 125 ml/well. Therefore 1 ul of 50mg/ml stock substrate per 125 ml.
- 6 Add 125ul substrate (400ug/ml substrate (400 mg/ml in OPTIMEM) per well. Put in 37 °C incubator for 1 h. 1h
- 7 Aspirate substrate. Wash cells three times with 250 ml PBS (37 °C). 10m
- 8 ADD 125UL ml OPTIMEM (37 °C). Measure t=0 minutes on plate reader: Ex, 488nm, Em, 520 nm. Return cells to incubator.
- 9 Measure fluorescence every 20-30 minutes in green channel (excitation/emission maxima ~492/516 nm) for up to 3 hours. Keep checking cells are alive with microscope. 3h

10 Optional: At end of experiment carefully aspirate OPTIMEM. Lyse cells in wells with 200 μ l 1% TX-100 in PBS on ice for 15 minutes. Pellet debris at 17,000 $\times g$, 10 minutes, 4°C. Measure protein concentration in supernatant with BCA assay.

11