



MAR 05, 2024

Transforming pJC8 into HB101

Bonnie Evans^{1,2}

¹MRC Laboratory of Medical Sciences; ²Imperial College London



Bonnie Evans

MRC Laboratory of Medical Sciences, Imperial College London

ABSTRACT

Taken from *Mix and Go E. coli* Transformation Kit (see attachment).
Transforming *E. coli* HB101 with pJC8-empty cosmid to use as a control strain.

ATTACHMENTS

[_t3001_t3002_mix_go_e_c
oli_transformation_kit_buff
er_set.pdf](#)

MATERIALS

Mix and Go E. coli Transformation Kit (Zymo Research, T3001/T3002)

OPEN ACCESS



DOI:

dx.doi.org/10.17504/protocols.io.n2bvj3kkwlk5/v1

Protocol Citation: Bonnie Evans 2024. Transforming pJC8 into HB101. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.n2bvj3kkwlk5/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Mar 05, 2024

Contents of SOB medium:

20g Bacto Tryptone

5g Yeast Extract

0.58g Sodium Chloride

0.19g Potassium Chloride

1L sterile H2O

Before starting

1 Prepare tetracycline LB agar plates

Protocol



NAME

Making tetracycline LB agar plates

CREATED BY

Bonnie Evans

PREVIEW

2 Warm plates in 37 °C incubator

Note

Chilled plates will decrease *Mix & Go* cell transformation efficiency.

3 Prepare SOC medium

3.1 Get SOB medium and 20 % glucose from media kitchen

3.2 Add 2 mL 20 % glucose to 100 ml SOB medium

Transformation

4 Isolate pJC8 cosmid from *E. coli* DH5α

Protocol



NAME

Qiagen QIAprep Spin Miniprep for cosmids from metagenomic library

CREATED BY

Bonnie Evans

PREVIEW

5 Thaw 100 uL aliquot of *E. coli* HB101 *Mix & Go* competent cells on ice

6 Add 5 uL pJC8 cosmid

7 Mix by tapping the tube and shaking downwards once

8 Place on ice for 5 minutes

9 Add 400 uL SOC medium

10 Incubate at 37 °C with shaking at 200 rpm for 1 hour

Note

An outgrowth in SOC medium is required for efficient transformation when selecting with tetracycline.

11 Pipette all (505 uL) of the transformation mix onto a tetracycline LB agar plate and spread using plating beads.

12 Dry under laminar flow for a few minutes

13 Incubate plate at 37 °C overnight

Expected result

If the transformation is successful, you will have individual colonies.